



K13902, K13992 CLEVELAND OPEN-CUP FLASH POINT TESTER

OPERATION AND INSTRUCTION MANUAL

REV E

Koehler Instrument Company, Inc.

1595 Sycamore Avenue • Bohemia, New York 11716-1796 • USA

Toll Free: 1-800-878-9070 (US only) • Tel: +1 631 589 3800 • Fax: +1 631 589 3815

<http://www.koehlerinstrument.com> • e-mail: info@koehlerinstrument.com

Petroleum Testing & Analysis Instrumentation • Custom Design & Manufacturing

CERTIFICATE OF CONFORMANCE

Manual Cleveland Open Cup Flash Point Tester K139XX

This certificate verifies that part number K139XX, Manual Cleveland Open Cup Flash Point Tester, was manufactured in conformance with the applicable standards set forth in this certification.

Specifications:

ASTM D92
ASTM D6074
AASHTO T48
ASN Z-11.6
IP 36
ISO 2592
DIN 51376
FTM 791-1103
FTM 141-4294

This unit is tested before it leaves the factory, to ensure total functionality and compliance to the above specifications and ASTM standards. Test and inspection records are on file for verification.



Jesse Kelly
Application Engineer
Koehler Instrument Company



EC Declaration of conformity

Koehler Instrument Company, Inc.
of 1595 Sycamore Av., Bohemia, New York USA

We declare that the product listed below meets all basic requirements in accordance with the following Directive(s) by design, type, and version placed upon the market by us.

2004/108/EC The Electromagnetic Compatibility Directive
2006/42/EC The Machinery Directive by way of the Low-Voltage directive 2006/95/EC

And hereby declare that:

Equipment : Manual Cleveland Open Cup Tester
Model Numbers: K13992
Date: 04/22/2013

Qualifications:


This product may only be used in a professional laboratory setting by authorized personnel following the instruction handbook.

and

This product declaration is valid for unmodified equipment when installed and operated by authorized personnel following the instruction handbook.

Conforms to the following standards:

Safety	Low-Voltage directive 2006/95/EC
EN 61010-1:2010	Safety Requirements for electrical equipment for measurement, control and laboratory use; by engineering design and risk review and by meeting the requirements of Hi-Pot Test (1900 VAC, 60 sec.) as detailed in the product's technical documentation.
EMC	Meets the essential requirements of EMC Directive 2004/108/EC by engineering design review and by meeting the requirements of Conducted Emissions Test for Group 1 Class A as detailed in the product's technical documentation.
EN 55011:2007	



James R. Ball
Chief Technology Officer

1595 Sycamore Av.
Bohemia, NY 11716
United States of America
April 22, 2013

www.koehlerinstrument.com

631-589-3800

WEEE Directive Compliance Statement

Background

The goal of the WEEE Directive is to encourage design of environment-friendly products that increase reuse, recycling and other forms of recovery to reduce waste streams and applies to listed Electronic and Electrical Equipment (EEE) and Koehler's equipment falls broadly into Appendix 1A; Section 9 Monitoring and Control Equipment: Measuring, weighing or adjusting appliances for household or as laboratory equipment.

Any associated non-embedded equipment such as Lighting (Saybolt Color) and PCs/Printers also fall under WEEE. If provided with an order these ancillary items must be WEEE compliant. For these and other reasons (printer cartridges are regionalized) the equipment must be supplied through a third party supplier in Europe.

The WEEE Directive applies to electrical and electronic equipment falling under the categories set out in Annex IA provided that the equipment concerned is not part of another type of equipment that does not fall within the scope of this Directive. Annex IB contains a list of products which fall under the categories set out in Annex IA.

<http://eur-lex.europa.eu/LexUriServ/LexUriServ.do?uri=OJ:L:2003:037:0024:0038:en:PDF>

We do not qualify for any of the 10 exemption categories.
<http://www.dpa-system.dk/en/WEEE/Products/Exemptions>

Professional use

For equipment defined for 'professional use' local authorities have no role to play. Producers and importers are basically responsible for collection of WEEE recyclables from the professional user and for subsequent management. A separate statement is given cataloging the items that require separation from the equipment along with basic information on subsequent processing or recycling prior to disposal of the equipment.

<http://www.dpa-system.dk/en/WEEE/Products/Private-or-professional-use>

Responsibility for Registration and Annual Reporting:

Koehler will not sell directly to end users in the EU and so has no responsibility to register within each EU state and to make annual reports. Koehler declares that this responsibility is born by the importer who is the first level of the distribution chain and is subject to producer responsibility. We will communicate this in writing to our distributor/importers in the EU stating they are responsible to satisfy WEEE registration and reporting requirements in the EU states where they conduct sales activities.

It is illegal to market electrical and electronic equipment covered by producer responsibility without being registered.

<http://www.dpa-system.dk/en/WEEE/Producers/Whoissubjecttoproducerresponsibility>

Product Design

Koehler's designs allow for complete disassembly to a modular level which usually allows for standard recycling. A qualified refrigeration system technician must be consulted when disassembling and de-commissioning any equipment with refrigeration systems.

Koehler's scientific testing equipment is robustly designed to function over a long service life and are typically repaired many times over the course of years rather than being replaced. We believe that re-use and refurbishment is the very best form of re-cycling.

All batteries must be readily removable not soldered in place.

Recycling instructions

In the event that replacement becomes necessary, we will include instructions, particularized to each instrument that informs the customer of their recycling responsibilities and giving them guidance in doing this. All Koehler equipment has been placed on the market since 13th August 2005 and so Koehler is defined as a "new WEEE producer". As such we must provide information on refurbishment, treatment, and re-use.

Our instrument manual will include this compliance statement and indicate that any collection of materials will be handled by their authorized distributor. In the event that the distributor is unreachable or is no longer a distributor for Koehler Instrument, Co., other arrangements may be made including accepting the materials directly.

Recycling is free of charge. Shipping is the responsibility of the end users. Whether shipping to a distributor or to Koehler directly, safe, properly declared, and labeled packaging and shipping expenses are the sole responsibility of the end user.

WEEE Marking



Since Koehler products are subject to the WEEE Directive we must display the WEEE symbol shown above in accordance with European Standard EN 50419 on the equipment. It must be indelible, at least 5mm in height, and clearly legible. If the equipment is too small the mark must be in the product literature, guarantee certificate, or on the packaging. Rules on marking are established in section 49 of the WEEE Order.

Koehler Instrument Company, Inc.
c/o RECYCLING
1595 Sycamore, Ave.
Bohemia, NY 11716

As a minimum the following substances, preparations and components have to be removed from any separately collected WEEE:

- Mercury containing components, such as switches or backlighting lamps (compact fluorescent lamps, CFL),
- Batteries
- Printed circuit boards if the surface of the printed circuit board is greater than 10 square centimeters (about 4 sq in.),
- Toner cartridges, liquid and pasty, as well as color toner,
- Chlorofluorocarbons (CFC), hydrochlorofluorocarbons (HCFC) or hydrofluorocarbons (HFC), hydrocarbons (HC)
- Liquid crystal displays (together with their casing where appropriate) of a surface greater than 100 square centimeters and all those back-lighted with gas discharge lamps,
- External electric cables
- Components containing refractory ceramic fibers as described in Commission Directive 97/69/EC of 5 December 1997 adapting to technical progress Council Directive 67/548/EEC relating to the classification, packaging and labeling of dangerous substances (2),
- Electrolyte capacitors containing substances of concern (height > 25 mm, diameter > 25 mm or proportionately similar volume)

2. The following components of WEEE that is separately collected have to be treated as indicated:

- Equipment containing gases that are ozone depleting or have a global warming potential (GWP) above 15, such as those contained in foams and refrigeration circuits: the gases must be properly extracted and properly treated. Ozone-depleting gases must be treated in accordance with Regulation (EC) No 2037/2000 of the European Parliament and of the Council of 29 June 2000 on substances that deplete the ozone layer (4).

Table of Contents

1	Introduction	6
1.1	<i>Koehler's Commitment to Our Customers</i>	6
1.2	<i>Recommended Resources and Publications</i>	6
1.3	<i>Instrument Specifications</i>	7
2	Safety Information and Warnings	7
3	Getting Started	8
3.1	<i>Packing List</i>	8
3.2	<i>Additional Accessories Required for Testing</i>	8
3.3	<i>Unpacking</i>	8
4	Setup	9
4.1	<i>Instrument Descriptions</i>	9
5	Operation	11
5.1	<i>Test Procedure</i>	11
6	Safety Features	12
6.1	<i>Over-Power Protection</i>	12
7	Maintenance	12
7.1	<i>Routine Maintenance</i>	12
7.2	<i>Instrument Cleaning</i>	12
7.3	<i>Replacement Parts</i>	12
8	Wiring Diagrams	13
8.1	<i>220-240V Unit Wiring</i>	13
9	Troubleshooting	14
9.1	<i>Unit does not power up</i>	14
9.2	<i>Unit is on and keeps resetting into start up routine</i>	14
10	Service	14
11	Storage	14
12	Warranty	14
13	Returned Goods Policy	15
	Notes	16

1 Introduction

The Koehler K13992 Cleveland Open Cup Flash Point Tester is the latest design for accurately determining flash and fire point temperatures of viscous petroleum products over an extended temperature range according to the ASTM D92 test method and related test specifications.

This manual provides important information regarding safety, technical reference, installation requirements, operating condition specifications, user facility resource requirements, and operating instructions for the Cleveland Open Cup Flash Point Tester. This manual should also be used in conjunction with applicable published laboratory procedures. Information on these procedures is given in section 1.2.

1.1 Koehler's Commitment to Our Customers

Providing quality testing instrumentation and technical support services for research and testing laboratories has been our specialty for more than 50 years. At Koehler, the primary focus of our business is providing you with the full support of your laboratory testing needs. Our products are backed by our staff of technically knowledgeable, trained specialists who are experienced in both petroleum products testing and instrument service to better understand your requirements and provide you with the best solutions. You can depend on Koehler for a full range of accurate and reliable instrumentation as well as support for your laboratory testing programs. Please do not hesitate to contact us at any time with your inquiries about equipment, tests, or technical support.

Toll Free: 1-800-878-9070 (US only)

Tel: +1 631 589 3800

Fax: +1 631 589 3815

Email: info@koehlerinstrument.com

<http://www.koehlerinstrument.com>

1.2 Recommended Resources and Publications

1. American Society for Testing and Materials (ASTM)
100 Barr Harbor Drive
West Conshohocken, Pennsylvania 19428-2959, USA
Tel: +1 610 832 9500
Fax: +1 610 832 9555
<http://www.astm.org>
email: service@astm.org

ASTM Publication:

- ASTM D92: Standard Test Method for Flash and Fire Points by Cleveland Open Cup Tester
- ASTM D6074: Standard Guide for Characterizing Hydrocarbon Lubricant Base Oils
- ASTM D6158: Standard Specification for Mineral Hydraulic Oils

2. International Organization for Standardization (ISO)
1, rue de Varembé
Case postale 56
CH-1211 Geneva 20, Switzerland
Tel: 41 22 749 01 11
Fax: 41 22 733 34 30
<http://www.iso.org>

ISO Publication:

- ISO 2592: Determination of Flash and Fire Points- Cleveland Open Cup Method

3. American Association of State Highway and Transport Officials (AASHTO)
444 North Capital Street N.W., Suite 249
Washington, D.C. 20001
Tel: (202) 624-5800
Fax: (202) 624-5806
E-mail: info@aaashto.org

AASHTO Publication:

- AASHTO T48: Flash and Fire Points by Cleveland Open Cup

4. Energy Institute (IP)
61 New Cavendish Street
London, W1M 8AR, United Kingdom
Tel: 44 (0)20 7467 7100
Fax: 44 (0)20 7255 1472
<http://www.energyinstpubs.org.uk/>

IP Publication:

- IP 36: Flash Point

5. Deutsche International Norm (DIN)
<http://www.din.de>

DIN Publication:

- DIN 51376: Flash and Fire Points by Cleveland Open Cup Tester

6. American National Standards Institute (ANSI)
1819 L Street, NW (between 18th and 19th Streets), 6th floor
Washington, DC 20036
Tel: 1.202.293.8020
Fax: 1.202.293.9287

ANSI Publication:

ANSI Z-11.6: Flash and Fire Points by Cleveland Open Cup Method

7. Federal Test Method (FTM)

FTM Publication:

- FTM 791-1103: Flashpoint by Cleveland Open Cup
- FTM 141-4294: Flash and Fire Points (Cleveland Open Cup)

1.3 Instrument Specifications

Models: K13992

Electrical Requirements: 220-240V 50/60Hz

Power: 1250W

Dimensions (wxdxh,in.(cm)): 10 x 17.5 x 11.5
(25.4 x 44.5 x 29.2)

Net Weight: 20 lbs (9.1kg)

Gross Weight: 30 lbs (13.6kg)

Altitude: Rated for use below 2000m

Environmental Conditions: As per section 1.4.1 of IEC 61010

2 Safety Information and Warnings

Safety Considerations. The use of this equipment may involve *hazardous* materials and operations. This manual does not purport to address all of the safety problems associated with the use of this equipment. It is the responsibility of any user of this equipment to investigate, research, and establish appropriate safety and health practices and determine the applicability of regulatory limitations prior to use.

Equipment Modifications and Replacement Parts. Any modification or alteration of this equipment from that of factory specifications is **NOT** recommended because it voids the manufacturer warranty, product safety, performance specifications, and/or certifications whether specified or implied, and may result in personal injury and/or property loss. Replacement parts must be O.E.M. exact replacement equipment.

Unit Design. This equipment is specifically designed for use in accordance with the applicable standard test methods listed in section 1.2 of this manual. The use of this equipment in accordance with any other test procedures, or for any other purpose, is not recommended and may be extremely hazardous.

Chemical Reagents Information. Chemicals and reagents used in performing the test may exhibit potential hazards. Any user must be familiarized with the possible dangers before use. We also recommend consulting the Material Data and Safety Sheet (MSDS) on each chemical reagent for additional information. MSDS information can be easily located on the internet at <http://siri.uvm.edu> or <http://www.sigma-aldrich.com>.

3 Getting Started

The instructions for preparing the equipment assume that the user is aware of the contents of this document, which lists the warranty conditions and important precautions.

3.1 Packing List


- K13992 Cleveland Open Cup Flash Point Tester, 220-240V 50/60Hz
- Cleveland Open Flash Cup
Precision machined brass Cup with heat resistant handle
- K13992-Manual K13992 Cleveland Open Cup Flash Point Tester Operation and Instruction Manual


3.2 Additional Accessories Required for Testing

- ASTM 11F Thermometer
Range: 20 to 760°F
- ASTM 11C Thermometer
Range: -6 to +400°C

3.3 Unpacking

1. Check Shock Watch Label on Cardboard Box for indication of rough handling and possible damage.
2. Check labeling for correct orientation of instrument. (e.g. This Side Up)
3. Carefully open top of box with box cutter and remove packing foam insert.
4. Extract instrument and place on suitable cart for transportation to work area / lab bench.

 **WARNING:** Be sure two or more individuals are available for extracting and lifting instrument from box to cart and from cart to bench. Individuals must lift in accordance to proper technique. See Figure below.



5. Lift instrument from cart and place on bench.
6. Carefully unpack and place the instrument and accessories in a secure location. Ensure that all parts listed on the packing list are present. Inspect the unit and all accessories for damage. If any damage is found, keep all packing materials and immediately report the damage to the carrier. We will assist you with your claim, if requested. When submitting a claim for shipping damage, request that the carrier inspect the shipping container and equipment. Do not return goods to Koehler without written authorization.


4 Setup

Equipment Placement. Make sure the instrument is placed on a firm, level table in an area with adequate ventilation or in a hood. The unit may be leveled by making minor turning adjustments to the feet located at the base of the unit. Please note that Koehler does not supply a level with this equipment.


Environmental Conditions: The instrument environment must comply with the following conditions for proper setup:

- No / Low Dust
- No direct sunlight
- Not near heating or AC ventilation ducts
- No Vibrations
- Clearance from other instruments
- Temperature Range: 5 to 40°C
- Elevation to 2000 meters
- Relative Humidity: < 80%


Ventilation. A fume hood or exhaust system is required for expending any fumes or vapors that have been generated while operating the unit. Flammable vapors and/or steam are generated during operation and must not be permitted to accumulate. A canopy-style hood may be used if the height from the top of the unit to the canopy is 5 feet or less. The exhaust blower should have a rating of 1000 C.F.M. or greater.

 **NOTE:** Adjust the exhaust blower rate so that it removes generated vapors / steam, however be sure that the exhaust system does not create a draft large enough to remove vapors from the test cup area or blow out the test flame as this will affect results of the test.

Gas Supply. Connect the gas inlet to a regulated low pressure gas supply (0.5 to 1.0 psi). Use only **propane, LPG, or natural gas**. Do not use direct unregulated pressure from an **LPG** tank.

 **NOTE:** Be sure Flame Control Switch (Figure 1, Item 3) is in the OFF position when Test Flame Applicator is not in use.

Power. Connect the line cords to properly fused and grounded receptacles with the correct voltage as indicated in section 1.3 or on the back of the unit.

 **WARNING:** For safety, disconnect the power when performing any maintenance and/or cleaning..

4.1 Instrument Descriptions





Figure 1: Instrument Descriptions_Front

1. Power Switch: Controls power to the entire unit. Pressing the switch to the ON position will energize the instrument. Pressing the switch to the OFF position will de-energize the instrument. Switch will illuminate when in the ON position

 **WARNING:** Be sure to completely Power Off instrument prior to performing any service of the instrument. This can only be done by switching the Power Switch to the **OFF** position. Turning the heating control dial to the off position **WILL NOT** de-energize the instrument. Only clean instrument or perform maintenance when power indicating light is off.

2. Temperature Control Dial. The analog dial is used to control the heating rate of the unit during the test procedure. This dial is non-linear, therefore, the numbers indicated on the dial plate **DO NOT** refer to specific temperatures or heating rates. The control dial

can be switched to an **OFF** position however, please **NOTE** that this **DOES NOT** power off the instrument.

3. **Flame Control ON/OFF Switch.** Controls cut off valve from gas line. Turn switch to OFF position when Flame Applicator is not in use to ensure gas is not flowing.
4. **Gas Adjustment Knob.** Turn to adjust the size of the test flame prior to testing. Adjustment knob designed for only minor adjustment of gas for proper sizing of test flame. Gas must be regulated prior to entry to instrument.
5. **Test Flame Applicator:** This is used to apply the test flame to the specimen. The flame should be applied across the top of the test cup to initiate a flash as per test method specifications.
6. **Thermometer Holder:** This holder is used to place the thermometer in the test cup to measure the temperature of the test specimen. The holder can be tilted back to a position away from the test cup when instrument is not performing a test.
7. **Cleveland Open Flash Cup:** The brass test cup contains the test specimen and is built with a heat resistant handle. An indicator line is engraved inside the test cup for proper sample level.
8. **Insulating Plate:** The insulating plate protects to operator from an extremely hot surface that will result from the heating of the test cup.
9. **Flame Guide:** Metal bead used for reference to determine the proper size of the test flame.



Figure 2: Instrument Descriptions_Back

10. **Gas Inlet:** For connection to source gas. Connect the gas inlet to a regulated low pressure gas supply (0.5 to 1.0 psi). Use only **propane, LPG, or natural gas**. Do not use direct unregulated pressure from an LPG tank.

5 Operation

5.1 Test Procedure

1. Be sure all tests be performed in a room free of excessive drafts. The room or compartment should be darkened sufficiently to allow the flash to be readily seen. Avoid breathing over the surface of the test oil.
2. Suspend the thermometer in a vertical position by the thermometer holder (**Figure 1, Item 5**) so that the bottom of the bulb is $\frac{1}{4}$ " from the bottom of the cup, and a point half way between the center and the back of the cup.
3. Fill the cup with the sample oil (approximately 70mL) so that the top of the meniscus is exactly at the filling line at ambient temperature. There should not be any test specimen on the outside of the test cup. If there is, empty, clean and refill the cup.



NOTE: Viscous samples should be heated until they are reasonably fluid before being poured into the cup, however the temperature during heating must not exceed 100°F (56°C) below the probable flash point.

4. Position the test cup on the center of the heating plate.
5. Press Flame Control Switch (**Figure 1, Item 3**) to the ON position.
6. Light Test Flame Applicator (**Figure 1, Item 5**) using external heat source.
7. Use a test flame approximately 5/23 of an inch in diameter, the same size as the test flame guide (**Figure 1, Item 8**). Turn the gas adjustment valve (**Figure 2, Item 10**) if needed, until the flame compares to the flame guide.
8. Press Power Switch (**Figure 1, Item 1**) to the ON position.
9. Turn Temperature Control Dial (**Figure 1, Item 2**) clockwise and apply heat initially at such a rate that the temperature on the thermometer increases 5 to 17°C (9 to 30°F) / min.
10. When the test specimen temperature is approximately 56°C (100°F) below the

expected flash point, decrease the rate so that the temperature rise during the last 28°C (50°F) before the flash point is 5 to 6°C (9 to 11°F) / min.

11. Apply the test flame when the temperature of the test specimen is approximately 28°C below the expected flash point and each time thereafter at a temperature reading that is a multiple of 2°C (5°F).
12. At the time of the next test flame application, pass the test flame in the opposite direction of the preceding application. The time consumed in passing the test flame across the test cup in each case shall be approximately 1 ± 0.1 s.
13. During the last 28°C (50°F) rise in temperature prior to the expected flash point, care shall be taken to avoid disturbing the vapors in the test cup with rapid movements or drafts near the test cup.
14. If foam persists on the top of the test specimen before the expected flash point, terminate the test and disregard any results.
15. Record, as the observed flash point, the reading on the temperature measuring device when a large flame appears and flash appears and instantaneously propagates itself over the entire surface of the test specimen.



NOTE: The application of the test flame can cause a blue halo or an enlarged flame prior to the actual flash point. This is not a flash point and shall be ignored.

16. To determine the Fire Point, after the flash point has been established and recorded, continue heating the oil at the specified rate of 5 to 6°C (9 to 11°F)/min and apply the test flame at the same intervals (every 2°C or 5°F) until the oil ignites and continues to burn for a period of at least 5 seconds. The temperature reading at the time of the flame application that caused burning for a period of 5 seconds or more is the fire point.

6 Safety Features

The Koehler K13992 Cleveland Open Cup Flash Point Tester is equipped with several safety and protection features, which are described in the following sections.

6.1 Over-Power Protection

The Koehler K13992 Cleveland Open Cup Flash Point Tester is equipped with Over-Power Protection circuitry, which prevents the unit from unsafe electrical conditions. If power to the unit is lost, then turn off the main power and turn it back on again. The main power switch also functions as a circuit breaker.



WARNING: Disconnect power to the unit before servicing and accessing any internal portion of the instrument to avoid exposure to high voltages and/or temperatures which may result in personal injury or death. If you have any questions about maintaining your equipment, then please do not hesitate to contact the Koehler technical service department.

7 Maintenance



WARNING: Disconnect power to the unit before servicing to avoid exposure to high voltages and/or temperatures which may result in personal injury or death. If you have any questions about maintaining your equipment, please do not hesitate to contact the Koehler technical service department.

7.1 Routine Maintenance

The K13992 Cleveland Open Cup Flash Point Tester requires little routine maintenance to provide many years of continuous service. However, over the course of time, some instrument parts may need to be replaced. When ordering replacement part(s), please provide the model number, serial number, and product shipment date of your equipment so that we can ensure you will receive the proper replacement part(s).

7.2 Instrument Cleaning

- To clean the instrument's exterior, which includes all painted surfaces, either a solution of soap and water or laboratory grade detergent may be used.

- Apply cleaner to clean wipe or cloth, not to the instrument directly. Wipe surface clean.
- **Do Not** clean bath exterior with organic chemicals such as Acetone, Toluene, Hexane, etc.
- For more difficult cleaning of finished surfaces, a dilute solution of isopropanol in water may be used.
- It is not recommended that more aggressive solvents be used on painted surfaces since paint color will tarnish or be stripped from the instrument.
- Stainless Steel surfaces, such as on the top plate, may be cleaned using a more aggressive solvent such as a stainless steel cleaner.




SHOCK AND BURN HAZARD: Only clean inside the bath when equipment is de-energized and unplugged from the mains power supply. Allow adequate time for heating coils to completely cool before cleaning.

7.3 Replacement Parts

Part Number	Description
225-230-002	Brick Heater, 1000W, 230V
K14000	Brass Test Cup
K13900-03120	Flame Application Assembly
036-108-01B	Needle Valve, Brass, 1/8 FPT
K13900-03117	Gas Inlet Tubing
271-010-014	Circuit Breaker, 2 Pole, 10A
K138-0-9	Burner Handle
K13900-03115	Thermometer Holder
K13900-03104	Heating Plate
K138-1-17	Insulation Plate
010-230-004	Wattstat, 230V
K160-1-2A	Flame Guide
K162-0-2	Heater Clips

8 Wiring Diagrams

8.1 220-240V Unit Wiring



KOEHLER INSTRUMENT CO. INC.	
1555 SYCAMORE AVE BOHEMIA NY 11716	
Title SCHEMATIC, POWERTRON HEATER 230V	
Size A	Document Number K13990-14000
Rev E	Date: Tuesday, April 30, 2013
Sheet 1 of 1	

REVISION HISTORY			
REV	ECO#	CHANGE DESCRIPTION	DATE
E	8188	CE GROUNDING PER ECO	4/30/13

9 Troubleshooting



WARNING: Troubleshooting procedures involve working with high voltages and/or temperatures which may result in personal injury or death, and should only be performed by trained personnel. Please do not hesitate to contact Koehler for assistance.

9.1 Unit does not power up

- Establish that the socket outlet is providing proper and adequate voltage.
- Check if Overpower Protection circuitry located directly behind the temperature controller inside the front tray has been activated.
- Check if line switch is in the **ON** position.
- Check fuse on wattstat.
- If problem persists, please call the Koehler technical service department for assistance.

9.2 Unit is on and keeps resetting into start up routine

- For 230V units, make sure that the socket outlet is greater than 215V.
- Check if there is a steady and reliable power source.
- Make sure the connector plug on the rear panel is firmly plugged in.

10 Service

Under normal operating conditions and with routine maintenance, the K13992 Cleveland Open Cup Flash Point Tester does not require service. Any service problem can be quickly resolved by contacting Koehler's technical service department either by letter, phone, fax, or email. In order to assure the fastest possible service, please provide us with the following information.

Model Number: _____

Serial Number: _____

Date of Shipment: _____

11 Storage

This laboratory test instrument is equipped with electrical components. Storage facilities should be consistent with an indoor laboratory environment. This testing equipment should not be subjected to extremes of temperature and/or moisture. This equipment was shipped from the factory in a corrugated cardboard container. If long term storage is anticipated, re-packing the instrument in a water-resistant container is recommended to ensure equipment safety and longevity.

12 Warranty

We, at Koehler, would like to thank you for your equipment purchase, which is protected by the following warranty. If within one (1) year from the date of receipt, but no longer than fifteen (15) months from the date of shipment, Koehler equipment fails to perform properly because of defects in materials or workmanship, Koehler Instrument Company, Inc. will repair or, at its sole discretion, replace the equipment without charge F.O.B. its plant, provided the equipment has been properly installed, operated, and maintained. Koehler Instrument Company must be advised in writing of the malfunction and authorize the return of the product to the factory. The sole responsibility of Koehler Instrument Company and the purchaser's exclusive remedy for any claim arising out of the purchase of any product is the repair or replacement of the product. In no event shall the cost of the purchaser's remedy exceed the purchase price, nor shall Koehler Instrument Company be liable for any special, indirect, incidental, consequential, or exemplary damages. KOEHLER INSTRUMENT COMPANY, INC. DISCLAIMS ALL OTHER WARRANTIES, EXPRESSED OR IMPLIED, INCLUDING ANY IMPLIED WARRANTIES OF FITNESS FOR A PARTICULAR PURPOSE. Please save the shipping carton in the event the equipment needs to be returned to the factory for warranty repair. If the carton is discarded, it will be the purchaser's responsibility to provide an appropriate shipping carton.

13 Returned Goods Policy

To return products for credit or replacement, please contact Koehler Customer Service with your purchase order number, our packing list/invoice number, the item(s) to be returned and the reason for the return. You will be issued a Returned Authorization (RA) number, which must be prominently displayed on the shipping container when you return the material to our plant. Shipping containers without an RA number prominently displayed will be returned to the sender. Goods must be returned freight prepaid. Returns will be subject to a restocking charge, the application of which will depend upon the circumstances necessitating the return. Some returns cannot be authorized, including certain products purchased from outside vendors for the convenience of the customer, products manufactured on special order, products shipped from the factory past ninety (90) days, and products which have been used or modified in such a way that they cannot be returned to stock for future sale.

