

Vapor Pressure of Petroleum Products



Automatic Microscale Vapor Pressure Analyzer

The design of vapor pressure determination in the instrument is based on Triple-expansion principle. Inject the sample of known volume into a temperature controlled measuring room with a piston in it, seal the measuring room, expand the volume of the sample by X times in three steps, measure the total pressure after every expansion step and calculate the dissolved partial pressure of the air and dissolved pressure of air in the sample. And then raise the temperature to a certain value and measure the total pressure at the temperature.

Dimensions wxdxh,in.(cm)

8.9x11.5x15.4 (22.5x29x39)

Net Weight: 26.5 lb (12kg)

Test Method

To determine the total vapor pressure exerted in vacuum by air-containing, volatile, liquid petroleum products, including automotive spark-ignition fuels with or without oxygenates. This test method is suitable for testing samples with boiling points above 0 °C (32 °F) that exert a vapor pressure between 7 kPa and 130 kPa (1.0 psi and 18.6 psi) at 37.8 °C (100 °F) at a vapor-to-liquid ratio of 4:1. Measurements are made on liquid sample sizes in the range from 1 mL to 10 mL. No account is made for dissolved water in the sample.

To determine the vapor pressure exerted in vacuum by volatile, liquid petroleum products, hydrocarbons, and hydrocarbon-oxygenate mixtures. This test method is suitable for testing samples with boiling points above 0°C (32°F) that exert a vapor pressure between 7 and 150 kPa (1.0 and 21 psi) at 37.8°C (100°F) at a vapor-to-liquid ratio of 4:1. The liquid sample volume size required for analysis is dependent upon the vapor-to-liquid ratio chosen and the measuring chamber volume capacity of the instrument.

Specifications

Conforms to the specifications of:

ASTM D5191, D6378;

Excellent Correlation to:

ASTM D323, D4953; EN 13016-1; IP394, JIS K2258-2; SH/T 0794; SH/T 0769; GB/T 8017; SN/T 2932

Temperature Range: 0 to 120°C (30 to 250°F)

Temperature Stability +/- 0.1°C

Pressure Range: 0 to 1000 kPa (0 to 145 psi)

Pressure Resolution: 0.01 kPa

Automated Calibration Routine

Fast & Easy On-site Calibration – No PC required

Sample Volume: 1 mL (2.5 mL rinse - 3x per cycle)

Vapor/Liquid Ratio: 0.02/1 to 4/1

Interfaces:

USB, RS 232, Ethernet

Power Supply:

100-240VAC,50/60Hz, 150W

Ordering Information

Catalog No.

K24870 Automatic Microscale Vapor Pressure Analyzer 100-240V 50/60Hz



1595 SYCAMORE AVENUE • BOHEMIA, NEW YORK 11716-1796
1-800-878-9070 (IN U.S. ONLY) • TEL: +1 631 589 3800 • FAX: +1 631 589 3815
www.koehlerinstrument.com • Email: sales@koehlerinstrument.com

©2019 Koehler Instrument Company, Inc.