

FREND™ System

COVID-19 IgG/IgM Duo

Microfluidic Qualitative Immunoassay

The FREND™ COVID-19 IgG/IgM Duo is a **point-of-care testing (POCT)** which can be used to check whether patients have developed immune responses to SARS-CoV-2 using human serum or plasma.

For COVID-19 IgG/IgM Duo, it provides results **qualitatively** but the principle of detection is **quantitative**.

The FREND™ System shows the test result as negative if the Cut-Off-Index (COI) is < 1.0 . If the Cut-Off-Index (COI) is ≥ 1.0 , then the test result will be provided as positive. The Cut-Off-Index (COI) is determined quantitatively by testing the specimens collected after day 8 from the onset of symptoms using the FREND™ System.

Fast




Easy



Monitoring




Steps




STEP 01.

Dilute sample.



STEP 02.


Add **35 µL** of the diluted sample.



STEP 03.

Insert the cartridge and get results.

Laboratory Information System Connectivity



In the event of a pandemic, the management of a vast amount of clinical results is important. As many laboratories face challenges in effective arrangement of the essential information⁽¹⁾ the **FREND™ System** provides the following advantages:

- **Paperless**
- **Time-saving**
- **Automatic Data Transfer & Collection**
- **Convenient Data Management**
- **Can help monitoring (serosurveillance)**

Clinical Performance

The total of 59 clinical specimens (43 positive and 16 negative) confirmed with RT-PCR were tested with the FRENDS™ COVID-19 IgG/IgM Duo.

The total result of all samples

Method			Clinically Confirmed Specimen	
			Positive	Negative
FRENDS™ COVID-19 IgG/IgM Duo	Positive	IgG(+), IgM(+)	14	0
		IgG(+), IgM(-)	26	0
		IgG(-), IgM(+)	0	0
	Negative		3	16
Total			43	16



*Positive Percent Agreement (PPA), Negative Percent Agreement (NPA)

The following table includes 16 negative specimens from the previous table and new 325 negative specimens, making a total of 341 negative specimens.

The total result of all samples

Method		Clinically Confirmed Specimen	
		Positive	Negative
FRENDS™ COVID-19 IgG/IgM Duo	Positive	40	10
	Negative	3	331
Total		43	341



The test result collected specimens after 8 days from symptom onset

Method		Clinically Confirmed Specimen	
		Positive	Negative
FRENDS™ COVID-19 IgG/IgM Duo	Positive	40	10
	Negative	0	331
Total		40	341



* The table shows result of specimens collected after 8 days from symptom onset. Three specimens were excluded in the original data.

Reference:

- (1) Division of Laboratory Systems (DLS), Clinical Laboratory COVID-19 Response Call
- (2) CDC Tests for COVID-19 from Centers for Disease Control and Prevention.
- (3) Profiling Early Humoral Response to Diagnose Novel Coronavirus Disease (COVID-19), Li Guo.
- (4) Purpose and Options for Testing for SARS-CoV2 (the COVID-19 Virus): Considerations for World Bank Task Teams Managing COVID-19 Fast Track Facility Operations, The World Bank
- (5) Cellular immune responses to severe acute respiratory syndrome coronavirus infection in senescent BALB/c Mice: CD4+ T cells are important in control of SARS-CoV infection. Jun Chen.
- (6) Chronological evolution of IgM, IgA, IgG and neutralisation antibodies after infection with SARS-associated coronavirus, P.-R. Hsueh.
- (7) SARS-CoV-2 Viral Load in Upper Respiratory Specimens of Infected Patients, Lirong Zou.
- (8) Clinical course and risk factors for mortality of adult inpatients with COVID-19 in Wuhan, China: a retrospective cohort study, Fei Zhou.

Why serological test is necessary?

The serological test using blood samples is required to detect whether symptoms developed from infection or the infection was asymptomatic⁽²⁾. After the onset of symptoms, the IgM starts to increase and then the IgG production begins from the next 7 days. Therefore, it is recommended to use the FRENDS™ System at a minimum day 8 (Fig. 1)⁽³⁻⁷⁾. After the onset of symptoms, the sensitivity of positive PCR decreased below 80% and 50% on day 6 and day 14, respectively. The combination of PCR with serological test increased the positive sensitivity by 47% (Fig. 2)⁽⁸⁾.




Fig. 1




Fig. 2

Specification

Item	Specification
Assay method	Fluorescence immunoassay
Specimen	Serum or Plasma
Sample volume	35 µL
Reaction time	< 3 min
Quantity	20 tests/kit
Storage condition	2 ~ 8 °C (35 ~ 46 °F)



Ordering Information

Product	Cat. No.
FRENDS™ System	F10
FRENDS™ COVID-19 IgG/IgM Duo	FRCOD 020

Related Products

Cat. No.	Product	Quantity (tests/kit)
EPCO19 100	COVID-19 RT-PCR Kit	100
FRCOG 020P	COVID-19 Ag	20
FRCGFL 020	COVID-19 Ag+FLU A&B	20
FRCOA 020	COVID-19 total Ab**	20
FRCOS 020	COVID-19 SP	20

**FDA EUA Approved



Scan for more information.

NanoEntek, Inc.

Head Office

12F, 5, Digital-ro 26-gil, Guro-gu, Seoul, 08389, Korea
Tel: +82-2-6220-7940/Fax: +82-2-6220-7999

NanoEntek America, Inc.

Waltham, MA, USA Tel: +1-781-472-2558

NanoEntek Bio-Technology(Beijing) Ltd.

Fengtai District, Beijing Tel: +86-10-5920-7980

website

www.nanoentek.com

e-mail

ivdst@nanoentek.com