

# FREND™ COVID-19 SP

Qualitative and semi-quantitative assay for COVID-19  
IgG and IgM of Spike Glycoprotein (S1)

---

**REF** FRCOS 020

**IVD** For *in vitro* diagnostic use only

## Intended use

The FREND™ COVID-19 SP is a microfluidic fluorescence immunoassay (FIA) using the FREND™ System intended for the qualitative and semi-quantitative detection of IgG and IgM antibodies to SARS-CoV-2 in human serum and plasma (EDTA). The FREND™ COVID-19 SP is intended for use as an aid in identifying individuals with an adaptive immune response to SARS-CoV-2, indicating recent or prior infection. At this time, it is unknown for how long antibodies persist following infection and if the presence of antibodies confers protective immunity.

The FREND™ cartridge utilizes micro-fluidics lateral flow technology where the analyte of interest in the sample forms immune complexes while moving through the fluidics pathway in the cartridge. Spike glycoprotein (SP) which is used to detect SP-antibodies has receptor-binding domain (RBD) fragment. Testing is limited to laboratories certified under the Clinical Laboratory Improvement Amendments of 1988 (CLIA), 42 U.S.C §263a, to perform moderate or high complexity tests and as applicable, point-of-care (POC) testing.

Results are for the detection of SARS-CoV-2 antibodies. IgG and IgM antibodies to SARS-CoV-2 are generally detectable in blood several days after initial infection, although the duration of time antibodies are present post-infection is not well characterized. Individuals may have detectable virus present for several weeks following seroconversion.

The sensitivity of FRENDS™ COVID-19 SP early after infection is unknown. Negative results do not preclude acute SARS-CoV-2 infection. If acute infection is suspected, direct testing for SARS-CoV-2 is necessary.

False positive results for FRENDS™ COVID-19 SP may occur due to cross-reactivity from pre-existing antibodies or other possible causes. Due to the risk of false positive results, confirmation of positive results should be considered using second, different IgG and IgM assay.

## **Summary and explanation of test**

COVID-19 is a respiratory disease caused by a novel coronavirus called SARS-CoV-2. The symptoms for patients have included mild to severe respiratory illness with fever, cough, and difficulty breathing.<sup>1-2</sup> The test result shows the presence of SARS-CoV-2 antibodies, IgG or IgM, are generally detectable in the blood several days after symptom onset.<sup>3-4</sup> The antibody test provides information specific to an adaptive immune response to SARS-CoV-2, indicating recent or prior infection. Presence of antibody may not be an indicator of immunity. Reactive or positive results should be confirmed by a second serology test. Final recommendations and patient management should be determined by healthcare professional based on clinical symptoms and other diagnostic tests.

## **Principle of the assay**

The FRENDS™ cartridge utilizes micro-fluidics lateral flow technology where the analyte of interest in the sample forms immune complexes while moving through the fluidics pathway in the cartridge.

A specimen is added to sample dilution tube and mixed. A well-mixed sample of 35 µL is transferred to the sample inlet of a single use FRENDS™ COVID-19 SP cartridge. The cartridge is then placed into the FRENDS™ System, which is programmed to begin analysis once the sample has reacted with the reagents. The reaction and analysis time is approximately 3-4 minutes. The anti-coronavirus IgG and/or IgM qualitative and semi-quantitative measurement is based on the ratio of fluorescence detected by the FRENDS™ System at the FRENDS™ Test and Reference zones. The magnitude of the fluorescent ratio is proportional to the presence and absence of IgG and IgM in the sample.

The FREND™ System is a bench-top fluorescence reader containing a touchscreen user interface. The System has a slot that accepts the FREND™ COVID-19 SP test cartridge (which contains the reagents and sample) and is programmed to analyze the test when the sample has fully reacted with the on-board cartridge reagents. Results of the test are displayed on the screen and can be printed on an optional printer.

The FREND™ System is a lateral flow immunofluorescence analyzer that can provide an semi-quantitative value through a function obtained from the ratios of test and reference values. The concentration of immunoglobulins against the spike glycoprotein of SARS-CoV-2 is determined quantitatively using FREND™ System. The FREND™ System shows the test result as non-reactive or reactive depending on the predetermined cutoff value saved in the calibration Code chip.

### Material provided

Q'ty	Contents	Catalogue number
20	Cartridges	FRCOS 020
20	Dilution tubes	
40	Disposable pipette tips	
01	Code chip	
01	Package insert	

### Components required but not included with the test

The following materials are not provided with the reagent but are required to perform COVID-19 IgG/IgM analysis using the FREND™ COVID-19 SP on the FREND™ System.

- FREND™ System including calibrated pipette,
- QC Cartridge and QC Code chip manufactured by NanoEntek.

### Warning and Precautions

- The FREND™ COVID-19 SP cartridges are intended for *in vitro* diagnostic use only.
- The FREND™ COVID-19 SP cartridges are only to be used on the NanoEntek FREND™ System.

- The FREND™ COVID-19 SP cartridges and dilution tubes are disposable, single use devices. Do not reuse!
- Allow sealed cartridges to come to room temperature for 15-30 minutes prior to use.
- Cartridges and dilution tubes should not be frozen.
- The humidity in the laboratory must be 10-80% range when running tests.
- Avoid cross-contamination between samples by using a new pipette tip for each new specimen when transferring the sample to the dilution tube. Use another new pipette tip when transferring are diluted, sample to the cartridge.
- Avoid high humidity, direct sunlight or heat when storing cartridges and dilution tubes.
- Testing of contaminated samples may cause erroneous results.
- Using specimens containing clotted fibrin could result in erroneous results.
- Over or under loading the cartridge with sample may result in inaccurate results.
- Do not use the cartridges beyond the expiration date on the pouch.
- Do not use the cartridge if the pouch is damaged or the seal is broken.
- Perform testing as specified in the Package insert and User manual
- Keep the cartridge sealed in the pouch until ready for use.
- Use the cartridge immediately after opening the pouch.
- Handle specimens in accordance with the OSHA Standard on Bloodborne Pathogens.
- Human specimens are not used in the preparation of this product, however, since human specimens will be used for samples and other quality control products in the laboratory may be derived from human materials. Use Universal Precautions when handling all specimens and controls. Wear disposable gloves when handling the cartridges and the samples.
- Wash hands thoroughly and often after handling reagent cartridges or samples.
- Do not ingest the silica gel packet found in the cartridge pouch.
- Used cartridges and diluent tubes should be disposed in accordance with local or national regulations.

## **Storage and Stability**

All unopened materials are stable until the expiration date on the label when stored at refrigerator temperature storage (2-8°C). Reagent stability has been demonstrated for twelve months from the date of manufacture. The expiration date is clearly indicated on the product box and the cartridges.

## **Specimen collection and handling**

Serum or plasma (EDTA) is required for the assay. No special patient preparation is necessary. Collect the appropriate venous blood sample in accordance with standard laboratory procedures. For serum, allow the sample to clot for 30 minutes at room temperature. For plasma (EDTA), centrifuge after collection. Centrifuge the sample for 10 minutes at 3,000 rpm within 2 hours of collection and immediately separate the serum or plasma from the packed cells.

It is recommended to use samples immediately. However, if testing is not done immediately, samples may be stored at 2-8°C for up to 7 days prior to testing. If testing is not performed within 7 days, store at -20°C or below. Samples can be stored frozen for up to 3 months prior to testing.

Repeated freeze-thaw cycles should be avoided. Turbid samples containing particulate matter such as fibrin clots or visible strands should be centrifuged before being tested. Prior to assay, slowly bring frozen samples to room temperature and mix gently but thoroughly before testing.

For optimal results, avoid grossly hemolytic, lipemic, or turbid specimens. Specimens should be free of aggregated fibrin, red blood cells, or other particulate matter. When pipetting into the cartridge sample inlet, ensure that bubbles in the sample are avoided. Bubbles may restrict flow and result in an incomplete or erroneous test result.

## Procedure

- Reagent preparation

There is no reagent preparation required to run the FREND™ COVID-19 SP cartridge on the FREND™ System. However, the cartridge and dilution tube should be incubated at room temperature for 15-30 minutes prior to testing.

- Code chip installation


A lot-specific Code chip is supplied with each kit of FREND™ COVID-19 SP. When using a new lot of reagent, the Code chip of the same lot must be installed in the FREND™ System. Please refer to the FREND™ System User manual for more detailed instructions relative to the Code chip installation. Abbreviated instructions are as follows:

- (1) Insert the FREND™ System electrical cord into an appropriate outlet.
- (2) Insert the Code chip into the Code chip slot at the rear of the system following the arrows.
- (3) Press the **'Setup'** button on the **'Main'** screen.
- (4) Press the **'Code chip'** button on the **'Setup'** screen.
- (5) The information embedded on the FREND™ COVID-19 SP Code chip is automatically saved on the FREND™ System.
- (6) When the Code chip installation is completed, press the **'OK'** button to go to the **'Setup'** screen.
- (7) Press the **'Item'** button on the **'Setup'** screen.
- (8) Check the FREND™ COVID-19 SP cartridge lot number and the installation date of the Code chip.
- (9) Press the **'Home'** button to go to the **'Main'** screen to begin running external quality control and patient samples.

- Specimen processing

Allow the tubes and the sealed pouches containing the FREND™ COVID-19 SP cartridges and dilution tubes to come to room temperature for 15-30 minutes prior to use. If using refrigerated patient samples, remove those from the refrigerator and allow to them to come to room temperature prior to testing. If frozen samples will be utilized, be sure these are removed from the freezer, thawed naturally and then mixed gently but thoroughly prior to testing. Testing should not begin on frozen samples until they have reached room temperature.

- Assay procedure



- (1) Have the FREND™ COVID-19 SP cartridge and samples ready for testing.
- (2) Record the Sample ID on the cartridge in the designated area.
- (3) Using the micropipette, transfer 35 µL of the sample to a sample dilution tube and mix by inverting the sample gently for 3-5 times. Using the mixed sample and a new pipette tip, transfer 35 µL into the sample inlet.
- (4) Press the 'Test' button on the 'Main' screen of the FREND™ System.
- (5) The system moves to the Patient ID screen automatically.
- (6) Type the Patient ID and press the 'Enter' button to begin the test.
- (7) Insert the cartridge into the cartridge slot using the cartridge arrows as a guide.

⚠ **Caution:** Please check the direction of the cartridge before insertion and assure the insertion is complete.

⚠ **Caution:** Please insert the cartridge into the FREND™ System after loading the sample into the cartridge.

- (8) When the reaction in the cartridges is complete, the FREND™ System will automatically begin the reading.
- (9) When the reading has been completed, the cartridge will automatically be expelled and the results displayed.
- ⚠ **Caution:** Do not remove power from the FREND™ System while a cartridge is in the reading chamber. This may cause a system error.
- (10) If the FREND™ System is connected to the optional printer, press the 'Print' button and the results will be output on the printer paper.
- ⚠ **Note:** For more detailed instructions, please refer to the 'FREND™ System User manual'.

## Control materials

- FREND™ System check

It is recommended for operators to use the QC Cartridge daily for the maintenance of the FREND™ System. Install QC Code chip in the FREND™ System before using the QC Cartridge. The QC Cartridge confirms the proper function of the FREND™ System including:

- Step 1 Laser power
- Step 2 Laser alignment
- Step 3 Calculate ratio

Please use the QC Cartridge and QC Code chip provided with the FREND™ System.

- Internal control

FREND™ COVID-19 SP cartridge contains built in control features. Fluorescence signal in the Reference zone of each cartridge shows: (1) that enough sample volume is added, (2) that proper flow is obtained, and (3) that the antibody is reactive. If this Reference zone signal is missing or lower than the threshold, FREND™ System considers it as an incorrect or failed test, and produces an error message instead of a test result. In addition, with each cartridge run, the system monitors, in part, for (1) flow of sample, (2) speed of sample flow, (3) shelf-life of cartridge components, (4) function of internal barcode scanner, and (5) function of scanner's mechanical components.



- Control materials

Control materials are available for purchase from NanoEntek (product name: COVID-19 SP Control LQ, cat.# FIC-COSLQ) and are not supplied with the kit. Commercially available controls may be used. Good laboratory practice suggests that controls are run routinely to ensure that test reagents are working and that the test is correctly performed. External controls should be used in accordance with local, state, federal accrediting organizations, or your laboratory's standard quality control procedures, as applicable.

## **Interpretation of results**

### Reactive result

"Reactive" indicates that IgG and/or IgM antibodies to COVID-19 SP were detected.

### Non-reactive result

"Non-reactive" indicates that IgG and/or IgM antibodies to COVID-19 SP were not detected.

### FREND™ System results display

FREND™ System qualitatively and semi-quantitatively detects of IgG and IgM antibodies to the SARS-CoV-2 virus.

Report options are indicated below: Note the following precautions:

- Non-reactive or negative results do not rule out SARS-CoV-2 infection, particularly in those who have been in contact with the virus. Follow-up testing with a molecular assay should be considered to rule out infection in these individuals.
- Results from antibody testing are not an indication of SARS-CoV-2 infection. Presence of SARS-CoV-2 antibodies are an indication of potential immunity to the virus.
- Reactive or positive results may indicate past or present infection with non-SARS-CoV-2 coronavirus strains, such as coronavirus HKU1, NL63, OC43, or 229E. The reactive or positive results should be confirmed using a second serology test.
- Not for the screening of donated blood.

Please read the section 'Interpretation of results' and 'Limitation of the procedure' carefully before you use the FREND™ COVID-19 SP and the FREND™ System.

## Display of Results

Displayed result	Description
 <p> <b>Displayed result</b>            Result            Print    Send to LIS            DateTime: 2021-04-19 10:21:43            Patient ID: NP            User ID: LEE-01            Order#: 765            Lab ID: NANCENTER            COVID SP IgG : Non-reactive (&lt;1.50 BAU/mL)            COVID SP IgM : Non-reactive (&lt;0.10 BAU/mL)         </p>	<p> <b>Description</b>            COVID SP IgG : Non-reactive (&lt;1.50 BAU/mL)            COVID SP IgM : Non-reactive (&lt;0.10 BAU/mL)         </p>
 <p> <b>Displayed result</b>            Result            Print    Send to LIS            DateTime: 2021-04-19 10:21:43            Patient ID: NP            User ID: LEE-01            Order#: 765            Lab ID: NANCENTER            COVID SP IgG : Reactive (&gt;100.00 BAU/mL)            COVID SP IgM : Reactive (&gt;20.00 BAU/mL)         </p>	<p> <b>Description</b>            COVID SP IgG : Reactive (&gt;100.00 BAU/mL)            COVID SP IgM : Reactive (&gt;20.00 BAU/mL)         </p>
 <p> <b>Displayed result</b>            Result            Print    Send to LIS            DateTime: 2021-04-19 10:21:43            Patient ID: NP            User ID: LEE-01            Order#: 765            Lab ID: NANCENTER            COVID SP IgG : Reactive (30.00 BAU/mL)            COVID SP IgM : Reactive (10.00 BAU/mL)         </p>	<p> <b>Description</b>            COVID SP IgG : Reactive (30.00 BAU/mL)            COVID SP IgM : Reactive (10.00 BAU/mL)         </p>
 <p> <b>Displayed result</b>            Result            Print    Send to LIS            DateTime: 2021-04-19 10:21:43            Patient ID: NP            User ID: LEE-01            Order#: 765            Lab ID: NANCENTER            COVID SP IgG : Non-reactive (2.50 BAU/mL)            COVID SP IgM : Reactive (0.20 BAU/mL)         </p>	<p> <b>Description</b>            COVID SP IgG : Non-reactive (2.50 BAU/mL)            COVID SP IgM : Reactive (0.20 BAU/mL)         </p>
 <p> <b>Displayed result</b>            Result            Print    Send to LIS            DateTime: 2021-04-19 10:21:43            Patient ID: NP            User ID: LEE-01            Order#: 765            Lab ID: NANCENTER            COVID SP IgG : Reactive (3.50 BAU/mL)            COVID SP IgM : Non-reactive (0.13 BAU/mL)         </p>	<p> <b>Description</b>            COVID SP IgG : Reactive (3.50 BAU/mL)            COVID SP IgM : Non-reactive (0.13 BAU/mL)         </p>

## **Limitation of the procedure**

- This test is for clinical laboratory use only. Not for home use.
- This test cannot rule out diseases caused by other bacterial or viral pathogens.
- This test must be read using only the FREND™ System.
- Immunocompromised patients who have COVID-19 may have a delayed antibody response and produce levels of antibody which may not be detected as Reactive (Positive) by the assay.
- The assay procedure and results interpretation must be followed closely when testing for the presence of SARS-CoV-2 virus specific antibodies. For optimal test performance, proper sample collection is critical. Failure to follow the procedure may give inaccurate results.
- If symptoms persist and the result of the COVID-19 SP test is non-reactive or negative, it is recommended to re-sample the patient a few days later or test with an alternative test device.
- The results obtained with this test should only be interpreted in conjunction with clinical finding, and the results from other laboratory tests and evaluations.
- Heterophilic antibodies in serum specimens may cause interference in immunoassay. These antibodies may be present in blood samples from individuals regularly exposed to animals or who have been treated with animal serum products. Results that are inconsistent with clinical observations indicate the need for additional testing.
- This test should not be used for screening of donated blood.
- Non-reactive or negative result should not rule out COVID-19 infection. Molecular testing should be consider along with other clinical examinations to confirm or eliminate infection.
- In the early phase of infection, the test sensitivity rate can be lower than the sensitivity rates showed in the 'Performance evaluation' section in the Package insert.

## Performance evaluation

### Cross-reactivity

No cross reactivity with 35 cross-reactive specimen were observed in the testing of the FRENDS™ COVID-19 SP.

Category	Result	
	IgG	IgM
N/A	-	-
HIV positive plasma	-	-
Japanese Encephalitis positive plasma	-	-
Zika virus positive plasma	-	-
Chikungunya positive plasma	-	-
Dengue IgM positive plasma	-	-
Salmonella typhi IgM positive plasma	-	-
Rubella IgM positive plasma	-	-
CMV IgG/IgM positive plasma	-	-
Tick borne encephalitis IgM positive plasma	-	-
West Nile Virus positive plasma	-	-
Treponema palladium positive plasma	-	-
HAV IgM positive plasma	-	-
HAV IgG positive plasma	-	-
HBV Ab positive plasma	-	-
HCV Ab positive plasma	-	-
Leishmania positive plasma	-	-
Brucella IgM positive plasma	-	-
Chagas positive plasma	-	-
Toxoplasma positive plasma	-	-
RSV positive serum	-	-
Mycoplasma pneumonia IgM positive plasma	-	-
Mycoplasma pneumonia IgG positive plasma	-	-
Influenza A IgM positive plasma	-	-

Influenza B IgM positive plasma	-	-
Influenza A and B IgG, IgM positive plasma	-	-
Parainfluenza 1, 2, 3 positive plasma	-	-
Adenovirus positive plasma	-	-
Enterovirus positive plasma	-	-
Human Rhinovirus positive plasma	-	-
SARS-CoV positive plasma	-	-
MERS-CoV positive plasma	-	-
HCoV OC43 positive plasma	-	-
HCoV 229E positive plasma	-	-
HCoV HKU1 positive plasma	-	-
HCoV NL63 positive plasma	-	-

### **Interference**

The interference evaluation test of FREND™ COVID-19 SP was conducted according to CLSI guidelines EP7-A2 using one lot. No interference in the testing of the FREND™ COVID-19 SP with 12 interfering substances was observed.

<b>Substances</b>	<b>Concentration tested</b>
Hemoglobin	2g/L
Bilirubin (unconjugated)	20 mg/dL
Bilirubin (conjugated)	20 mg/dL
Glucose	1000 mg/dL
Ascorbic acid	3 mg/dL
Human serum albumin	15 g/L
Triglycerides (total)	500 mg/dL
Rheumatoid factor	100 IU/mL
Cholesterol (Total)	250 mg/dL
Ibuprofen	2425 µmol/L
Zanamivir	3.3 mg/mL
Tamiflu	25 mg/mL

### **Clinical performance**

Eighty eight (88) clinical specimens (serum and plasma) were collected at Kangwon National University Hospital, Department of diagnostic examination, 156 Baeknyeong-ro, Chuncheon-si, Gangwon 24289, Korea from April 7<sup>th</sup> to May 28<sup>th</sup> 2020. Patients were confirmed for SARS-CoV-2 by rtPCR with the Allplex™ 2019-nCoV Assay (Seegene, Inc.). The results were as follows.















Method			rtPCR		Total
			Positive	Negative	
FREND™ COVID-19 SP	Reactive (Positive)	IgG(+), IgM(+)	32	0	32
		IgG(+), IgM(-)	1	0	1
		IgG(-), IgM(+)	0	2	2
	Non-reactive (Negative)	1	52	53	
Total			34	54	88

- Positive percent agreement: 97.06% (33/34)
- Negative percent agreement: 96.30% (52/54)

## References

1. Coronavirus Disease 2019 (COVID-19) Frequently Asked Questions, FDA
2. Q&A Detail, Coronavirus disease (COVID-19), WHO
3. Yu Chen, Lanjuan Li, SARS-CoV-2: virus dynamics and host response, *Lancet Infectious Diseases*. Vol. 20, May 2020
4. Kelvin Kai-Wang To MD, et al., Temporal profiles of viral load in posterior oropharyngeal saliva samples and serum antibody responses during infection by SARS-CoV-2: an observational cohort study, *Lancet Infectious Diseases*, Vol. 20, May 2020

## Glossary of symbols

	Caution, warning, Consult accompanying documents
	Catalogue number/Reference number
	Lot number/Batch number
	Use by YYYY-MM-DD or YYYY-MM
	Manufacturer
	Authorized representative in the European Community
	CE marking
	<i>In vitro</i> diagnostic medical device
	Temperature limitation
	Contains sufficient for <n> tests
	Do not reuse
	Do not use if package is damaged
	For prescription use only
	Irritant






**NanoEntek**

ivdst@nanoentek.com

www.nanoentek.com

 **Manufactured by  
NanoEntek, Inc.**

851-14 Seohae-ro, Paltan-myeon, Hwaseong-si, Gyeonggi-do, 18531, Korea  
Tel.:+82-2-6220-7942, Fax.:+82-2-6220-7999

 **MT Promedt Consulting GmbH**

Altenhofstrasse 80, 66386 St. Ingbert, Germany