# Influenza A & B Rapid Card

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※ Product code : 11311 ※ Please read this Instructions before use.

## INTENDED PURPOSE

BioTracer™ Influenza A&B Rapid Card is intended for the qualitative determination of influenza A virus and influenza B virus from nasal / throat swab specimens of symptomatic patients.

## EXPLANATION OF THE TEST

Human influenza A and B viruses are major infectious pathogens of the respiratory tract that affect people of all ages, with the highest rates of mortality and severe morbidity in the elderly and the young. Patients who present with suspected influenza may benefit from treatment with an antiviral agent especially if given within the first 48 hrs of onset of illness. It is important to distinguish influenza A from influenza B in order to allow physicians a choice in selective antiviral intervention. Diagnostic tests available for influenza include rapid immunoassay, immunofluorescence assay, polymerase chain reaction(PCR), serology, and viral culture.

BioTracer™ Influenza A&B Rapid Card is a qualitative rapid test to detect influenza A or B antigens from respiratory specimens of symptomatic patients

## PRINCIPLE OF THE METHOD

BioTracer™ Influenza A&B Rapid Card is based on the principle of an immunochromatography in vitro test for the qualitative determination of influenza A or influenza B. When the test device including specimens dissolved in extraction solution, the specimens moves to the conjugate pad via sample pad and resuspends anti-influenza A or anti-influenza B conjugated gold complex that is dried on the conjugate pad. The mixture moves along the membrane by capillary action and reacts with anti-influenza A or anti-influenza B that is immobilized on membrane. If influenza A or influenza B is present enough in the sample, a colored band is appeared. If there is no influenza A or influenza B in the sample or not sufficient, the area will remain colorless. As the control band(C) which is red or purple color is appeared, it means the test is working properly and the result is valid.

## COMPONENTS


- ① Test device ② Specimen extraction tube ③ Filter cap ④ Specimen collection swab ⑤ Instructions for use  
\*Control swabs are optional.

## SPECIMEN COLLECTION AND PREPARATION

### 01. Specimen collection

#### (01) Nasal swab specimens

To collect nasal swab specimens, insert the sterile swab into the nostril that presents the most secretion under visual inspection. Using gentle rotation push the swab until resistance is met at the level of the turbinate(less than one inch into the nostril). Rotate the swab a few times against the nasal wall.

Nasal swab

Throat swab

#### (02) Throat swab specimens

To collect throat swab specimens, rub the sterile swab on both tonsillar surfaces and the posterior pharynx.


### 02. Specimen transport and storage

Testing should be performed immediately after the specimens have been collected. The specimens extracted from swab may be stored at 1-8°C for up to one day. For long term storage, specimens should be kept below -20°C.

## TEST PROCEDURE

01. Insert the patient sample swab into extraction tube and swirl the swab 10 times.
02. Remove the swab while pressing the swab head against tests tube to extract the liquid from the swab.
03. Discard the used swab.
04. Close the tube with filter cap.
05. Take out test device from pouch and place it on a flat surface.
06. Add 3 drops into the sample well(S) of the test device.
07. After adding sample, reaction will be finished within 8 minutes.

## Test Procedure

## READING AND INTERPRETATION OF RESULT

01. Control band(C) means that the test is working properly.
02. Test band(A,B) indicates the test results.

### 01. NEGATIVE


The presence of only control band indicates a negative result.

### 02. POSITIVE

(1) Positive for Influenza virus A

The presence of two color bands ("A" band and "C" band) indicates a positive for Influenza A.

(2) Positive for Influenza virus B

The presence of two color bands ("B" band and "C" band) indicates a positive for Influenza B.

(3) Positive for Influenza virus A&B

The presence of three color bands ("A", "B" band and "C" band) indicates a positive for Influenza A&B.

### 03. INVALID

If control band("C" band) is not appeared after performing test, the result is considered invalid. The directions may not have been followed correctly or the device may have been deteriorated. It is recommended that the specimen be re-tested with the new device.

## FOLLOW-UP ACTION

01. BioTracer™ Influenza A&B Rapid Card is a qualitative test kit. For knowing the quantitative concentration of Influenza, it should be performed with additional quantitative test.
02. As with all diagnostic tests, a definitive clinical diagnosis should not be based on the results of a single test, but should only be made by the physician after all clinical and laboratory findings have been evaluated.

## LIMITATIONS OF THE METHOD

01. BioTracer™ Influenza A&B Rapid Card is for the qualitative detection of influenza A and B viral antigens. The test performance depends on antigen load and may not correlate with cell culture performed on the same specimen. Negative test results are not intended to rule out other non-influenza viral infections.
02. Sensitivity can differ with various strains of influenza due to difference in antigen expression. Specimens might contain new, non-identified strains of influenza that express varying amounts of antigen.
03. As with all diagnostic tests, the results obtained with this test kit yield data that must be used only as an adjust to other information available to the physician.
04. Individuals who received nasally administered influenza vaccine may have positive test results for up to three days after vaccination.
05. Monoclonal antibodies may fail to detect(or with less sensitivity) influenza A viruses that have undergone minor amino acid changes in the target epitope region.

## PERFORMANCE CHARACTERISTICS

### 01. Sensitivity and Specificity

A study was performed using 400 positive and negative specimens. Each specimen was assayed with the BioTracer™ Influenza A&B Rapid Card and a commercially available DFA after cell culture according to the respective instructions for use. Compared to culture and confirmed for influenza A or influenza B by DFA, the BioTracer™ Influenza A&B Rapid Card correctly identified 95.5 % (191/200) positive specimens and 100 % (200/200) negative specimens.

		Influenza A&B Rapid Card			Total
		A	B	Negative	
Cell Culture Test	A	97	-	5	102
	B	-	94	4	98
	Neg.	-	-	200	200
Total		97	94	209	400

※ Sensitivity : Influenza A 95.10 % (97/102), Influenza B 95.92 % (94/98)  
※ Specificity : more than 99.9% (200/200)

### 02. Analytical specificity and cross reactivity

BioTracer™ Influenza A&B Rapid Card was evaluated with 8 bacterial and 8 viral isolates. These studies had been performed to demonstrate that have been no cross reactivity with 8 bacterial panels and 8 virus panels on the BioTracer™ Influenza A&B Rapid Card.

### 03. Analytical sensitivity

Analytical sensitivity was established using a total 9 strains of influenza virus : 6 influenza A and 3 influenza B.

### 04. Precision

- (1) Within run performance test was determined by one analyst with ten devices each of 3 different lots containing different concentration every influenza type. There was no variation between the test devices of each lot and between three different lots.
- (2) Between run performance test was determined by three analysts with 3 different lots containing different concentration every influenza type. There was no variation between different analysts.

### 05. Analytical Interference test

The following potential interfering substances were tested and were found to have no affect on the performance of the BioTracer™ Influenza A&B Rapid Card.

Potential interfering substances	Concentration	Potential interfering substances	Concentration
Acetyl Salicylic Acid	20 mg/mL	Phenylephrine HCl	1 mg/mL
Acetamidophenol	10 mg/mL	Fexofenadine	500 mg/mL
Chlorpheniramine	5 mg/mL	Amantadine	500 mg/mL
Diphenhydramine HCl	5 mg/mL	Ribavirin	500 mg/mL
Guaiacol Glyceryl Ether	20 mg/mL	Pseudoephedrine HCl	20 mg/mL
Oxymetazoline HCl	0.05 mg/mL	Ibuprofen	10 mg/mL

## WARNING AND PRECAUTIONS

01. For in vitro diagnostic use only.
02. Do not eat or smoke while handling specimens.
03. Wear protective gloves while handling specimens. Wash hands thoroughly afterwards.
04. Do not use test kit if the packings damaged or the seal is broken.
05. Avoid splashing or aerosol formation.
06. Clean up spills thoroughly using an appropriate disinfectant.
07. Decontaminate and dispose of all specimens, reaction kits and potentially contaminated materials, as if they were infectious waste, in a bio hazard container.
08. Do not re-use the test device.
09. Use the only swab provided (Do not use the swab not provided).

## STORAGE AND SHELF-LIFE


BioTracer™ Influenza A&B Rapid Card should be stored at 2-30 °C (36-86 °F). The test device is sensitive to humidity as well as to heat. Do not use it beyond the expiration date, 24 months from date of manufacture.

## LITERATURE REFERENCES

01. Barenfanger, J., C. Drake, T. Mueller, T. Troutt, J. O'Brien, and K. Guttman. R-Mix cells are faster, at least as sensitive and marginally more costly than conventional cell lines for the detection of respiratory viruses. *J. Clin. Virol.* 2001; 22:101-110.
02. Glezen W. Emerging infections: pandemic influenza. *Epidemiol. Rev.* 1996;18:64-72.
03. Hayden FG, Treanor JJ, Betts RF, Lobo M, Esinhart JD, Hussey EK. Safety and efficacy of the neuraminidase inhibitor GG167 in experimental human influenza. *J. Am. Med. Assoc.* 1996;275:295-9.

## Manufactured by

NanoEnTek Inc. | 851-14, Seohae-ro, Paltan-myeon, Hwaseong-si, Gyeonggi-do, 18531, Korea  
Tel : +82-2-6220-7942 | Fax : +82-2-6220-7999 | ivdst@nanoentek.com | www.nanoentek.com

MT Promedt Consulting GmbH Altenhofstrasse 80, 66386 St. Ingbert, Germany