

Iqon Low DNA Ladder

Range: 100 bp – 1000 bp

MADE IN DENMARK



Cat. No.: A610541

A610541

100 lanes of 5 µl

-	Iqon Low DNA Ladder
ID No.	5401600
Cap colour	Black
Content	1 x 0.5 ml

Features and General Description

Iqon Low DNA Ladder is a broad range dsDNA ladder with bands from 100 bp to 1000 bp (figure 1).

The ladder is supplied in loading buffer, ready-to-use on agarose and polyacrylamide gels. It is suitable with TBE, TAE, SB and LB electrophoresis systems.

The 500 bp band of Iqon Low DNA Ladder comes with a higher intensity and can be used as a reference point.

Suggestions for use of Iqon PCR Ladder

- **Important:** Mix ladder briefly before use. Do not heat the ladder.
- Load 5 µl of ladder per lane (for gel wells of ~5 – 10 mm width).
- If band intensity is too high, then load a smaller volume of the ladder per lane.
- Agarose gel electrophoresis: Prepare 1-1.5 % gel. The dye should migrate 60 – 70 % the length of the gel.
- Polyacrylamide gel electrophoresis: Prepare 8 % gel. The dye should migrate approx. 90 % the length of the gel.
- Ethidium bromide (0.5 µg/ml) is the recommended gel stain.

Storage, Stability and Shipment

Iqon Low DNA ladder may be kept safely at room temperature for at least 6 months from date of shipment. Iqon DNA Ladders are guaranteed for 12 months when stored at 4°C. Aliquot product if necessary, to avoid repeated freezing and thawing cycles.

Shipped at ambient temperature. For long term storage, keep the Iqon Low DNA Ladder at -20 °C for up to 3 years.

Quality control

Agarose gel analysis shows that all bands are present at the expected location and band intensity.




Figure 1: Iqon Low DNA Ladder. 5 µl Iqon PCR Ladder was loaded on a 1.5 % agarose in 1x TBE and stained with ethidium bromide.

Related Products

DNA Ladder	Cat. No.
Low Range DNA Ladder: 100 -1000 bp. 250 lanes	A610241
PCR DNA Ladder: 100 -3000 bp. 250 lanes	A610341
Iqon Mini DNA Ladder: 100 - 500 bp. 100 lanes	A610441
Iqon Low DNA Ladder: 100 - 1000 bp. 100 lanes	A610541
Iqon PCR Ladder: 100 - 3000 bp. 100 lanes	A610641

Loading Buffers	Cat. No.
5x Loading Buffer Red, 5 x 1 ml	A608104
5x Loading Buffer Blue, 5 x 1 ml	A608204
5x Loading Buffer Orange, 5 x 1 ml	A608304
5x Loading Buffer Cyan, 5 x 1 ml	A608404

Ultrapure dNTPs*	Cat. No.
dNTP Mix 40 mM (2 x 500 µl): 10 mM each dA, dC, dG, dT	A502004
dNTP Set, 100 mM each: 250 µl of each dA, dC, dG and dT	A511104

*Other concentrations and Single dNTPs are available.

PCR Grade Water	Cat. No.
6 x 5 ml	A360056

Taq Polymerase (500 units)	Cat. No.
Taq DNA Polymerase 5 U/µl	A110003
• with 10x Ammonium Buffer	A111103

AQ97 High Fidelity DNA Polymerase 2 U/µl*	Cat. No.
With 5x AQ97 Buffer – 500 U	A767503

* Available in different size ranges.

Hot Start Polymerase (500 units)	Cat. No.
TEMPase Hot Start DNA Polymerase, 5 U/µl	A220003
• with 10x Ammonium Buffer	A221103

AQ97 High Fidelity DNA Polymerase 2x Master Mix*	Cat. No.
100 reactions	A770101

* Available in different size ranges.

Taq Master Mixes (500 x 50 µl reactions)	Cat. No.
2x Master Mix, 1.5 mM MgCl ₂ final concentration	A140303
2x Master Mix RED, 1.5 mM MgCl ₂ final concentration	A180303

TEMPase Hot Start Master Mixes (500 x 50 µl reactions)	Cat. No.
2x Master Mix A**, 1.5 mM MgCl ₂ final concentration	A230303
2x Master Mix A**BLUE, 1.5 mM MgCl ₂ final concentration	A290403

Real-time PCR Master Mixes (400 x 25 µl reactions)	Cat. No.
RealQ Plus 2x Master Mix for probe,	
• without ROX TM	A313402
• with low ROX TM	A314402
• with high ROX TM	A315402
RealQ Plus 2x Master Mix Green	
• without ROX TM	A323402
• with low ROX TM	A324402
• with high ROX TM	A325402

For Research Use Only. Not for use in diagnostics procedures.

Other product sizes, combinations and customized solutions are available. Please look at www.ampliqon.com or ask for our complete product list for PCR Enzymes. For customized solutions please contact us.

Made in Denmark

Issued 08/2021

Iqon DNA Ladders

