***Instruction Manual
Manuale di istruzioni
Manuel d'instructions
Manual de instrucciones
Bedienungsanleitung***

SE Radiation Detector

R10800340

General Information / Informazioni Generali / Informations Générales / Información General / Allgemeine Hinweise

Before using the unit, please read the following instruction manual carefully.

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Do not dispose of this equipment as urban waste, in accordance with EEC directive 2002/96/CE.

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This unit must be used for laboratory applications only.

The manufacturer declines all responsibility for any use of the unit that does not comply with these instructions.

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This unit has been designed and manufactured in compliance with the following standards:

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Safety requirements for electrical equipment for measurement, control and for laboratory use

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IEC/EN 61010-1

Electrical equipment for laboratory use

UL 61010-1

General requirement - Canadian electrical code

CAN/CSA-C22.2 No.61010-1

VELP reserves the right to modify the characteristics of its products with the aim to constantly improving their quality.

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Avoid exposing the instrument to liquid, moisture, corrosive gas, high temperature and direct sunlight for extended periods.

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Mica window tube can be easily damaged if struck directly. *Do not insert anything through the screen!*

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To avoid contamination, do not touch the instrument to the surface being tested.

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This instrument may be sensitive to and may not operate in radio frequency, microwave, electrostatic, and magnetic fields.

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The unit is susceptible to EMP (electromagnetic pulse) and may be rendered inoperable by an atomic detonation. It has not been determined what distance from an atomic blast would be considered safe for semiconductor circuitry.

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1. Introduction

The meter senses the four main types of ionizing radiation: alpha, beta, gamma, X-rays.

Calibrated by Cesium 137, is suitable to reveal also other sources of ionizing radiations. Shows maximum accuracy when measuring Cesium 137 and isotopes of similar energies; some isotopes which can be revealed are Cobalt 60, Phosphorus 32, Strontium 90, and different forms of Radium, Plutonium, Uranium and Thorium.

Gamma and X radiations are measured in milli-Roentgens/hour (mR/hr), micro-Sieverts ($\mu\text{Sv/hr}$) or milli-Sieverts (mSv/hr). Alpha and Beta are measured in counts/minute (CPM) or counts/second (CPS).

The revelation of ionizing radiations is performed by means of a Geiger-Mueller tube having a thin mica window equipped with a screen. Some levels of alpha or beta, gamma and X rays, having a low energy, which cannot penetrate the plastic covering of the instrument or the GM tube walls, can be revealed through the window.

The GM tube is placed inside the instrument.

Fig.1

When a ray or particle of ionizing radiation enters the GM tube or passes through it, it is electronically revealed and visualized by means of a red count light.

If the switch is in AUDIO position, the instrument emits also an acoustic signal every time a radioactive event occurs.

5-25 counts/minute about, with random intervals, in relation to position and altitude, are due to natural background radiation.

2. Assembly and installation

Check the integrity of the unit after unpacking. The box includes:

- SE Radiation detector
- Instruction manual

2.1 Start-up

Before turning on the instrument, install a 9 V alkaline battery. If a battery is already installed, turn the instrument on and switch the range switch to the BATT position: battery condition will be indicated on the meter. If you don't think that high levels of radiation are present, set the range switch in the X1 position. If the meter goes off scale, move the range switch to the next higher setting, X10 or X100. For an audible signal, position the OFF/ON/ AUDIO switch to the AUDIO position.

NOTE: The flashes from the count light and the audible beeps are progressively shorter in X10, X100 positions. At the end of the work cycle set the stirring speed knob on "0" rpm, the timer knob on "OFF" and turn off the light of central column. If the instrument is not used for a long period it's better to disconnect the external power supply from the plug.

3. Operating controls

To determine whether the radiation detected is alpha, beta or gamma, move the top of the instrument toward the source.

Gamma: if there is an indication of radioactivity, it is most likely gamma or high energy beta.

Low energy gamma and x-rays (10-40 KeV) cannot penetrate the sidewall of the GM tube, but may be detected through the window.

Beta: place a piece of aluminum about 1/8" (3mm) thick between the instrument and the source.

If the indication stops, decreases or changes, it is most likely beta radiation.

Most common isotopes emit both beta and gamma radiation.

For this reason, the signal could decrease or change, but not stop totally.

Alfa: if there is no indication through the back of the case (the side of the tube GM), position the window close to but not touching the source. If there is an indication, it is alpha, beta or low energy gamma.

If a sheet of paper placed between the window and the source stops the indication, it is most likely alpha.

Do not hold the source above the window, to avoid particles falling into the instrument.

If you perform the alpha/beta test above and there is no change or only a very slight change in the indication, the source is emitting primarily gamma radiation.

4. Maintenance

No routine or extraordinary maintenance is necessary apart from periodically cleaning the unit as described in this manual. In compliance with the product guarantee law, repairs to our units must be carried out in our factory, unless previously agreed otherwise with local distributors. The instrument must be transported in its original packaging and any indications present on the original packaging must be followed (e.g. palletized).

It is the responsibility of the user, to properly decontaminate the unit in case of hazardous substances remaining on the surface or interior of the device. If in doubt about the compatibility of a cleaning or decontamination product, contact the manufacturer or distributor.

4.1 Cleaning

Use a cloth dampened with an non-inflammable non-aggressive detergent.

5. Technical data

Detector

Geiger-Mueller LND 712 tube, uncompensated, halogen quenched. Thin mica window, 1,5-2,0 mg/cm² thick. About 1000 CPM/mR/hr for Cesium 137.

Energy sensitivity

- detects alpha down to 2.5 MeV; typical detection efficiency at 3,6 MeV is greater than 80%.
 - detects beta at 50 KeV with typical 35% detection efficiency, at 150 KeV with 75% detection efficiency.
 - detects gamma and X rays down to 10 KeV through the window, 40 KeV minimum through the case.
- Normal background is 5-20 CPM.

Operating range

0-50 mR/hr and 0-50000 CPM or 0-500 µSv/hr and 0-50 mR/hr , analog meter with dual scale.

Range switch

X1, X10, X100; battery check. When the switch is in position X1, the instrument reveals from 0 to 0,5 mR/hr, corresponding to 0-500 CPM or to 0-5 µSv/hr. In the position X10 readings are to be multiplied X10, indicating values up to 5mR/hr or 5000 CPM or 50 µSv/hr. In the position X100, readings are to be multiplied X100, indicating values up to 50 mR/hr or 50000 CPM or 500 µSv/hr, which correspond to about 2500 times the natural background radiation.

Accuracy

±15% of full scale (referenced to Cesium 137).

Audio

built-in piezoelectric transducer gives audible beep (3 KHz frequency about) when switch is in the AUDIO position. Can be switched off for silent.

Anti-saturation

Meter will hold at full scale (X100), in fields as high as 100 times the maximum reading.

Operating voltage

7-11 V , CC.

High voltage

Fully regulated, having a variation less than 3V in the outlet for a variation of supply from 9 to 6V.3.10

Power requirements

1 alkaline battery 9V. Battery life is up to 2000 hours at normal background radiation levels.

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0.22 Kg, ohne Batterie; 70x210x48 mm.

Thank you for having chosen VELP!

Established in 1983, VELP is today one of the world's leading manufacturer of analytical instruments and laboratory equipment that has made an impact on the world-wide market with Italian products renowned for innovation, design and premium connectivity. VELP works according to **ISO 9001**, **ISO14001** and **OHSAS 18001** Quality System Certification.

Our instruments are manufactured in Italy according to the IEC 1010-1 and CE regulation.

Our product lines:

Analytical instruments

Elemental Analyzers
Digestion Units
Distillation Units
Solvent Extractors
Fiber Analyzers
Dietary Fiber Analyzers
Oxidation Stability Reactor
Consumables

Laboratory Equipment

Magnetic Stirrers
Heating Magnetic Stirrers
Heating Plates
Overhead stirrers
Vortex mixers
Homogenizers
COD Thermoreactors
BOD and Respirometers
Cooled Incubators
Flocculators
Overhead Shakers
Turbidimeter
Radiation Detector
Open Circulating Baths
Pumps

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VELP Scientifica Srl