


*insta***BAG**[®]

Filter blender bag with dehydrated media

Your diluent ready in minutes!



NEW

Available in BPW and Fraser 1/2

interscience

Ready-to-use media for easy preparation

instaBAG[®] is a filter bag which includes a dehydrated pre-dosed media, allowing a 2-in-1 process. Skip the time-consuming media preparation.



1 Place the sample in an *instaBAG*[®]



2 Add sterile and purified water



3 Blend the *instaBAG*[®] with the BagMixer[®]

➔ The sample is ready to be pipetted

All the advantages of a filter bag

Debris are blocked by the filter during blending

Long shelf life

Easy pouch opening

Zip pouch with desiccant: optimal conservation after opening

Light and humidity proof aluminum pouch



Just add water

Ready-to-use bag with pre-dosed dehydrated media

Complete dissolution during blending

Non-woven filter identical to the bag filter

Available in BPW 90/225 mL and Fraser 225 mL

PATENTED

SHELF LIFE
> 12 MONTHS

STORAGE
10°C to 30°C

Space saving



100
instaBAG



100
glass bottles

No waste



Less garbage



instaBAG
+ diluent + sample



Media bottle
+ bag + diluent + sample

Easy preparation



with *instaBAG*.



with media preparation

Green

instaBAG: no water transported




Diluent pouches



Long shelf life

instaBAG



Media bottle (BPW)



Technical specifications



	<i>instaBAG</i> [®] EPT 90	<i>instaBAG</i> [®] EPT 225	<i>instaBAG</i> [®] Fraser 1/2
Reference	114 090	114 225	114 220
Type of bag	Bag with lateral filter		
Composition of the bag	MultiLayer [®] : reinforced multicoated complex		
Type of filter	Non-woven		
Filter porosity	< 250 micrometers		
Type of dehydrated media	Buffered peptone water		Supplemented Fraser 1/2
Weight of dehydrated media	1.8 g	4.5 g	13 g
Sample weight	10 g (± 5 %)	25 g (± 5 %)	
Volume of water to be added (sterile and deionized)	90 mL (± 2 %)	225 mL (± 2 %)	
Blending time	1 mn		2 mn
Bag dimensions	190 x 300 mm		
Aluminum pouch of	10 bags		
Shelf life	24 months		12 months
Storage conditions	Temperature from + 6°C to + 30°C		Temperature from + 10°C to + 30°C
Treatment	Gamma ray (10 to 25 KGy)		
Standards	ISO 11133 - ISO 11290-2 - ISO 21528-1 - ISO 7218 - ISO 6579 - ISO 6887- FDA-BAM		ISO 11133 - ISO 11290-2 - ISO 7218 - ISO 6887- FDA-BAM
Box of	100 bags		80 bags

MSDS available on request. Each batch comes with a complete ISO 11133 certificate, controlling fertility, selectivity and sterility.

Ask for free samples

Certified production

Products made for INTERSCIENCE by Interlab, an ISO 9001 certified company.



The winning combination



DiluFlow[®] Gravimetric dilutors



instaBAG[®]



BagMixer[®] Lab blenders

Your local distributor

interscience

PARIS

Phone: +33 (0)1 34 62 62 61 - Email: info@interscience.com

FRANKFURT

Phone: +49 611 7238 7770 - Email: info@interscience.com

BOSTON

Phone: +1 781 937 0007 - Email: sales.usa@intersciencelab.com

SHANGHAI

Phone: +86 (0)21-64739390 - Email: sales.china@interscience.com

SINGAPORE

Phone: +65 6977 7232 - Email: sales.asia@interscience.com

TOKYO

Phone: +81 3 6712 9715 - Email: sales.japan@interscience.com

www.interscience.com