

Lovibond® Water Testing

Tintometer® Group



Spectrophotometer XD 7500 (UV-VIS)

Simplify your Work - Save Time and Money!



- Premium optical system with reference beam
- Automatic test recognition with internal barcode reader
- Automatic cuvette type detection
- Support of analytical quality assurance

Part Number: 71307500

Latest technology

The XD 7000 and XD 7500 spectrophotometers are equipped with VIS respectively UV / VIS reference beam technology.

Automatic test recognition

Each of the more than 150 pre-programmed Lovibond® methods is recognized by barcode as well as the cuvette type.

Analytical quality assurance

Standard procedures for analytical quality assurance support the verification of the photometer, the overall system (including the chemical methodology) and the verification of matrix effects.

Built-in security levels

The devices support the assignment of passwords and the assignment of up to three different levels of rights.

Additional functions

Measurement of transmission and absorption, spectral scan, kinetics analysis, as well as the creation of user-defined methods.

Interfaces for data processing

You want to process your data? With Ethernet, USB B, USB A for external memory, keyboard, barcode scanner and printer different versions are available.

Industry

Chemical Industry | Food and Beverage Industry | Industries Others | Marine Industry | Municipalities | NGO | Oil and Gas | Pharmaceutical Industry | Power and Energy

Application

Boiler Water | Cooling Water | Disinfection Control | Drinking Water Treatment | Food and Beverage | Galvanization | Others | Pool-Water Control | Pool Water Treatment | Raw Water Treatment | Waste Water Treatment

Spectrophotometer XD 7500 (UV-VIS)

The XD series instruments are ideal for routine and spectral analysis. As true all-rounders, the XD 7000 and XD 7500 facilitate work routines, especially with automatic method recognition using barcodes, automatic cell recognition and more than 150 pre-programmed methods. In addition to analytical quality control, they also support Good Labor Practice (GLP). The instruments are available together with bar-coded test tubes and a wide range of accessories from a single source - which is also easy on your wallet. High-quality reference beam optics, usability, easy handling and multilingual global applicability round off the versatility of the VIS and UV/VIS spectrophotometers, which are also equipped for flexible on-site use.

Measuring Range

| Test Name | Measuring Range | Chemical Method |
|------------------------------------|---|--|
| Alkalinity-m HR T | 5 - 500 mg/L CaCO ₃ | Acid / Indicator |
| Alkalinity-m T | 5 - 200 mg/L CaCO ₃ | Acid / Indicator |
| Alkalinity-p T | 5 - 500 mg/L CaCO ₃ | Acid / Indicator |
| Aluminium PP | 0.01 - 0.25 mg/L Al | Eriochrom Cyanine R |
| Aluminium T | 0.01 - 0.3 mg/L Al | Eriochrom Cyanine R |
| Ammonia HR TT | 1.0 - 50 mg/L N | Salicylate |
| Ammonia LR TT | 0.02 - 2.5 mg/L N | Salicylate |
| Ammonia PP | 0.01 - 0.8 mg/L N | Salicylate |
| Ammonia T | 0.02 - 1 mg/L N | Indophenole Blue |
| Arsenic | 0.02 - 0.6 mg/L As | Silver Diethyldithiocarbamate |
| Bromine 10 T | 0.1 - 3 mg/L Br ₂ | DPD |
| Bromine 50 T | 0.05 - 1 mg/L Br ₂ | DPD |
| Bromine PP | 0.05 - 4.5 mg/L Br ₂ | DPD |
| Bromine T | 0.05 - 13 mg/L Br ₂ | DPD |
| Cadmium M. TT | 0.025 - 0.75 mg/L Cd | Cadion |
| Chloramine (M) PP | 0.02 - 4.5 mg/L NH ₂ Cl as Cl ₂ | Indophenole method |
| Chloride L (A) | 5.00 - 60 mg/L Cl ⁻ | Iron(III)-thiocyanate |
| Chloride L (B) | 0.5 - 20 mg/L Cl ⁻ | Mercury Thiocyanate / Iron Nitrate |
| Chloride T | 0.5 - 25 mg/L Cl ⁻ | Silver Nitrate / Turbidity |
| Chlorine 10 T | 0.1 - 6 mg/L Cl ₂ | DPD |
| Chlorine 50 T | 0.02 - 0.5 mg/L Cl ₂ ^{a)} | DPD |
| Chlorine (free) and Monochloramine | 0.02 - 4.50 mg/L Cl ₂ | Indophenole method |
| Chlorine dioxide 50 T | 0.05 - 1 mg/l ClO ₂ | DPD / Glycine |
| Chlorine dioxide PP | 0.04 - 3.8 mg/l ClO ₂ | DPD |
| Chlorine dioxide T | 0.02 - 11 mg/l ClO ₂ | DPD / Glycine |
| Chlorine HR 10 T | 0.1 - 10 mg/L Cl ₂ ^{a)} | DPD |
| Chlorine HR (KI) T (105) | 5 - 200 mg/L Cl ₂ | KI / Acid |
| Chlorine L | 0.02 - 4.0 mg/L Cl ₂ ^{a)} | DPD |
| Chlorine MR PP | 0.02 - 3.5 mg/L Cl ₂ ^{a)} | DPD |
| Chlorine PP | 0.02 - 2 mg/L Cl ₂ ^{a)} | DPD |
| Chlorine T | 0.01 - 6.0 mg/L Cl ₂ ^{a)} | DPD |
| Chromium 50 PP | 0.005 - 0.5 mg/L Cr ^{b)} | Diphenylcarbazide |
| Chromium PP | 0.02 - 2 mg/L Cr ^{b)} | Diphenylcarbazide |
| COD HR TT | 200 - 15000 mg/L COD ^{b)} | Dichromate / H ₂ SO ₄ |
| COD LMR TT | 15 - 300 mg/L COD ^{b)} | Dichromate / H ₂ SO ₄ |
| COD LR TT | 3 - 150 mg/L COD ^{b)} | Dichromate / H ₂ SO ₄ |
| COD MR TT | 20 - 1500 mg/L COD ^{b)} | Dichromate / H ₂ SO ₄ |
| Copper 50 T | 0.05 - 1 mg/L Cu ^{a)} | Biquinoline |
| Copper L | 0.05 - 4 mg/L Cu ^{a)} | Bicinchoninate |
| Copper PP | 0.05 - 5 mg/L Cu | Bicinchoninate |
| Copper T | 0.05 - 5 mg/L Cu ^{a)} | Biquinoline |
| CyA HR T | 10 - 200 mg/L CyA | Melamine |
| Cyanide 50 L | 0.005 - 0.2 mg/L CN ⁻ | Pyridine-barbituric Acid |
| Cyanide L | 0.01 - 0.5 mg/L CN ⁻ | Pyridine-barbituric Acid |
| CyA T | 10 - 160 mg/L CyA | Melamine |
| DEHA PP | 0.02 - 0.5 mg/L DEHA | PPST |
| DEHA T (L) | 0.02 - 0.5 mg/L DEHA | PPST |
| Fluoride L | 0.05 - 2 mg/L F ⁻ | SPADNS |
| Formaldehyde 10 M. L | 1.00 - 5.00 mg/L HCHO | H ₂ SO ₄ / Chromotropic acid |
| Formaldehyde 50 M. L | 0.02 - 1.00 mg/L HCHO | H ₂ SO ₄ / Chromotropic acid |
| Formaldehyde M. TT | 0.1 - 5 mg/L HCHO | H ₂ SO ₄ / Chromotropic acid |

| Test Name | Measuring Range | Chemical Method |
|------------------------------------|---|--|
| H ₂ O ₂ 50 T | 0.01 - 0.5 mg/L H ₂ O ₂ | DPD / Catalyst |
| H ₂ O ₂ HR L | 40 - 500 mg/L H ₂ O ₂ | Titanium Tetrachloride / Acid |
| H ₂ O ₂ LR L | 1 - 50 mg/L H ₂ O ₂ | Titanium Tetrachloride / Acid |
| H ₂ O ₂ T | 0.03 - 3 mg/L H ₂ O ₂ | DPD / Catalyst |
| Hardness Ca and Mg L | 0.05 - 4 mg/L CaCO ₃ | Calmagite |
| Hardness Ca and Mg MR TT | 10 - 360 mg/L CaCO ₃ | Calmagite |
| Hardness Calcium (B) T | 20 - 500 mg/L CaCO ₃ | Murexide |
| Hardness Calcium (B) T | 50 - 900 mg/L CaCO ₃ | Murexide |
| Hardness total HR T | 20 - 500 mg/L CaCO ₃ ⁱ⁾ | Metallphthaleine |
| Hardness total T | 2 - 50 mg/L CaCO ₃ | Metallphthaleine |
| Hazen 24 | 10 - 500 mg/L Pt | (APHA) Platinum Cobalt Standard Method |
| Hazen 50 | 10 - 500 mg/L Pt | (APHA) Platinum Cobalt Standard Method |
| Hydrazine C | 0.01 - 0.7 mg/L N ₂ H ₄ ^{c)} | PDMAB |
| Hydrazine L | 5 - 600 µg/L N ₂ H ₄ | Dimethylaminobenzaldehyde |
| Hydrazine P | 0.05 - 0.5 mg/L N ₂ H ₄ | Dimethylaminobenzaldehyde |
| Hypochlorite T | 0.2 - 17 % NaOCl | Potassium Iodide |
| Iron 10 T | 0.05 - 1 mg/L Fe | Ferrozine / Thioglycolate |
| Iron 50 PP | 0.01 - 1.5 mg/L Fe ⁹⁾ | 1,10-Phenanthroline |
| Iron 50 T | 0.01 - 0.5 mg/L Fe | Ferrozine / Thioglycolate |
| Iron (TPTZ) PP | 0.02 - 1.8 mg/L Fe | TPTZ |
| Iron HR L | 0.1 - 10 mg/L Fe | Thioglycolate |
| Iron in Mo PP (224) | 0.01 - 1.8 mg/L Fe | TPTZ |
| Iron LR L (A) | 0.03 - 2 mg/L Fe | Ferrozine / Thioglycolate |
| Iron LR L (B) | 0.03 - 2 mg/L Fe | Ferrozine / Thioglycolate |
| Iron PP | 0.02 - 3 mg/L Fe ⁹⁾ | 1,10-Phenanthroline |
| Iron T | 0.02 - 1 mg/L Fe | Ferrozine / Thioglycolate |
| K _{S4.3} T | 0.1 - 4 mmol/L K _{S4.3} | Acid / Indicator |
| Lead 10 | 0.1 - 5 mg/L Pb | 4-(2-Pyridylazo)-resorcin |
| Lead (A) TT | 0.1 - 5 mg/L Pb | 4-(2-Pyridylazo)-resorcin |
| Lead (B) TT | 0.1 - 5 mg/L Pb | 4-(2-Pyridylazo)-resorcin |
| Iodine T | 0.05 - 3.6 mg/L I | DPD |
| Manganese HR PP | 0.1 - 18 mg/L Mn | Periodate Oxidation |
| Manganese L | 0.05 - 5 mg/L Mn | Formaloxime |
| Manganese LR PP | 0.01 - 0.7 mg/L Mn | PAN |
| Manganese T | 0.2 - 4 mg/L Mn | Formaloxime |
| Molybdate HR L | 1 - 100 mg/L MoO ₄ | Thioglycolate |
| Molybdate HR PP | 0.3 - 40 mg/L Mo | Mercaptoacetic Acid |
| Molybdate LR PP | 0.03 - 3 mg/L Mo | Ternary Complex |
| Molybdate T | 1 - 50 mg/L MoO ₄ | Thioglycolate |
| Nickel 50 L | 0.02 - 1 mg/L Ni | Dimethylglyoxime |
| Nickel L | 0.2 - 7 mg/L Ni | Dimethylglyoxime |
| Nitrate HR | 1.2 - 35 mg/L N | 2,6-Dimethylphenole |
| Nitrate LR TT | 0.5 - 14 mg/L N | 2,6-Dimethylphenole |
| Nitrate T | 0.08 - 1 mg/L N | Zinc Reduction / NED |
| Nitrate TT | 1 - 30 mg/L N | Chromotropic Acid |
| Nitrite HR PP | 2 - 250 mg/L NO ₂ ⁻ | Ferrous Sulfate Method |
| Nitrite HR TT | 0.3 - 3 mg/L N | Sulfanilic / Naphthylamine |
| Nitrite LR TT | 0.03 - 0.6 mg/L N | Sulfanilic / Naphthylamine |
| Nitrite PP | 0.01 - 0.3 mg/L N | Diazotation |
| Nitrite T | 0.01 - 0.5 mg/L N | N-(1-Naphthyl)-ethylenediamine |
| Nitrite VHR L | 25 - 2500 mg/L NO ₂ ⁻ | Ferrous Sulfate Method |

| Test Name | Measuring Range | Chemical Method |
|-----------------------------|---|-------------------------------------|
| Oxygen active T | 0.1 - 10 mg/L O ₂ | DPD |
| Oxygen dissolved C | 10 - 1100 µg/L O ₂ ^{e)} | Rhodazine D TM |
| Ozone 50 T | 0.02 - 0.5 mg/L O ₃ | DPD / Glycine |
| Ozone PP | 0.015 - 2 mg/L O ₃ | DPD / Glycine |
| Ozone T | 0.02 - 2 mg/L O ₃ | DPD / Glycine |
| Phenol T | 0.1 - 5 mg/L C ₆ H ₅ OH | 4-Aminoantipyrine |
| PHMB T | 2 - 60 mg/L/PHMB | Buffer / Indicator |
| Phosphate h. TT | 0.02 - 1.6 mg/L P ^{b)} | Phosphomolybdenum Blue |
| Phosphate HR C | 1.6 - 13 mg/L P ^{c)} | Vanadomolybdate |
| Phosphate HR L | 5 - 80 mg/L PO ₄ | Vanadomolybdate |
| Phosphate HR T | 0.33 - 26.09 mg/L P | Vanadomolybdate |
| Phosphate HR TT | 0.98 - 19.57 mg/L P | Vanadomolybdate |
| Phosphate LR C | 0.016 - 1.6 mg/L P ^{c)} | Stannous Chloride |
| Phosphate LR L | 0.1 - 10 mg/L PO ₄ | Phosphomolybic Acid / Ascorbic Acid |
| Phosphate LR T | 0.016 - 1.305 mg/L P | Phosphomolybdenum Blue |
| Phosphate PP | 0.02 - 0.815 mg/L P | Phosphomolybdenum Blue |
| Phosphate t. TT | 0.02 - 1.1 mg/L P ^{b)} | Phosphomolybdenum Blue |
| Phosphate total HR TT | 1.5 - 20 mg/L P ^{b)} | Phosphomolybdenum Blue |
| Phosphate total LR TT | 0.07 - 3 mg/L P ^{b)} | Phosphomolybdenum Blue |
| Phosphate TT | 0.02 - 1.63 mg/L P | Phosphomolybdenum Blue |
| Phosphonate PP | 0.02 - 125 mg/L PO ₄ | Persulfate UV Oxidation Method |
| pH-value HR T | 8.0 - 9.6 pH | Thymol Blue |
| pH value L | 6.5 - 8.4 pH | Phenol Red |
| pH-value LR T | 5.2 - 6.8 pH | Bromocresolpurple |
| pH-value T | 6.5 - 8.4 pH | Phenol Red |
| Polyacrylate L | 1 - 30 mg/L Polyacryl | Turbidity |
| Potassium T | 0.7 - 16 mg/L K | Tetraphenylborat Turbidity |
| SAC 254 nm (344) | 0.25 - 50 m ⁻¹ | Direct Reading EN ISO 7887:1994 |
| SAC 436 nm | 0.5 - 50 m ⁻¹ | Direct Reading EN ISO 7887:1994 |
| SAC 525 nm | 0.5 - 50 m ⁻¹ | Direct Reading EN ISO 7887:1994 |
| SAC 620 nm | 0.5 - 50 m ⁻¹ | Direct Reading EN ISO 7887:1994 |
| Selenium | 0.05 - 2 mg/L Se | 3,3'-Diaminobenzidine in Toluene |
| Silcate T | 0.05 - 4 mg/L SiO ₂ | Silicomolybdenum Blue |
| Silicate HR PP | 1 - 100 mg/L SiO ₂ | Silicomolybdate |
| Silicate L | 0.1 - 8 mg/L SiO ₂ | Heteropolyblue |
| Silicate LR PP | 0.05 - 1.6 mg/L SiO ₂ | Heteropolyblue |
| Silica VLR PP | 0.005 - 0.5 mg/L SiO ₂ | Heteropolyblue |
| Sulphate HR PP | 50 - 1000 | Bariumsulphate Turbidity |
| Sulphate PP | 5 - 100 mg/L SO ₄ ²⁻ | Bariumsulphate Turbidity |
| Sulphate T | 5 - 100 mg/L SO ₄ ²⁻ | Bariumsulphate Turbidity |
| Sulphide L | 8 - 1400 µg/L S ²⁻ | Methylene Blue |
| Sulphide T | 0.04 - 0.5 mg/L S ²⁻ | DPD / Catalyst |
| Sulphite 10 T | 0.1 - 12 mg/L SO ₃ | DTNB |
| Sulphite T | 0.1 - 6 mg/L SO ₃ | DTNB |
| Surfactants M. (anion.) TT | 0.05 - 2 mg/L SDSA | Methylene Blue |
| Surfactants M. (cation.) TT | 0.05 - 1.5 mg/L CTAB | Disulphine Blue |

| Test Name | Measuring Range | Chemical Method |
|-------------------------------|---|--|
| Surfactants M. (not ionic) TT | 0.1 - 7.5 mg/L Triton X-100 | TBPE |
| Suspended solids 24 | 10 - 750 mg/L TSS | Turbidity / Attenuated Radiation Method |
| Suspended solids 50 | 10 - 750 mg/L TSS | Turbidity / Attenuated Radiation Method |
| Tannin L | 0.5 - 20 mg/L Tannin | |
| TN HR 2 TT | 5 - 140 mg/L N ^{b)} | 2,6-Dimethylphenole |
| TN HR TT | 5 - 150 mg/L N ^{b)} | Persulphate Digestion |
| TN LR 2 TT | 0.5 - 14 mg/L N ^{b)} | 2,6-Dimethylphenole |
| TN LR TT | 0.5 - 25 mg/L N ^{b)} | Persulphate Digestion |
| TOC HR M. TT | 50 - 800 mg/L TOC ^{b)} | H ₂ SO ₄ / Persulphate / Indicator |
| TOC LR M. TT | 5 - 80 mg/L TOC ^{b)} | H ₂ SO ₄ / Persulphate / Indicator |
| Triazole PP | 1 - 16 mg/L Benzotriazole or Tolytriazole | Catalyzed UV Digestion |
| Turbidity 50 | 5 - 500 FAU | Attenuated Radiation Method |
| Turbidity 24 | 10 - 1000 FAU | Attenuated Radiation Method |
| Urea T | 0.1 - 2.5 mg/L Urea | Indophenol / Urease |
| Zinc L | 0.1 - 2.5 mg/L Zn | Zincon / EDTA |
| Zinc T | 0.02 - 1 mg/L Zn | Zincon |

Technical Data

| | |
|--------------------------------------|--|
| Optics | Grid monochromator with reference beam and beam splitter after exit slit |
| Measurement | Concentration, single and multi-wavelength measurement of absorbance and % transmission, kinetics, spectra |
| Wavelength Range | 190 - 1100 nm (nm) |
| Wavelength Resolution | 1 nm |
| Wavelength Accuracy | ± 1 nm on all Holmium peaks |
| Wavelength Reproducibility | better than 0,5 nm |
| Spectral Scope | 4 nm |
| Photometric Range | -3.3 - +3.3 Abs |
| Photometric Resolution | Absorption: 0.001 ; Transmission: 0.1 % |
| Photometric Accuracy | 0.003 Abs below 0.6 Abs ; 0.5 % from 0.6 to 2.0 Abs |
| Photometric Reproducibility | 0.003 Abs below 0.6 Abs ; 0.5 % from 0.6 to 2.0 Abs |
| Photometric Linearity | < 1 % up to 2.0 Abs between 340 to 900 nm |
| Scan Speed | 700 - 2000 nm/min. |
| Drift | < 0.005 Abs per hour after 15 minutes heat up time |
| Stray Light | < 0.05 % Transmission at 340 and 408 nm |
| Suitable Vials | Rectangular Cuvettes 10 mm Rectangular Cuvettes 20 mm Rectangular Cuvettes 50 mm Round Cuvettes 13 mm Round Cuvettes 16 mm Round Cuvettes 24 mm |
| Display | 7" high contrast colour graphic-display |
| Interfaces | Ethernet USB B USB A for External Memory Keypad Barcode-Scanner PCL Compatible Printer |
| Operation | Membrane Keyboard |
| Automatic Cuvette Recognition | Round cuvettes: 13, 16 and 24 mm ; Rectangular cuvettes: 10, 20 and 50mm |
| Auto – OFF | Yes |
| Test Recognition | via internal barcode reader |

| | |
|------------------------------------|---|
| Auto-Check | Self-test at each switch-on: test of memory, processor, internal interface, filter lamp and additional calibration of each wavelength |
| LIMS Compatibility | ASCII, .csv-files |
| Internal Storage | approx. 5000 data sets (method, user ID, date, result), autostorage function / manual storage function |
| Security | Password protection possible: 3 different user levels (guest, user, admin) |
| Power Drain | 100 - 240 V, 50/60 Hz |
| Power Supply | Buffer batteries (4 x AA), power supply unit with cable |
| Portability | Benchtop |
| Environmental Conditions | +10 °C to 35 °C (41 °F to 95 °F), ≤ 75 % average humidity within one year 95 %, max. 30 days/year, 85 % all other days |
| Stock Conditions | -25 °C to +65 °C (-13 °F to 268 °F) |
| Protection Class | IP 30 |
| Compliance | CE |
| IP Protection Class | EN 60529 |
| Interference Emission | Class B |
| Interference Immunity | IEC 61000-4-3 |
| Tolerance Extension | 0.008 E |
| Meter Safety | EC Directive 2014/35/EC EN 61010-1:2010 |
| Languages User Interface | German, English, French, Spanish, Italian, Portuguese, Polish, Indonesian, Russian, Chinese, Japanese, Dutch, Swedish, Norwegian, Czech, Romanian, Macedonian, Slovenian, Hungarian, Turkish, Korean, Vietnamese, Thai, Serbian, Malaysian, Danish, Bulgarian |
| Languages Quick Start Guide | German, English, French, Spanish, Italian, Portuguese, Polish, Indonesian, Russian, Chinese, Japanese, Dutch, Swedish, Norwegian, Czech, Romanian, Macedonian, Slovenian, Hungarian, Turkish, Korean, Vietnamese, Thai, Serbian, Malaysian, Danish, Bulgarian |
| Languages Full User Manual | German, English, Spanish, French, Italian, Portuguese, Chinese, Japanese |
| Dimensions | 422 x 195 x 323 mm |

Delivery Scope

- 4 batteries (AA)
- 1 power supply cable
- 4 round cuvettes with lid and 1 zero cuvette XD 7x00 (ø 24 mm)
- 1 zero cuvette (ø 16 mm)
- Quickstart guide in 24 languages
- Full user manual in 8 languages
- Warranty information

Accessories

| Title | Part Number |
|---|-------------|
| Adapter (13 mm) MultiDirect for Vacu-vial | 192075 |
| Batteries (AA), set of 4 | 1950025 |
| Round cuvette 24 mm, set of 12 | 197620 |
| Round cuvette 24 mm, set of 5 | 197629 |
| Cleaning cloth | 197635 |
| Round cuvette 16 mm, set of 10 | 197665 |
| Adapter for round cuvettes 13 mm | 19802192 |
| Mixing cylinder, 25 ml | 19802650 |
| Inspection and calibration package for XD7500 | 19802708 |
| Zero cuvette ø 16 mm for XD 7000/7500 | 215661 |
| Zero cuvette ø 24 mm for XD 7000/7500 | 215662 |
| Thermoreactor RD 125 | 2418940 |
| USB cable 3 m | 2444482 |
| Pipette, 1000 µl | 365045 |
| Measuring spoon, 1 g | 384930 |
| UV Pen Lamp, 254 nm | 400740 |
| UV protection glasses, orange | 400755 |
| Cuvette stand for 6 round cuvettes Ø 24 mm | 418951 |
| Cuvette stand for 10 round cuvettes Ø 16 mm | 418957 |
| Pipette tips, 1-5 ml (white) 100 pc. | 419066 |
| Pipette tips, 0,1-1 ml (white), 1000 pc. | 419073 |
| Automatic pipette, 1-5 ml | 419076 |
| Automatic pipette, 0,1-1 ml | 419077 |
| Screw caps TOC | 420757 |
| Measuring spoon no. 8, black | 424513 |
| Universal Container - Cap | 424648 |
| Plastic funnel with handle (white) | 471007 |
| ValidCheck Chlorine 1,5 mg/l | 48105510 |
| Stirring rod and spoon | 56A006601 |
| W100/OG/10MM Rectangular cell, optical glass | 601040 |
| W100/OG/20MM Rectangular cell, optical glass for determination of arsenic | 601050 |
| W100/OG/50MM Rectangular cell, optical glass | 601070 |
| W110/UV/10MM Rectangular cell, Quartz UV | 661130 |
| W110/UV/20MM Rectangular cell, Quartz UV | 661140 |
| W110/UV/50MM Rectangular cell, Quartz UV | 661160 |
| Energy station XD series/SpectroDirect | 711051 |
| Secondary standard set VIS with DAKkS calibration certificate | 711160 |
| 12 V-plug connector for XD 7000/7500 | 71310020 |
| Barcode Hand-held Scanner | 71310030 |
| Semimicro cell, 50 mm with lid | 71310045 |
| Factory calibration certificate ISO 9001 for XD7500 | 999755 |

Tintometer GmbH

Lovibond® Water Testing
Schleefstraße 8-12
44287 Dortmund
Tel.: +49 (0)231/94510-0
sales@lovibond.com
www.lovibond.com
Germany

The Tintometer Limited

Lovibond House
Sun Rise Way
Amesbury, SP4 7GR
Tel.: +44 (0)1980 664800
Fax: +44 (0)1980 625412
sales@lovibond.uk
www.lovibond.com
UK

Tintometer China

Room 1001, China Life Tower
16 Chaoyangmenwai Avenue,
Beijing, 100020
Tel.: +86 10 85251111 App. 330
Fax: +86 10 85251001
chinaoffice@tintometer.com
www.lovibond.com
China

Tintometer South East Asia

Unit B-3-12, BBT One Boulevard,
Lebuh Nilam 2, Bandar Bukit Tinggi,
Klang, 41200, Selangor D.E
Tel.: +60 (0)3 3325 2285/6
Fax: +60 (0)3 3325 2287
lovibond.asia@tintometer.com
www.lovibond.com
Malaysia

Tintometer Brazil

Caixa Postal: 271
CEP: 13201-970
Jundiaí – SP
Tel.: +55 (11) 3230-6410
sales@tintometer.com.br
www.lovibond.com.br
Brazil

Tintometer Inc.

6456 Parkland Drive
Sarasota, FL 34243
Tel: 941.756.6410
Fax: 941.727.9654
sales@lovibond.us
www.lovibond.com
USA

Tintometer India Pvt. Ltd.

Door No: 7-2-C-14, 2nd, 3rd & 4th Floor
Sanathnagar Industrial Estate,
Hyderabad, 500018
Telangana
Tel: +91 (0) 40 23883300
Toll Free: 1 800 599 3891/ 3892
indiaoffice@lovibond.in
www.lovibondwater.in
India

Tintometer Spain

Postbox: 24047
08080 Barcelona
Tel.: +34 661 606 770
sales@tintometer.es
www.lovibond.com
Spain

Technical changes without notice
Printed in Germany

Lovibond® and Tintometer® are Trademarks of the Tintometer Group of Companies