OPTIKA

M I C R O S C O P E S
I T A L Y

Ver. 4.1.0

 **OPERATION MANUAL**

 **GUIDA UTENTE**

 **MANUAL DE INSTRUCCIONES**

 **MANUEL D'INSTRUCTIONS**

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OPTIKA MICROSCOPES - ITALY

www.optikamicroscopes.com - info@optikamicroscopes.com

1.0 DESCRIPTION	page 3
2.0 INTRODUCTION	page 4
3.0 UNPACKING AND ASSEMBLY	page 4
4.0 USING THE MICROSCOPE	page 4
5.0 MAINTENANCE	page 5
6.0 RECOVERY AND RECYCLING	page 6

2.0 INTRODUCTION

This microscope is a scientific precision instrument designed to last for many years with a minimum of maintenance. It is built to high optical and mechanical standards and to withstand daily laboratory use.

Optika reminds you that this manual contains important information on safety and maintenance, and that it must therefore be made accessible to the instrument users.

Optika declines any responsibility deriving from instrument uses that do not comply with this manual.

3.0 UNPACKING AND ASSEMBLY

Unpack the microscope from the box and put it on a stable surface. Unscrew the bottom part of the penlight and insert two AAA batteries. Screw the lamp holder onto the microscope and fixate the lamp to the lamp holder with the screws. The microscope is now ready to use.

4.0 USING THE MICROSCOPE

- 4.1 Put the microscope on top of the object or surface you want to examine.
- 4.2 Switch on the light on the bottom of the lamp.
- 4.3 Rotate the focusing knob to bring the object into focus.
- 4.4 To avoid that the objective hits against the surface that are to be examined and in order to always use the same focusing distance, the minimum height of the objective can be set. When the image is in focus, lower the focus-stop ring on the microscope tube as much as is possible and fixate it using the screw.

5.1 Always think about

- The following environment is required: Indoor temperature: 0-40°C, Maximum relative humidity: 85 % (non condensing)
- Keep the microscope away from dust and shocks while in use.
- Turn off the light immediately after use.
- Use a soft lens tissue to clean the optics after use.
- Only if needed, use a cloth moistened with water and a mild detergent, rinsing with water and drying immediately with a lint-free cloth.
- After use, cover the microscope with the included dust-cover, and keep it in a dry and clean place.

5.2 Do not!

- Wipe the surface of any optical items with your hands. Fingerprints can damage the optics.
- Use solvents, neither on the microscope, nor on the optics.
- Disassemble objective or eyepieces to attempt to clean them.
- Mishandle or impose unnecessary force on the microscope.
- Clean the unit with volatile solvents or abrasive cleaners.
- Attempt to service the microscope yourself.

5.3 If you need to send the microscope to Optika for maintenance, please use the original packaging.

Art.13 Dlsg 25 july 2005 N°151. “According to directives 2002/95/EC, 2002/96/EC and 2003/108/EC relating to the reduction in the use of hazardous substances in electrical and electronic equipment and waste disposal.”

The basket symbol on equipment or on its box indicates that the product at the end of its useful life should be collected separately from other waste.

The separate collection of this equipment at the end of its lifetime is organized and managed by the producer. The user will have to contact the manufacturer and follow the rules that he adopted for end-of-life equipment collection. The collection of the equipment for recycling, treatment and environmentally compatible disposal, helps to prevent possible adverse effects on the environment and health and promotes reuse and/or recycling of materials of the equipment. Improper disposal of the product involves the application of administrative penalties as provided by the laws in force.

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