


# Instruction Manual

---

## HI8424

### Portable pH/mV/°C Meter



[www.hannainst.com](http://www.hannainst.com)

Dear Customer,

Thank you for choosing a Hanna Instruments Product.

Please read this instruction manual carefully before using this instrument. This manual will provide you with the necessary information for the correct use of this instrument, as well as a precise idea of its versatility.

If you need additional technical information, do not hesitate to e-mail us at [tech@hannainst.com](mailto:tech@hannainst.com) or view our worldwide contact list at [www.hannainst.com](http://www.hannainst.com).

## TABLE OF CONTENTS

PRELIMINARY EXAMINATION .....	3
GENERAL DESCRIPTION .....	3
FUNCTIONAL DESCRIPTION .....	4
SPECIFICATIONS .....	5
OPERATIONAL GUIDE .....	6
pH CALIBRATION .....	9
pH BUFFER TEMPERATURE DEPENDENCE .....	11
MENU SELECTION .....	12
mV CALIBRATION .....	13
TEMPERATURE CALIBRATION .....	13
BATTERY REPLACEMENT .....	14
LCD MESSAGES AND TROUBLESHOOTING .....	15
ELECTRODE CONDITIONING & MAINTENANCE .....	17
ACCESSORIES .....	19
pH ELECTRODE APPLICATION REFERENCE GUIDE .....	25

## PRELIMINARY EXAMINATION

Remove the instrument from the packing material and examine it carefully to make sure that no damage has occurred during shipment. If there is any damage, please contact your local Hanna Instruments Office.

The meter is supplied complete with:

- **HI1230B** pH electrode
- **HI7662** temperature probe
- pH 4.01 & pH 7.01 buffer solutions, 20 mL each
- **HI700661** cleaning solution, 20 mL sachet (2 pcs.)
- 9V battery
- Instruction manual

**Note:** Save all packing material until you are sure that the instrument functions correctly. Any defective items must be returned in the original packing with the supplied accessories.

## GENERAL DESCRIPTION

**HI8424** is a portable microprocessor-based pH/mV/temperature meter.

It features an enhanced user interface, rainproof casing, battery percentage indication, low battery detection, automatic shut-off, automatic calibration and error codes to guide the user in calibration and troubleshooting.


The pH calibration procedure automatically recognizes three memorized buffer values: pH 4.01, 7.01 and 10.01.

This instrument can also measure ORP. The resolution automatically switches from 0.1 mV to 1 mV when readings exceed  $\pm 699.9$  mV.

The user can set the following parameters by entering the setup mode through the keyboard:

- enabling or disabling the auto-off feature
- enabling or disabling the acoustic signal
- selecting the temperature unit, °C or °F

## FUNCTIONAL DESCRIPTION



- 1) BNC connector for pH or ORP electrodes
- 2) Liquid Crystal Display (LCD)
- 3) **RANGE** key, to select pH, mV or temperature range
- 4) **Up Arrow** key, to manually set the temperature value when no temperature probe is connected
- 5) **CAL/SELECT** key, to enter calibration mode or select menu options
- 6) **CFM/SETUP** key, to confirm data or enter/exit setup menu
- 7) **Down Arrow** key, to manually set the temperature value when no temperature probe is connected
- 8) **HOLD** key, to freeze the reading on display
- 9) **ON/OFF** key, to switch the instrument ON or OFF
- 10) RCA socket for temperature probe

## SPECIFICATIONS

Range	-2.00 to 16.00 pH ±699.9 mV/ ±1999 mV -20.0 to 120.0 °C/ -4.0 to 248.0 °F
Resolution	0.01 pH/ 0.1 mV/ 1 mV/ 0.1 °C/ 0.1 °F
Accuracy (@20 °C/68 °F)	±0.01 pH/ ±0.2 mV/ ±1 mV/ ±0.4 °C/ ±0.8 °F
Typical EMC Deviation	±0.02 pH/ ±0.2 mV/ ±1 mV/ ±0.4 °C/ ±0.8 °F
pH Calibration	Automatic, 1 or 2 point, with 3 memorized buffer values (pH 4.01, 7.01, 10.01) Offset: ±1 pH; Slope: from 75 to 110%
Temperature Compensation	Automatic, -20 to 120 °C (-4 to 248 °F) or manual without temperature probe
Probes (included)	H11230B double junction, gel-filled pH electrode H17662 temperature probe
Battery Type	9V (1 pc.)
Battery Life	Approx. 150 hours of continuous use
Auto-off	After 20 minutes of non-use or disabled (user-selectable)
Environment	0 to 50 °C (32 to 122 °F); RH max 100%
Dimensions	164 x 76 x 45 mm (6.5 x 3.0 x 1.8")
Weight	180 g (6.3 oz.)

# OPERATIONAL GUIDE


## INITIAL PREPARATION

- Remove the electrode protective cap before taking any measurements. If the electrode has been left dry, soak the tip in **H170300** storage solution for a few hours or overnight to reactivate it.
- Connect the pH electrode to the BNC connector on the top of the instrument.
- Connect the temperature probe to the RCA connector. The temperature probe can be used independently to take temperature measurements, or in conjunction with the pH electrode to utilize the ATC capability of the meter.
- Turn the meter ON by pressing the ON/OFF key. The display shows all the used segments for a few seconds (or as long as the button is held), followed by the percentage indication of the remaining battery life, and then enters normal measurement mode.



## pH MEASUREMENTS

- To take a pH measurement simply submerge the electrode tip (at least 4 cm/1½") and the temperature probe into the sample to be tested.
- Select the pH mode by pressing the RANGE key until the display changes to pH.
- Stir gently and wait for the stability symbol (hourglass) to turn off. The display will show the pH value automatically compensated for temperature.



### Notes:

- In order to take accurate pH measurements, make sure that the instrument has been calibrated before use (see page 9).
- If measurements are taken in different samples successively, it is recommended to rinse the electrode thoroughly to avoid cross-contamination. After cleaning, it is recommended to rinse the electrode with some of the sample to be measured.

## TEMPERATURE COMPENSATION

The meter is designed to compensate for temperature, as the response of the pH electrode is directly affected by temperature.


### Automatic Temperature Compensation (ATC shown on LCD)

To use the ATC feature, submerge the temperature probe into the sample as close as possible to the electrode and wait for a few minutes. The displayed pH reading is compensated for the temperature of the sample.

### Manual Temperature Compensation (MTC shown on LCD)

If the temperature probe is not connected, it is possible to enter the temperature value manually.

- Record the sample temperature by using a **ChecktempC** (if you are measuring temperature in °C, or **ChecktempF** for °F readings) or another accurate thermometer.
- Press RANGE to select the temperature mode. The “°C” (or “°F”) symbol will blink to indicate that the temperature probe is not connected.




- Use the UP and DOWN keys to display and set the sample temperature (e.g. 25 °C).
- Press RANGE to select the pH measurement mode and immerse the electrode into the sample. The displayed pH reading will be temperature compensated at the set value (in this case at 25 °C).



## ORP MEASUREMENTS

Oxidation Reduction Potential (ORP) measurements provide a quantification of the oxidizing or reducing power of the sample tested.

- Connect the ORP electrode (optional) to the BNC connector.
- To enter the “mV” mode turn the instrument ON and press the RANGE key until the display changes to mV.
- Submerge the ORP electrode tip (at least 4 cm / 1½”) into the sample to be tested and allow time for the reading to stabilize (hourglass symbol turns off).
- Measurements within the  $\pm 699.9$  mV range are displayed with 0.1 mV resolution, while outside this range the resolution is 1 mV.



### Notes:


- To perform correct ORP measurements, the surface of the ORP electrode must be clean and smooth.
- When not in use, the tip of the electrode should be kept moist (use **HI70300** storage solution) and safe from any mechanical stress which might cause damage to the glass/platinum junction.

## TEMPERATURE MEASUREMENTS

- Turn the instrument ON and press the RANGE key to select the temperature mode.
- Make sure the temperature probe is connected to the meter.
- Dip the temperature probe into the sample, allow the reading to stabilize (hourglass symbol turns off) and read the temperature value.
- Temperature measurements can be displayed in °C or °F units (see “Menu selection” for details).

### Notes:

- A blinking full scale value means that the reading is out of range.
- To freeze a reading on display while in measurement mode, press the HOLD key. The “HOLD” tag will blink. The pH, mV and temperature values are held, and the RANGE key can be used to view the values. Press HOLD again to return to normal mode.
- If enabled, keypresses are followed with an acoustic signal. A lower note indicates that the key is not currently active.
- To save battery life, the meter is provided with an **auto-off feature**, which turns the instrument off after 20 minutes of non-use. This feature can be disabled by the user (see “Menu selection” for details).





## pH CALIBRATION

For better accuracy, frequent calibration of the instrument is recommended. The instrument should be recalibrated for pH:

- Whenever the pH electrode or temperature probe is replaced
- At least once a week
- After testing aggressive chemicals
- When extreme accuracy is required

### PREPARATION

Pour small quantities of pH 7.01 (HI7007) and pH 4.01 (HI7004) or pH 10.01 (HI7010) buffer solutions into two clean beakers.

For accurate calibration use two beakers for each buffer solution, one for rinsing the electrode tip, and one for calibration. In this way contamination of the buffers is minimized.


For measurements in acidic samples, it is recommended to calibrate the meter by using pH 7.01 (HI7007) and pH 4.01 (HI7004) buffers, while for alkaline measurements use pH 7.01 (HI7007) and pH 10.01 (HI7010) buffers.

### PROCEDURE



- Connect the pH electrode and the temperature probe, then switch the meter ON.
- Remove the electrode protective cap, rinse the electrode tip with pH 7.01 solution, then immerse the pH electrode and temperature probe into pH 7.01 buffer solution; stir gently and wait a few minutes for the electrode to stabilize and reach thermal equilibrium.

**Note:** The electrode should be submerged approximately 4 cm (1½") into the solution. The temperature probe should be located close to the pH electrode.

- Press RANGE to display pH measurement.
- Press CAL to enter the calibration mode. The buffer value @25 °C (77 °F) and the "pH" symbol will blink on the display.
- The meter expects a pH 7.01 buffer.



- When the buffer value is recognized and the reading is stable, an acoustic signal (if enabled) advises the user, the “pH” symbol stops flashing, the hourglass indicator disappears and the “CFM” tag starts blinking to indicate that the value can be confirmed.
- Press CFM to store the first calibration point.
- The meter expects a pH 4.01 or 10.01 buffer.
- Rinse and immerse the pH electrode and the temperature probe in pH 4.01 or pH 10.01 buffer (2<sup>nd</sup> calibration point) and stir gently.
- When the buffer value is recognized and the reading is stable, an acoustic signal (if enabled) advises the user, the “pH” symbol stops flashing, the hourglass indicator disappears and the “CFM” tag starts blinking to indicate that the value can be confirmed.
- Press CFM to store the second calibration point.
- The meter returns to normal mode.




The pH calibration is now complete; “CAL” and the pH tags corresponding to the buffers used for calibration are lit on the LCD.



### Notes:

- If the buffer value is not recognized, after 12 seconds the meter will display blinking dashes together with the “WRONG” tag. Either the buffer solution is wrong or out of specification and needs to be replaced or the electrode is damaged.
- The meter will retain the calibration if the battery is removed.
- To quit calibration and keep previous data: press CAL after entering the calibration mode and before the first point is accepted.
- To perform a single-point calibration: press CAL after the first point has been confirmed.



If the temperature probe is not connected and manual temperature compensation is required, follow the procedure below:

- Press RANGE to select the temperature mode.
- Rinse the pH electrode and place it into the pH 7.01 buffer, stir briefly and wait a few minutes to reach thermal equilibrium.

- Rinse the temperature probe of a **ChecktempC** (or **ChecktempF**) or another accurate thermometer, and place it close to the pH electrode.
- Use the UP and DOWN arrow keys to manually adjust the temperature to match the reference thermometer.



- Follow the pH calibration procedure explained in the previous pages.

## pH BUFFER TEMPERATURE DEPENDENCE

The temperature has an effect on pH. The calibration buffer solutions are affected by temperature changes to a lesser degree than normal solutions. During calibration the instrument will automatically calibrate to the pH value corresponding to the measured or set temperature.

During calibration the instrument will display the pH buffer value at 25 °C.

TEMP		pH VALUES		
°C	°F	4.01	7.01	10.01
0	32	4.01	7.13	10.32
5	41	4.00	7.10	10.24
10	50	4.00	7.07	10.18
15	59	4.00	7.04	10.12
20	68	4.00	7.03	10.06
25	77	4.01	7.01	10.01
30	86	4.02	7.00	9.96
35	95	4.03	6.99	9.92
40	104	4.04	6.98	9.88
45	113	4.05	6.98	9.85
50	122	4.06	6.98	9.82
55	131	4.07	6.98	9.79
60	140	4.09	6.98	9.77
65	149	4.11	6.99	9.76
70	158	4.12	6.99	9.75
75	167	4.14	7.00	9.74
80	176	4.16	7.01	9.73
85	185	4.17	7.02	9.74
90	194	4.19	7.03	9.75
95	203	4.20	7.04	9.76

## MENU SELECTION

While in normal measurement mode, press and hold the CFM/SETUP key for about 5 seconds until the meter enters the menu selection mode.

The following parameters can be set from the menu:


1. Auto-off feature: 20 minutes (default setting) or disabled;
2. Acoustic signal: enabled (default setting) or disabled;
3. Temperature unit: °C (default setting) or °F.

When entering the menu mode, the auto-off selection is entered. The LCD shows in three subsequent screens "Auto", "OFF" and "20" to indicate that the 20 minutes selection is active, or "Auto", "OFF" and "no" if the feature is disabled.




The user can toggle the selection by pressing CAL/SELECT or move to the next step with CFM/SETUP.

The following selection is the acoustic signal, which is displayed on two subsequent screens: "bEEP", "OFF" when the feature is disabled, and "bEEP", "On" when the feature is enabled.



Press CAL/SELECT to toggle the selection and CFM/SETUP to move to next step.

At this point it is possible to set the temperature unit, by selecting "Unit", "°C" or "Unit", "°F".



Press CAL/SELECT to toggle the selection and CFM/SETUP to exit the menu selection mode and return to normal measurement mode.

## **mV CALIBRATION**

**HI8424** has been accurately precalibrated for mV range at the factory.

For optimum accuracy, it is recommended to recalibrate the meter for mV readings at least once a year.

Contact your local Hanna Instruments Office for more information.

## **TEMPERATURE CALIBRATION**


**HI8424** has been accurately precalibrated for temperature at the factory.

For optimum accuracy, it is recommended to recalibrate the meter for temperature at least once a year.

Contact your local Hanna Instruments Office for more information.


## BATTERY REPLACEMENT

The meter displays the remaining battery percentage when turned on. When the level is below 5%, the battery symbol on the bottom left of the LCD blinks to indicate a low battery condition.





If the battery level is low enough to cause erroneous readings, the Battery Error Prevention System (BEPS) turns the meter off. When the battery needs to be replaced, remove the cover on the rear of the meter and replace the rundown battery with a new one, while paying attention to the correct polarity. Reattach the back making sure that the gasket is in place and tighten the 3 screws to ensure a good seal.

Replacement should take place in a non-hazardous area using a 9V alkaline battery.



## LCD MESSAGES & TROUBLESHOOTING

### TAGS & SYMBOLS

- pH, mV, °C, °F Measurement unit of the selected mode
- ATC Indicates Automatic Temperature Compensation (in pH or temperature mode)
- MTC Indicates Manual Temperature Compensation (in pH or temperature mode)
- HOLD Blinks when in Hold mode. Reading frozen on LCD. The user can scroll through the three ranges by pressing RANGE
- CAL In pH calibration mode, or in pH mode when the meter is calibrated
- CFM Blinks in pH calibration mode when the meter is ready to confirm a value
- WRONG During pH calibration, when the meter does not recognize the pH buffer
- % At startup, when showing the percentage of the remaining battery life
- **7** In pH mode, when meter was calibrated with pH 7.01 buffer
- **4** In pH mode, when meter was calibrated with pH 4.01 buffer
- **10** In pH mode, when meter was calibrated with pH 10.01 buffer
-   
(hourglass symbol) When reading is not stable
-   
(battery symbol) At startup, if remaining battery life is below 5%


## TROUBLESHOOTING

Symptom	Problem	Solution
Slow response or excessive drift	Dirty pH electrode	Soak the electrode tip in HI7061 solution for 30 minutes
Reading fluctuates up and down (noise)	Clogged/dirty junction or low electrolyte level (refillable electrodes)	Soak the electrode tip in warm HI7082 solution for one hour, then rinse it with distilled water (refill with fresh electrolyte if necessary)
Blinking full scale value	Reading is out of range	
Blinking "°C" (or "°F")	Temperature probe is not connected or broken	
"WRONG" & blinking dashes	Calibration error	Check buffer solution or replace pH electrode
Blinking battery symbol	Low battery level	Replace battery
Meter shuts off	Auto-off enabled or dead battery	Replace battery
"Clr" message	Loaded default pH calibration values	Perform pH calibration
"Er1" and "Er2" messages	EPROM error	Contact your local Hanna Instruments Office

**Note:** For field applications, it is always recommended to keep a conditioned spare electrode handy. When anomalies cannot be resolved with simple maintenance, change the electrode and recalibrate the meter.



# ELECTRODE CONDITIONING & MAINTENANCE



\* Only for refillable electrodes; must be open while measuring.

## PREPARATION PROCEDURE

Remove the electrode protective cap.

DO NOT BE ALARMED IF ANY SALT DEPOSITS ARE PRESENT. This is normal with electrodes and they will disappear when rinsed with water.

During transport tiny bubbles of air may have formed inside the glass bulb. The electrode cannot function properly under these conditions. These bubbles can be removed by "shaking down" the electrode as you would do with a glass thermometer.

If the bulb and/or junction are dry, soak the electrode in **HI70300** storage solution for at least one hour.

For refillable electrodes, if the refill solution (electrolyte) is more than 2.5 cm (1") below the fill hole, add the appropriate electrolyte solution.

## MEASUREMENT

Rinse the electrode tip with distilled water, immerse it (4 cm / 1½") in the sample and stir gently for a few seconds.

For a faster response and to avoid cross contamination of the samples, rinse the electrode tip with the solution to be tested, before taking any measurements.

## STORAGE PROCEDURE

To minimize clogging and ensure a quick response time, the glass bulb and the junction should always be kept moist.

When not in use, store it with a few drops of **HI70300** storage solution in the protective cap.

**NEVER STORE THE ELECTRODE IN DISTILLED OR DEIONIZED WATER.**

## PERIODIC MAINTENANCE

Inspect electrode and cable. The cable used for the connection to the meter must be intact and there must be no points of broken insulation on the cable or cracks on the electrode stem or bulb. If any scratches or cracks are present, replace the electrode. Rinse off any salt deposits with water.

Connectors must be perfectly clean and dry.

### For refillable electrodes:

Refill the electrode with fresh electrolyte (see the electrode's specifications to select the correct refilling solution). Allow the electrode to stand upright for 1 hour. Follow the Storage Procedure above.

## CLEANING PROCEDURE

- *General* Soak in **HI7061** general cleaning solution for approximately 30 minutes.
- *Protein* Soak in **HI7073** protein cleaning solution for 15 min.
- *Inorganic* Soak in **HI7074** inorganic cleaning solution for 15 minutes.
- *Oil/grease* Rinse with **HI7077** Oil & Fat cleaning solution for 1 minute.

**IMPORTANT:** After performing any of the cleaning procedures, rinse the electrode thoroughly with distilled water and soak it in **HI70300** storage solution for at least 1 hour before taking measurements.

## ACCESSORIES

### pH CALIBRATION SOLUTIONS

HI70004P	pH 4.01 buffer solution, 20 mL sachet (25 pcs)
HI7004M	pH 4.01 buffer solution, 230 mL bottle
HI7004L	pH 4.01 buffer solution, 500 mL bottle
HI8004L	pH 4.01 buffer solution, 500 mL FDA bottle
HI70007P	pH 7.01 buffer solution, 20 mL sachet (25 pcs)
HI7007M	pH 7.01 buffer solution, 230 mL bottle
HI7007L	pH 7.01 buffer solution, 500 mL bottle
HI8007L	pH 7.01 buffer solution, 500 mL FDA bottle
HI70010P	pH 10.01 buffer solution, 20 mL sachet (25 pcs)
HI7010M	pH 10.01 buffer solution, 230 mL bottle
HI7010L	pH 10.01 buffer solution, 500 mL bottle
HI8010L	pH 10.01 buffer solution, 500 mL FDA bottle

### STORAGE & CLEANING SOLUTIONS

HI70300M	Storage solution, 230 mL bottle
HI80300M	Storage solution, 230 mL FDA bottle
HI70300L	Storage solution, 500 mL bottle
HI80300L	Storage solution, 500 mL FDA bottle
HI70000P	Electrode rinsing solution, 20 mL sachet (25 pcs.)
HI7061M	General cleaning solution, 230 mL bottle
HI8061M	General cleaning solution, 230 mL FDA bottle
HI7061L	General cleaning solution, 500 mL bottle
HI8061L	General cleaning solution, 500 mL FDA bottle
HI7073M	Protein cleaning solution, 230 mL bottle
HI8073M	Protein cleaning solution, 230 mL FDA bottle
HI7073L	Protein cleaning solution, 500 mL bottle
HI8073L	Protein cleaning solution, 230 mL FDA bottle
HI7074M	Inorganic cleaning solution, 230 mL bottle
HI7074L	Inorganic cleaning solution, 500 mL bottle
HI7077M	Oil & Fat cleaning solution, 230 mL bottle
HI8077M	Oil & Fat cleaning solution, 230 mL FDA bottle
HI7077L	Oil & Fat cleaning solution, 500 mL bottle
HI8077L	Oil & Fat cleaning solution, 500 mL FDA bottle

### REFILLING ELECTROLYTE SOLUTIONS

HI7071	3.5 M KCl + AgCl electrolyte solution, 4 x 30 mL bottle, for single junction electrodes
HI8071	3.5 M KCl + AgCl electrolyte solution, 4 x 30 mL FDA bottle, for single junction electrodes
HI7072	1 M KNO <sub>3</sub> electrolyte solution, 4 x 30 mL bottle
HI8072	1M KNO <sub>3</sub> Electrolyte, 4x30 mL (FDA approved bottle)


- HI7082 3.5 M KCl electrolyte solution, 4 x 30 mL bottle, for double junction electrodes
- HI8082 3.5 M KCl electrolyte solution, 4 x 30 mL FDA bottle, for double junction electrodes
- HI8093 1 M KCl+AgCl electrolyte solution, 4 x 30 mL FDA bottle, for double junction electrodes

## ORP SOLUTIONS


- HI7091L Reducing pretreatment solution, 500 mL bottle
- HI7092M Oxidizing pretreatment solution, 230 mL bottle
- HI7092L Oxidizing pretreatment solution, 500 mL bottle
- HI7021M Test solution @240 mV, 230 mL bottle
- HI7021L Test solution @240 mV, 500 mL bottle
- HI7022M Test solution @470 mV, 230 mL bottle
- HI7022L Test solution @470 mV, 500 mL bottle

## pH ELECTRODES

B = BNC CONNECTION PLUG + 1 m (3.3') CABLE




S = SCREW-TYPE CONNECTOR



### HI1043B / HI1040S

Glass body, double junction, refillable, combination pH electrode.


Use: strong acid/alkali.



### HI1053B / HI1050S

Glass body, triple ceramic, conic shape, refillable, combination pH electrode.


Use: emulsions.



### HI1083B

Glass body, micro, Viscolene, non-refillable, combination pH electrode.


Use: biotechnology, micro titration.



### HI1131B / HI1111S

Glass body, single junction, refillable, combination pH electrode.


Use: general purpose.



### HI1230B / HI1210S

Plastic body, double junction, gel-filled, combination pH electrode.


Use: general purpose.



### HI1330B / HI1310S

Glass body, semimicro, single junction, refillable, combination pH electrode.


Use: laboratory.



### HI1331B / HI1311S

Glass body, semimicro, single junction, refillable, combination pH electrode.


Use: flasks.



### HI1332B / HI1312S

Plastic body, double junction, refillable, combination pH electrode.


Use: general purpose.



### HI2031B / HI2020S

Glass body, semimicro, conic, refillable, combination pH electrode.


Use: semisolid products.



### FC100B


PVDF body, double junction, refillable, combination pH electrode.

Use: general purpose for food industry.




### FC200B / FC200S

PVDF body, single junction, conic, Viscolene, non-refillable, combination pH electrode. Use: meat & cheese.




### FC210B

Glass body, double junction, conic, Viscolene, non-refillable combination pH electrode. Use: milk, yogurt.




### FC220B

Glass body, triple ceramic, single junction, refillable, combination pH electrode.  
Use: food processing.




### FC911B

PVDF body, double junction, refillable with built-in amplifier, combination pH electrode. Use: very high humidity.



### HI1413B


Glass body, single junction, flat tip, Viscolene, non-refillable combination pH electrode. Use: surface measurement.



## ORP ELECTRODES


### HI3131B / HI3111S

Glass body, refillable, combination platinum ORP electrode.  
Use: titration.



### HI3230B / HI3210S


Plastic body, gel-filled, combination platinum ORP electrode.  
Use: general purpose.



## HI4430B / HI4410S

Plastic body, gel-filled, combination gold ORP electrode.

Use: general purpose.




Consult the HANNA Instruments General Catalog for a complete and wide selection of electrodes.

## EXTENSION CABLES FOR SCREW-TYPE ELECTRODES (SCREW TO BNC CONNECTOR)

HI 7855 SERIES CABLE CONNECTORS  
CONNECTOR AND 3.0 mm (0.12") CABLE WITH BNC

CONNECT TO  
SCREW TYPE  
ELECTRODES



CONNECT TO THE  
BNC SOCKET  
OF THE METER

- HI7855/1 Extension cable 1 m (3.3') long
- HI7855/3 Extension cable 3 m (9.9') long
- HI7855/5 Extension cable 5 m (16.5') long
- HI7855/10 Extension cable 10 m (33') long
- HI7855/15 Extension cable 15 m (49.5') long

## OTHER ACCESSORIES

- HI98501 **ChecktempC** pocket-size thermometer (-50.0 to 150.0 °C)
- HI98502 **ChecktempF** pocket-size thermometer (-58.0 to 302.0 °F)
- HI710015 Shockproof rubber boot, blue
- HI710016 Shockproof rubber boot, orange
- HI710022 Spare protective case
- HI76405 Electrode holder
- HI7662 Temperature probe with 1 m (3.3') screened cable
- HI8427 pH/ORP electrode simulator with 1 m (3.3') coaxial cable and BNC connector
- HI931001 pH/ORP electrode simulator with LCD, 1 m (3.3') coaxial cable and BNC connector



## pH ELECTRODE APPLICATION REFERENCE GUIDE

Application	Electrodes *
1. Aquarium	HI1332B, HI1312S
2. Bath-water	HI1130B, HI1110S
3. Beer	HI1131B, HI1111S
4. Bread	HI2031B, FC200B, HI2020S, FC200S
5. Cheese	FC 200B, FC 200S
6. Dairy products	FC 911B, FC100B
7. Dirty water	HI1230B, HI1210S
8. Emulsions	HI1053B, HI1050S
9. Environment	HI1230B, HI1210S
10. Flasks	HI1331B, HI1310S
11. Food industry general use	FC 911B, FC100B
12. Fruit	FC200B, FC 220B, FC200S
13. Fruit juices, organic	FC210B
14. Galvanizing waste solution	HI1130B, HI1110S
15. High purity water	HI1053B, HI1050S
16. Horticulture	HI1053B, FC200B, HI1050S, FC 200S
17. Laboratory general use	HI1131B, HI1230B, HI1332B, HI1330B HI1111S, HI1210S, HI1312S, HI1310S
18. Leather	HI1413B, HI1410S
19. Lemon juice	FC100B
20. Meat	FC200B, HI2031B, FC200S, HI2020S
21. Micro plate sampling of less than 100 ml	HI1083B
22. Milk and Yogurt	FC 210B
23. Paints	HI1053B, HI1050S
24. Paper	HI1413B, HI1410S
25. Photographic chemicals	HI1230B, HI1210S
26. Quality control	HI1332B, HI1312S
27. Sausages	FC 200B, HI2031B, FC 200S, HI2020S
28. Semi-solid products	HI2031B, HI2020S
29. Skin	HI1413B, HI1410S
30. Soil samples	HI1230B, HI1210S
31. Solvents	HI1043B, HI1040S
32. Strong acid	HI1043B, HI1040S
33. Submersion application	HI1130B, HI1110S
34. Surface measurements	HI1413B, HI1410S
35. Swimming pool	HI1130B, HI2114P/2
36. Titrations with constant temperature range	HI1131B, HI1111S
37. Titrations with wide temperature range	HI1131B, HI1111S
38. Very high humidity	FC 911B
39. Vials and test tube	HI1330B, HI1310S
40. Wine processing	FC220B

B = BNC-type connector

S = Screw-type connector

\* All electrodes ending with "B" are supplied with 1 m (3.3') cable and BNC Connector

## **RECOMMENDATIONS FOR USERS**

Before using this product, make sure that it is entirely suitable for the environment in which they are used.

Operation of this instrument in residential area could cause unacceptable interferences to radio and TV equipments, requiring the operator to take all necessary steps to correct interferences.

The metal band at the end of the sensor is sensitive to electrostatic discharges. Avoid touching this metal band at all times. During calibration, ESD wrist straps should be worn to avoid possible damage to the sensor by electrostatic discharge.

Any variation introduced by the user to the supplied equipment may degrade the instrument's EMC performance. To avoid electrical shock, do not use this instrument when voltages at the measurement surface exceed 24 Vac or 60 Vdc. Use plastic beakers to minimize any EMC interferences.

To avoid damages or burns, do not perform any measurement in microwave ovens.

Hanna Instruments reserves the right to modify the design, construction or appearance of its products without advance notice.

## **WARRANTY**

**HI8424** is guaranteed for two years against defects in workmanship and materials when used for their intended purpose and maintained according to instructions. Electrodes and probes are guaranteed for a period of six months. This warranty is limited to repair or replacement free of charge.

Damages due to accident, misuse, tampering or lack of prescribed maintenance are not covered.

If service is required, contact your local Hanna Instruments Office. If under warranty, report the model number, date of purchase, serial number and the nature of the problem. If the repair is not covered by the warranty, you will be notified of the charges incurred. If the instrument is to be returned to Hanna Instruments, first obtain a Returned Goods Authorization number from the Technical Service department and then send it with shipping costs prepaid.

When shipping any instrument, make sure it is properly packed for complete protection.

*All rights are reserved. Reproduction in whole or in part is prohibited without the written consent of the copyright owner.*



**Hanna Instruments Inc.**

Highland Industrial Park

584 Park East Drive

Woonsocket, RI 02895 USA

**Technical Support for Customers**

Tel. (800) 426 6287

Fax (401) 765 7575

E-mail [tech@hannainst.com](mailto:tech@hannainst.com)

[www.hannainst.com](http://www.hannainst.com)