# Lovibond® Colour Measurement

# Tintometer® Group



## **RCMSI**

Ensures Correct Operation and Compliance with International Standards



- Reduced Cost & Carbon Footprint
- Guarantees Instrument is Fully Calibrated to ISO 17025 Certified Standard
- Results Compared to Master Data at UKAS Accredited, ISO 9001 : 2008 Qualified Centre
- Certificate of Calibration Provided on Successful Completion

Part Number: RCMSI

# Remote Calibration & Maintenance Service via internet (**RCMSi**)

RCMSi is the process by which a number of tightly controlled calibration liquids are measured in a 50 mm cell\* via the calibration function of a PFXi. The data is communicated via the internet to a secure server and compared to Master Data. If the data passes within tolerance, then a Certificate of Calibration is made electronically available to the user. This allows The Tintometer Ltd to ensure a PFXi instrument is operating correctly and fully calibrated using ISO 17025 certified standards.

## On-Site Calibration - No Need To Return Instrument to Service Centre

Colour measurement instruments need regular checking with calibration standards to ensure correct operation and compliance to established international standards. Being able to do this on-site, saves costs, carbon footprint and means you are never without your instrument.

### **Utilizes Secure Communication Procedures**

The data is communicated via the internet to a secure server and compared to Master Data at a UKAS accredited, ISO 9001: 2008 qualified centre. No other data is sent, only calibration data. This allows Lovibond Tintometer to ensure a PFXi instrument is operating correctly and is traceable to ISO17025 standards.

## Certificate of Calibration Supplied on Successful Completion

If the data passes within tolerance, then a Certificate of Calibration can be downloaded by logging into the Lovibond® website, giving you written confirmation that your instrument is performing correctly,

has been calibrated and is compliant to international standards - everything you need for internal and external audits.

### Take Control of Your Processes

You can determine how often you carry out the procedure in accordance with your Good Laboratory Practice (GLP). There is no limit as to how often you can run RCMSi.

## RCMSi is Not Mandatory

It is simply an additional, purely optional, level of performance guarantee. There are a number of PFXi users worldwide who, for security or infrastructure reasons, do not use RCMSi. This does not have any negative effect on their day-to-day operations. You can simply use the Glass Filter that ships as standard with the instrument or purchase additional Liquid Reference Standards or Glass Filters.

\*(Please note 50mm cell must be used when performing RCMSi calibration.)

### **RCMSI**

Colour measurement instruments need regular checking with calibration standards to ensure correct operation and compliance to established international standards. Being able to do this on-site, saves costs and carbon footprint and means you are never without your instrument.

## **Delivery Scope**

The Following Components Are Required To Carry Out RCMSi

- A Lovibond® PFXi instrument with LAN connection
- An internet connection
- A valid email address to register and access the servers & calibration data at the Lovibond<sup>®</sup> secure calibration server in Amesbury
- A valid Calibration Pack

RCMSi Starter Pack

- 2 x PFXi CalSol A (non haz)
- 2 x PFXi CalSol B (non haz)
- 2 x PFXi CalSol C (non haz)
- 2 x PFXi CalSol D (non haz)
- 2 x PFXi CalSol E (non haz)
- 2 x PFXi CalSol F (non haz)
- 2 x PFXi CalSol G (non haz)
- 1 x RCMSi MSDS Sheet

RCMSi Four Pack

- 2 x PFXi CalSol A (non haz)
- 2 x PFXi CalSol B (non haz)
- 2 x PFXi CalSol C (non haz)
- 2 x PFXi CalSol D (non haz)

### Tintometer GmbH

Lovibond® Water Testing Schleefstraße 8-12 44287 Dortmund Tel.: +49 (0)231/94510-0 sales@lovibond.com www.lovibond.com Germany

### **Tintometer Brazil**

Brazil

Caixa Postal: 271 CEP: 13201-970 Jundiaí – SP Tel.: +55 (11) 3230-6410 sales@tintometer.com.br www.lovibond.com.br

### The Tintometer Limited

Lovibond House Sun Rise Way Amesbury, SP4 7GR Tel.: +44 (0)1980 664800 Fax: +44 (0)1980 625412 sales@lovibond.uk www.lovibond.com UK

## Tintometer Inc.

6456 Parkland Drive Sarasota, FL 34243 Tel: 941.756.6410 Fax: 941.727.9654 sales@lovibond.us www.lovibond.com USA

### **Tintometer China**

Room 1001, China Life Tower 16 Chaoyangmenwai Avenue, Beijing, 100020 Tel.: +86 10 85251111 App. 330 Fax: +86 10 85251001 chinaoffice@tintometer.com www.lovibond.com China

## Tintometer India Pvt. Ltd.

Door No: 7-2-C-14, 2nd, 3rd & 4th Floor Sanathnagar Industrial Estate, Hyderabad, 500018 Telangana Tel: +91 (0) 40 23883300 Toll Free: 1 800 599 3891/3892 indiaoffice@lovibond.in www.lovibondwater.in India

#### **Tintometer South East Asia** Unit B-3-12, BBT One Boulevard,

Unit B-3-12, BB1 One Bottlevard, Lebuh Nilam 2, Bandar Bukit Tinggi, Klang, 41200, Selangor D.E Tel.: +60 (0)3 3325 2285/6 Fax: +60 (0)3 3325 2287 lovibond.asia@tintometer.com www.lovibond.com Malaysia

### **Tintometer Spain**

Postbox: 24047 08080 Barcelona Tel.: +34 661 606 770 sales@tintometer.es www.lovibond.com Spain

Technical changes without notice Printed in Germany Lovibond® and Tintometer® are Trademarks of the Tintometer Group of Companies