# Elettrodi pH e redox

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# pH and Redox Electrodes

## **Commissioning the electrode**

Remove the protective cap and rinse the electrode tip with some water. Do not touch the sensitive membrane by hands.

## **Calibration**

Always use fresh and clean standard solutions. Follow the instructions given in the manual of the instrument.

**Caution!** Do not rub the probe while drying it, to avoid wrong readings due to the static electricity.

## **Cleaning and storage**

Typically, it is recommended to clean the electrodes when the response is slow or measurements are not reliable, and when they have been used for a long time, especially in aggressive solutions, pollutants, very acidic or very alkaline environments.

To have your electrode at top performance, the sensing bulb must always be kept moist. Never store the electrode in distilled water, but use the appropriate solution.

A kit of solutions is available, for the periodic cleaning and storage of pH & redox electrodes.

## **Electrode life**

Measurement probes are subject to a natural ageing even if they are correctly handled. Please notice that the average life of an electrode varies from 6 to 18 months, depending on the application.

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