

SCI-VS/SCI-FS Vortex Mixers

User Manual



SCI-VS

Adjustable Speed Vortex

SCI-FS

Fixed Speed Vortex

Please read the User Manual and the related Video of this instrument on our website carefully before use, and follow all operating and safety instructions!

The website www.scilogex.com

SCILOGEX

Contents

Contents	1
Preface	2
Service	2
Warranty	2
1 Safety Instructions	3
2 Proper use	4
3 Inspection	4
3.1 Receiving Inspection	4
3.2 Listing of Items	4
4 Trial run	4
5 Operation	5
5.1 Continuous Operating	5
5.2 Touch Operating	5
6 Maintenance and Cleaning	6
7 Associated Standards and Regulations	6
8 Technical Data	7
9 Products and Accessories	7

Preface

Welcome to the “Vortex Mixers (Adjustable and Fixed Speed Models) User Manual” . Users should read this manual carefully, follow the instructions and procedures, and beware of all the cautions.

Service

When help needed, you can always contact the service department of manufacturer for technical support in the following ways:

SCIOGEX, LLC
1275 Cromwell Ave
Suite C6
Rocky Hill, CT 06067
USA
Tel: 1- (860) 828-5614
Fax: 1- (860) 828-5389

info@scilogex.com
www.scilogex.com

- Description of problem
- Methods and procedures adopted once to resolve the problems
- Your contact information.

Warranty

You have purchased a Scilogex instrument. This instrument is warranted to be free from defects in materials and workmanship under normal use and service, for a period of 24 months from the date of invoice. The warranty is extended only to the original purchaser. It shall not apply to any product or parts which have been damaged on account of improper installation, improper connections, misuse, accident or abnormal conditions of operation.

For claims under the warranty please contact your local dealer. You may also send the instrument direct to our works, enclosing the invoice copy and by giving reasons for the claim. You would be solely liable for freight costs.

1 Safety Instructions



Warning!

- Read the operating instructions carefully before use.
- Ensure that only trained staff use the instrument.



Protective ground contact!

- Make sure the socket is earthed (protective ground contact) before use.

- When work , wear the personal guard to avoid the risk from

- Splashing of liquids.
- Mechanical vibration may leads to glass breakage.
- Being trapped of parts of the body, hair or pieces of clothing in the moving parts.
- During manual manipulation:
 - Keep medium vessels vertical as much as possible.
 - Mix the medium in the minimum pressure required.
 - Use attachment supplied when mixing a mass of different medium.
- Set up the instrument in a spacious area on an stable, clean, non-slip, dry and fireproof surface. Do not operate the instrument in explosive atmospheres, with hazardous substances or under water.

- Put the speed control knob to the lowest position. Gradually increase the speed (MX-S), reduce the speed if :
 - The medium splashes out of the vessel because of too high speed.
 - The instrument is not running smoothly or moves on the stage.
- Safe operation is only guaranteed with the accessories described in the “Accessories” chapter.
- Accessories must be securely attached to the instrument and can not come off by themselves.
- Check the instrument and accessories before hand for damage each time you use them. Do not use damaged components.
- Place a vessel at the centre of the mixing head and spread out others uniformly.
- A sharp edged vessel can wear-out the mixing head.
- Always disconnect the plug before fitting accessories.
- The instrument may heat up when in use.
- The operating voltage required on the label and voltage supply network must match.
- Protect the instrument and accessories from bumps and impacts.
- The instrument may only be opened by qualified and trained technician only.

2 Proper use

The instrument is designed for mixing liquids in schools, laboratories or factories. This device is not suitable for use in residential areas or other areas that may cause danger to the user or instrument as mentioned in Chapter 1.

3 Inspection

3.1 Receiving Inspection

Unpack the equipment carefully and check for any damages which may have arisen during transport. If it happens, please contact manufacturer for technical support.



Note:

If there is any apparent damage to the system, please do not plug it into the power line.

3.2 Listing of Items

The packing includes the following items:

Items	Qty
Main unit	1
User manual	1

Table 1

Please retain all packing materials. Additional shipping charges will be assessed if packing is unavailable for the system warranty service.

4 Trial run




Figure 1



Note:

Ensure the two-direction switch must be in the middle (off state) before powering on.

If these operations above are normal, the instrument is ready to operate following the “operation” chapter. If these operations are not normal, the instrument may be damaged. please contact manufacturer for technical support.

5 Operation




Figure 2

5.1 Continuous operating



Figure 3

- Ensure the speed control knob to the lowest position Put the instrument on the stable and safe place and plug in the main power.

- Push the two-direction switch leftward to ON and power LED is lit. Then the instrument begins to work.
- Turn the speed control knob to set the rated speed (MX-S).
- If use the Nonstandard accessories, please ensure the speed in | indicating range (MX-S)
- Push the two-direction switch rightward gently to OFF and LED power is closed and the instrument is turned off.

5.2 Touch operating




Figure 4

- Ensure the speed control knob to the lowest position Put the instrument on the stable and safe place and plug in the main power.
- Push the two-direction switch rightward to the touch mode.
- Turn the speed control knob to set the rated speed (MX-S).
- If a test tube is pressed into the mixing head vertically, the instrument will begin to work and the power LED will be lit.
- Push the two-direction switch leftward gently to OFF and the instrument is turned off.



Note:

- Place vessels such as a test tube as vertical as possible into the centre of the mixing head in order to mixing smoothly.
- Adjust the speed control knob slowly in order to have the instrument work smoothly.

6 Maintenance and cleaning

- Proper maintenance can keep instruments working in a good state and lengthen its life time.
- Be careful not spray the cleanser into the instrument when cleaning.
- Unplug the power line when cleaning.
- Only use cleanser that we advised as below:

Dyes	Isopropyl alcohol
Construction materials	Water containing tenside / isopropyl alcohol
Cosmetics	Water containing tenside / isopropyl alcohol
Foodstuffs	Water containing tenside
Fuels	Water containing tenside

Table 2

For materials which are not listed, please request information

from manufacturer.

- Before using other method for cleaning or decontamination, the user must ascertain with the manufacturer that this method does not destroy the instrument.
- Wear the proper protective gloves during cleaning of the instruments.

NOTE:

- Electrical instruments may not be placed in the cleansing agent for the purpose of cleaning.
- The instrument is a sophisticated electronic products, avoiding severe shake in work and transportation
- The instrument must be cleaned and put it into the initial packaging carton before sending to service for repair, avoiding the contamination of hazardous.
- Switch OFF the instrument and put it in the dry, clean, stable place at room temperature in long-term disuse.

7 Associated standards and regulations

Construction in accordance with the following safety standards:

EN 61010-1

Construction in accordance with the following EMC standards:

EN 61326-1

Associated EU guidelines:

EMC-guidelines: 89/336/EWG

Instrument guidelines: 73/023/EWG

Items	Parameters
Permitted relative humidity	80%
Protection class acc. to DIN 60529	IP21

Table 3

8 Technical data

Items	Parameters
Voltage [VAC]	220-230 / 110-120
Frequency [Hz]	50 / 60
Power [W]	60
Shaking movement	Orbital
Orbital diameter [mm]	4
Motor type	Shaded-pole motor
Motor rating input [W]	58
Motor rating output [W]	10
Permissible ON time	100%
Speed range [rpm]	0-2500(MX-S)/ 2500(MX-F)
Speed display	Scale
Run type	Continuous/touch operation
Dimensions (mm)	127 × 130 × 160
Weight [kg]	3.5
Permitted ambient temperature[° C]	5-40







9 Products and Accessories

Cat.No.	Description
821100019999	SCI-FS fixed speed Vortex, with VT 1.1 Standard Attachment, Cn plug, 220v/50Hz
821100029999	SCI-FS fixed speed Vortex, with VT 1.1 Standard Attachment, Cn plug, 220v/60Hz
821100039999	SCI-FS fixed speed Vortex, with VT 1.1 Standard Attachment, USA plug, 110v/50Hz
821100049999	SCI-FS fixed speed Vortex, with VT 1.1 Standard Attachment, USA plug, 110v/60Hz
821100059999	SCI-FS fixed speed Vortex, with VT 1.1 Standard Attachment, Euro plug, 220v/50Hz

Vortex Mixers (Adjustable and Fixed Speed Models)

821100069999	SCI-FS fixed speed Vortex, with VT 1.1 Standard Attachment, UK plug, 220v/50Hz	821200049999	SCI-VS adjustable speed Vortex, with VT 1.1 Standard Attachment, USA plug, 110v/60Hz
821100079999	SCI-FS fixed speed Vortex, with VT 1.1 Standard Attachment, Euro plug, 220v/60Hz	821200059999	SCI-VS adjustable speed Vortex, with VT 1.1 Standard Attachment, Euro plug, 220v/50Hz
821100089999	SCI-FS fixed speed Vortex, with VT 1.1 Standard Attachment, UK plug, 220v/60Hz	821200069999	SCI-VS adjustable speed Vortex, with VT 1.1 Standard Attachment, UK plug, 220v/50Hz
821100159999	SCI-FS fixed speed Vortex, with VT 1.1 Standard Attachment, USA plug, 220v/60Hz	821200099999	SCI-VS adjustable speed Vortex, with VT 1.1 Standard Attachment, Euro plug, 220v/60Hz
821100169999	SCI-FS fixed speed Vortex, with VT 1.1 Standard Attachment, USA plug, 220v/50Hz	821200109999	SCI-VS adjustable speed Vortex, with VT 1.1 Standard Attachment, UK plug, 220v/60Hz
821200019999	SCI-VS adjustable speed Vortex, with VT 1.1 Standard Attachment, Cn plug, 220v/50Hz	821200159999	SCI-VS adjustable speed Vortex, with VT 1.1 Standard Attachment, USA plug, 220v/60Hz
821200029999	SCI-VS adjustable speed Vortex, with VT 1.1 Standard Attachment, Cn plug, 220v/60Hz	821200169999	SCI-VS adjustable speed Vortex, with VT 1.1 Standard Attachment, USA plug, 220v/50Hz
821200039999	SCI-VS adjustable speed Vortex, with VT 1.1 Standard Attachment, USA plug, 110v/50Hz		

Vortex Mixers (Adjustable and Fixed Speed Models)

Accessories			
	Cat.No and Designation	Serviceable range	Operation mode
	18900034 VT1.1 standard top	For test tubes and small vessels upto diameter 30 mm	Continuous /touch operation
	18900035 VT1.3 universal top plate	100 mm-plate attachment with rubber pad	Continuous operation
	18900020 VT1.3.1 tube adapte	For 48 Eppendorf tubes	Continuous operation
	18900021 VT1.3.2 tube adapter	For 18 test tubes Ø10mm for VT 1.3	Continuous operation
	18900022 VT1.3.3 tube adapter	For 12 test tubes Ø12mm , for VT 1.3	Continuous operation
	18900023 VT1.3.4 tube adapter	For 8 test tubes Ø16mm, for VT 1.3	Continuous operation



	18900024 VT1.3.5 tube adapter	For 8 test tubes Ø20mm, for VT 1.3	Continuous operation
	18900043 VT1.3.6, platform pad	For < Ø99 mm tubes and small vessels	Touch operation

Table 4

NOTE: VT1.3 ~ VT1.3.5 adapter and pad are only suitable to MX-S Adjustable speed Vortex Serial Number and pattern number and accessories names on the type plate must be supplied when ordering other accessories.

Changing attachment




Figure 5