



SYMMETRIC OCHRATOXIN ES

LATERAL FLOW TEST KIT

for the quantitative determination of Ochratoxin in grains, cereals and nuts.

ProGnosis Biotech S.A. is ISO 9001:2015 certified by TÜV Hellas (TÜV NORD).

Use only the current version of Product Data Sheet enclosed with the kit.

Symmetric OCHRATOXIN ES, S6424/S6448, is a Lateral Flow Test kit for the quantitative determination of Ochratoxin A in grains, cereals and nuts.

This kit contains all reagents required for 24 or 48 reactions.

Matrices:

Type I: Barley, Brown Rice, Buckwheat, Buckwheat Flour*, Corn, Corn Flour, Dried Dates, Dried Figs, Dried Prunes, Malt, Oats, Raisins, Sesame, Soy, Soybeans, Triticale, Wheat, Wheat Flour, White Rice,

Rye, Spelt, Pea flour, Pistachio, Pistachio paste, Hazelnut paste

Type II: Almond, Cashews, Sunflower seeds, Corn Gluten Meal, DDGS

- Sample preparation: extraction
- Test time (reaction time after samples and reagents preparation): 10min
- Range: 0 - 20ppb
- Shelf life: 12 months
- Storage: 2-8°C
- Testing temperature: 18-30°C

**Contact our company for sample preparation procedure*

This is an electronic version, please verify always the last one included in the kit.

Specifications

- The LOD of the method is: 0,45ppb (Type I) & 0.7 ppb (Type II)
- The LOQ of the method is: 0.7ppb (Type I) & 1 ppb (Type II)
- Cross-reactivity: The cross-reaction of the anti-Ochratoxin antibody with Ochratoxin A and B is 100 and <0.1% respectively.

1. Description

Symmetric OCHRATOXIN ES is an innovative Lateral Flow test, utilizing state-of-the-art features for the quantitative detection of Ochratoxin A in grains, cereals and nuts.

2. General Information

Ochratoxins are a group of mycotoxins produced by some Aspergillus species (mainly *A. ochraceus*, but also by 33% of *A. niger* industrial strains) and some Penicillium species, especially *P. verrucosum* and *P. carbonarius*. Ochratoxin A (OTA) is the most prevalent and relevant fungal toxin of this group, while ochratoxins B and C are of lesser importance. OTA is a potent nephrotoxin and causes both acute and chronic effects in the kidneys of all mammalian species tested. It is also genotoxic (damages DNA) and teratogenic (damages the foetus) and is considered a probable carcinogen, causing renal carcinoma and other cancers in a number of animal species, Most controlling government agencies worldwide have regulations regarding the amount of aflatoxins allowable in human and animal foodstuffs. Accurate and rapid determination of the presence of OTA in commodities is of paramount importance.

3. Principle of the Method

The quantitative lateral flow test is based on the immunochromatography assay principles. The wells of the microtiter strips contain OTA specific antibodies conjugated to colloidal gold. Diluted extract is added into the well. A dipstick with two capture lines, test and control, is dipped into the well. The suspended mixture starts flowing vertically on the dipstick and passes through the two lines. While running, OTA (if it is present) binds to the antibodies. A valid test should always have the upper control line red. If the sample is free of OTA, a color development occurs at the test line, indicating the absence of OTA in the sample. On the contrary, the presence of OTA in the sample will cause a reduced colored signal at the test line. The test line color intensity is indirectly proportionate to the concentration of OTA present in the samples. By utilizing S-Flow software and the symmetric quantification technology, OTA is accurately quantified.

4. Reagents Provided

Symmetric Ochratoxin ES kit contains sufficient reagents and materials for 24/48 reactions.

Reagents (Store at 2-8°C)	Quantity for 24 reactions	Quantity for 48 reactions
Pots each with 1 strip of 8 reagent microwells and 8 dipsticks	3	6
Sample Diluent Tubes	24	48
High Range Solution	1	1

5. Materials required but not provided

- A grinder sufficient to render sample to particle size of fine instant coffee
- Balance with 0 - 50g measuring capability and Graduated cylinder - 50ml
- Ethanol (9.75ml reagent grade per sample) and Deionized water
- Filter Paper Whatman #1 or equivalent, Filter Funnel and Miscellaneous laboratory plastic or glass tubes 5 - 15ml or
- Mini centrifuge (spin) and plastic tubes 1,5 or 2ml
- Tube roller or Vortex mixer
- 100 - 300µl adjustable micropipettes (single or multi channel) with disposable tips
- **S-Flow** software along with matching scanner device

6. Storage Instructions

Store kit components between 2 - 8°C. Do not freeze any components provided. Reseal the unused strips in the storing tube together with the desiccant bag provided. The expiry date of the kit and reagents is stated on their labels and no quality guarantee is accepted after the expiration date. The expiry of the kit components can only be guaranteed if the components are stored properly and the reagent is not contaminated due to prior handling. Do not interchange individual components between kits of different lot numbers.

7. Safety and Precautions for use

Let the reagents warm to room temperature (21 - 25°C) before the analysis (at least half an hour) and cover them when not in use. Use a clean disposable plastic pipette tip for each reagent, to avoid cross contamination.

8. Preparation of the extraction solvent

Prepare the Extraction Solution (65% Ethanol) by adding 35ml of distilled or deionized water to 65ml of ethanol (reagent grade) and transfer it into a glass bottle.

The quantity of extraction solvent (100ml of 65% ethanol) is sufficient for 6 samples (5gr each). A user may prepare smaller or larger amounts of the extraction solvent maintaining the ratio of 6.5 parts ethanol to 3.5 parts distilled or deionized water by volume. Prepare a new extraction solvent prior to use

9. Sample preparation

1. The sample must be collected according to established sampling techniques. Grind a representative sample to the particle size of fine instant coffee (85-95% passes through a 20 mesh screen).
2. Weigh out a 5g ground portion of the sample and add 15ml of the **65% ethanol** for each sample to be tested. Mix using a tube roller for 5 minutes (or vortex for 2min). **The ratio of sample to Extraction Solution is 1:3 (w/v). To achieve good homogenization, ensuring that any portion of the sample will be representative of the whole, weight at least 20gr of the sample.**
3. Allow the particulate matter to settle. Filter the extract through a Whatman #1 filter paper (or equivalent) and collect the filtrate. Alternatively, centrifuge 1ml of the extract for 2min using a mini centrifuge (spin).
4. **Type I - Type II (quantification range 0.7-20ppb):** add **100µl** of extract (supernatant) into the Sample diluent tube provided and mix well. (Use this **diluted extract** within **30 minutes**).
5. **If sample exceeds 20ppb** first prepare the diluted extract (9.4.). All the samples should be diluted as described in 9.4.step. The diluted extract should be diluted 5 or 10-folds with High Range Solution for re-analysis. Use the second dilution within 30 minutes.
6. Choose 5X Dilution Type and set the suitable dilution factor type to multiply the results by **5or 10**.

5x Dilution Matrix Type	
Dilution Factor 5 (Quantify from 3 up to 50ppb)	100µl of the diluted sample (9.4) +400µl of High Range Solution
Dilution Factor 10 (Quantify from 6 up to 100ppb)	100µl of the diluted sample (9.4) +900µl of High Range Solution

10. Method Procedure

1. Before opening the reagents, take the kit out of the fridge and wait until the temperature of the reagents reaches the ambient temperature.
2. Download and/or set the kit's **lot number**, as provided in the Quality Assurance Certificate and then set the suitable **type** and **Dilution Factor**.
3. Open one plastic pot and take out as many test strips and microwells as samples to be tested .
4. The pot with dipsticks should **always be well closed** after reagents have been taken out.
5. Dispense **100µl of diluted filtrate** into the microwell and pipette **up and down 4 times** to completely mix the lyophilized gold particles in the sample, while avoiding bubbles. The sample should turn into a **uniform pink color**. In case of more than 2 samples, an 8 channel multipipette should be used.
6. Place the appropriate number of sticks into microwells **immediately** and set timer for 10 minutes.
7. When the 10 minutes are over, take the dipsticks out of the microwells and remove the white cotton sample-pad of the stick **immediately**. Touch the stick with your hand from the colorful pad and remove the white pad with your hands. Do not use a paper towel or any other material.
8. Place the stick inside the plastic holder in order to be scanned. In case of EPSON scanner, the sticks must be facing down (inverted) and the colored side must be facing the orange sticker. **NOTE:** The sticks should be scanned within 2 minutes after the sample-pad removal.
9. The software will use a Lot specific curve to calculate the results (ppb) according to the matrix sample type. A simple visual interpretation of the stick is NOT possible

11. Performance Evaluation

11.1 Reference Materials


Several reference materials are being used for the evaluation of each product of ProGnosis Biotech S.A. in the context of Quality Control performed by Quality Control Department. Please request a validation report, including the results, at exports@prognosis-biotech.com.

11.2 Proficiency Tests

All products participate frequently in Proficiency Tests. For more information, visit the individual product page in our website: www.prognosis-biotech.com

12. Method Summary

Total method time: 10 minutes



All immune assays supplied by ProGnosis Biotech S.A., are warranted to meet or exceed our published specification when used under normal conditions in your laboratory. If the product fails during the stated period, a replacement product will be issued.

ProGnosis Biotech S.A. makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. There is no warranty of merchantability of this product, or of the fitness of the product for any purpose. ProGnosis Biotech S.A. shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product. This method is considered to be a screening method, before a legal action, samples detected as positives must be confirmed with a confirmation method. This product is meant to be used only For Research or Manufacturing use and by qualified technicians.

S6424-S6448 Manual_Symmetric_Ochratoxin ES_v6_en

VERSION N6

CAT.NUMBER: S6424/S6448

STORAGE: 2-8°C



LATERAL FLOW TEST KIT

for the quantitative determination of Ochratoxin in grains, cereals and nuts.

ProGnosis Biotech S.A. is ISO 9001:2015 certified by TÜV Hellas (TÜV NORD).

Use only the current version of Product Data Sheet enclosed with the kit.

Symmetric Ochratoxin ES, S6424/S6448, is a Lateral Flow Test kit for the quantitative determination of Ochratoxin in grains, cereals and nuts.

This kit contains all reagents required for 24 or 48 reactions.

Matrices:

Type I: Barley, Brown Rice, Buckwheat, Buckwheat Flour*, Corn, Corn Flour, Dried Dates, Dried Figs, Dried Prunes, Malt, Oats, Raisins, Sesame, Soy, Soybeans, Triticale, Wheat, Wheat Flour, White Rice, Rye, Spelt, Pea flour, Pistachio, Pistachio paste, Hazelnut paste
Type II: Almond, Cashews, Sunflower seeds, Corn Gluten Meal, DDGS

- Sample preparation: extraction
- Test time (reaction time after samples and reagents preparation): 10min
- Range: 0 - 20ppb
- Shelf life: 12 months
- Storage: 2-8°C
- Testing temperature: 18-30°C

**Contact our company for sample preparation procedure*



www.prognosis-biotech.com
e: exports@prognosis-biotech.com
t: +30 2410 623922
Farsalon 153 | 41335 Larissa, Greece

