# **RapidChek**<sup>®</sup> SELECT™ *Salmonella* Enteritidis Quick Guide Chicken Carcass Rinse Samples

# RapidChek® SELECT™ Salmonella Enteritidis Test Kit, 10001396

Test Kit Includes:

- 250 test strips in a desiccated canister
- 50 transfer pipettes
- 50 sample tubes
- Package Insert

Does not include media.

## RapidChek® SELECT™ Salmonella Enteritidis Media System, 10001397

- 500 g bottle of RapidChek® SELECT™ Primary Media
- 250 mL bottle of RapidChek<sup>®</sup> SELECT™ Primary Supplement
- 4 x 10 g bottles of RapidChek® SELECT™ Secondary Media Secondary Media

Media system is sufficient for 25 L of primary media and 500 mL of secondary media.

# RapidChek® SELECT™ Salmonella Enteritidis Test System, 10001398

- 250 test strips in a desiccated canister
- 250 transfer pipettes
- 250 sample tubes
- 500 g bottle of RapidChek® SELECT™ Primary Media
- 250 mL bottle of RapidChek® SELECT™ Primary Supplement
- 4 x 10 g bottles of RapidChek® SELECT™ Secondary Media
- Package Insert

Test system is sufficient for 25 L of primary media and 500 mL of secondary media.



#### Also Available:

- RapidChek® CONFIRM™ Salmonella Enteritidis
  IMS
- RapidChek® SELECT™ Salmonella

- RapidChek® Listeria
- RapidChek® *Listeria* NextDay™
- RapidChek® E. coli O157 (including H7)
- RapidChek® CONFIRM™ non-O157 STEC IMS

RapidChek®. Simply Accurate.



# RapidChek® SELECT™ Salmonella Enteritidis Quick Guide Chicken Carcass Rinse Samples



#### **Primary Media Preparation**

Prepare RapidChek® SELECT™ Salmonella primary media with supplement.\* Note this media should be prepared with 2 x concentration of both the powder base and supplement.

Number of Samples	Sterilized or Distilled** Water (mL)	RapidChek® SELECT™ Primary Media (Grams)	RapidChek® SELECT™ Primary Supplement (mL)
50	1500 (1.5 L)	60	30
30	900	36	18
20	600	24	12
10	300	12	6

<sup>\*</sup>The supplement must be stored at 2 - 8 °C when not in use.

# 2 Options for Media **Preparation**

Non-Autoclaved Option: Add dehydrated media and supplement to sterile water.

Shelf Life: Three hours within preparation.


Autoclaved Option: Add dehydrated media to DI water and autoclave. Just prior to use add supplement.

Shelf Life: Four weeks at 2 - 8 °C or RT without supplement.



#### **Primary Sample Enrichment**

Add 30 mL of prepared 2 x RapidChek® SELECT™ Salmonella Primary Media to 30 mL of raw chicken rinse. Incubate 16 - 22 hours at 42 °C.



#### Secondary Media Preparation\*

Prepare RapidChek® SELECT™ Salmonella Secondary Media and aliquot 1 mL of the prepared media into supplied tubes, 1 tube per each enriched sample.

Number of Samples	Sterilized Water (mL)	RapidChek® SELECT™ Secondary Media (Grams)
50	50	3.7
30	30	2.2
20	20	1.5
10	10	0.74

<sup>\*</sup>Prepared secondary media can be stored two weeks at 2 - 8 °C




### **Secondary Sample Enrichment**

Transfer 0.1 mL of the primary enriched sample into the prepared tubes of 1 mL of RapidChek® SELECT™ Salmonella Secondary Media and incubate an additional 16 - 22 hours at 42 °C.



Test

Place test strips into tubes with arrows facing down.



Read Result After 10 minutes.

read result.

1 Line = Negative 2 Lines = Positive

<sup>\*\*</sup>If you are using DI water, the re-hydrated media must be autoclaved and cooled prior to adding the supplement.