RapidChek[®] SELECT™ *Salmonella* Enteritidis Quick Guide Drag Swab Samples

RapidChek® SELECT™ Salmonella Enteritidis Test Kit, 10001396

Test Kit Includes:

- 250 test strips in a desiccated canister
- 50 transfer pipettes
- 50 sample tubes
- Package Insert

Does not include media.

RapidChek® SELECT™ Salmonella Enteritidis Media System, 10001397

- 500 g bottle of RapidChek® SELECT™ Primary Media
- 250 mL bottle of RapidChek® SELECT™ Primary Supplement
- 4 x 10 g bottles of RapidChek® SELECT™ Secondary Media Secondary Media

Media system is sufficient for 25 L of primary media and 500 mL of secondary media.

RapidChek® SELECT™ Salmonella Enteritidis Test System, 10001398

- 250 test strips in a desiccated canister
- 250 transfer pipettes
- 250 sample tubes
- 500 g bottle of RapidChek® SELECT™ Primary Media
- 250 mL bottle of RapidChek® SELECT™ Primary Supplement
- 4 x 10 g bottles of RapidChek® SELECT™ Secondary Media
- Package Insert

Test system is sufficient for 25 L of primary media and 500 mL of secondary media.




Also Available:

- RapidChek® CONFIRM™ Salmonella Enteritidis IMS
- RapidChek® SELECT™ Salmonella
- RapidChek® Listeria
- RapidChek® Listeria NextDay™
- RapidChek® E. coli O157 (including H7)
- RapidChek® CONFIRM™ non-O157 STEC IMS

RapidChek®. Simply Accurate.



RapidChek® SELECT™ Salmonella Enteritidis Quick Guide **Drag Swab Samples**



Primary Media Preparation

Prepare RapidChek® SELECT™ Salmonella primary media with supplement.*

Number of Samples	Sterilized or Distilled** Water (Liters)	Rapid Chek® SELECT™ Primary Media (Grams)	Rapid Chek® SELECT™ Primary Supplement (mL)
100	10	200	100
20	2	40	20
10	1	20	10
1	0.1 (100 mL)	2	1

^{*}The supplement must be stored at 2 - 8 °C when not in use.

2 Options for Media **Preparation**

Non-Autoclaved Option: Add dehydrated media and supplement to sterile water.

Shelf Life: Three hours within preparation.


Autoclaved Option: Add dehydrated media to DI water and autoclave. Just prior to use add supplement.

Shelf Life: Four weeks at 2 - 8 °C or RT without supplement.



Primary Sample Enrichment

Add 100mL of prepared RapidChek® SELECT™ Salmonella Primary Media with supplement to the drag swab sample. Incubate 16 - 22 hours at 42 °C.



Secondary Media Preparation*

Prepare RapidChek® SELECT Secondary Media and aliquot 2 mL of the prepared media into supplied tubes, 1 tube per each enriched sample.

Number of Samples	Sterilized Water (mL)	Rapid Chek® SELECT™ Secondary Media (Grams)
100	200	14.8
20	40	3.0
10	20	1.5
1	2	0.15

^{*}Prepared secondary media can be stored two weeks at 2 - 8 °C




Secondary Sample Enrichment

Transfer 0.2 mL of the primary enriched sample into the prepared tubes of 2 mL of RapidChek® SELECT™ Salmonella Secondary Media and incubate an additional 16 - 22 hours at 42 °C.



Test

Place test strips into tubes with arrows facing down.



Read Result After 10 minutes.

read result. 1 Line = Negative 2 Lines = Positive

^{**}If you are using DI water, the rehydrated media must be autoclaved and cooled prior to adding the supplement.