

Catalog number	Unit Size	Precision	Description	Price (NGN)
4421120	Each	2 μ L	Range: 0.2-1.0mL, \pm 6.0 μ L	Contact Us >

Full specifications ▼

Certifications/Compliance	Individual calibration certificate with serial number
Description	Range: 0.2-1.0mL, \pm 6.0 μ L
Type	Dispenser
Precision	2 μ L
Unit Size	Each

4421130	Each	4 μ L	Range: 0.4-2.0mL, \pm 12.0 μ L	Contact Us >
---------	------	-----------	--------------------------------------	---------------------------------

Full specifications ▼

Certifications/Compliance	Individual calibration certificate with serial number
Description	Range: 0.4-2.0mL, \pm 12.0 μ L

Catalog number	Unit Size	Precision	Description	Price (NGN)
----------------	-----------	-----------	-------------	-------------

Type	Dispenser
Precision	4 µL
Unit Size	Each

4421140	Each	10 µL	Range: 1-5mL, ±30.0µL	Contact Us >
---------	------	-------	-----------------------	---------------------------------

Full specifications ▼

Certifications/Compliance	Individual calibration certificate with serial number
Description	Range: 1-5mL, ±30.0µL
Type	Dispenser
Precision	10 µL
Unit Size	Each

4421150	Each	20µL	Range: 2-10mL, ±60.0µL	Contact Us >
---------	------	------	------------------------	---------------------------------

Full specifications ▼

Certifications/Compliance	Individual calibration certificate with serial number
----------------------------------	---

Etcon Analytical and Environmental Systems & Services Limited

Office Address: No 7, Abeokuta Street, Off Murtala Mohammed Way, Ebutte-Metta, Yaba, Lagos, Nigeria.
 Tel: +234-1-4536335, +234-8021422237 | E-mail: info@etconanalytical.com. | Website: www.etconanalytical.com

Catalog number	Unit Size	Precision	Description	Price (NGN)
----------------	-----------	-----------	-------------	-------------

Description	Range: 2-10mL, ±60.0µL
Type	Dispenser
Precision	20µL
Unit Size	Each

4421160	Each	60 µL	Range: 5-30mL, ±180µL	Contact Us >
---------	------	-------	-----------------------	---------------------------------

Full specifications ▼

Certifications/Compliance	Individual calibration certificate with serial number
Description	Range: 5-30mL, ±180µL
Type	Dispenser
Precision	60 µL
Unit Size	Each

4421170	Each	120µL	Range: 10-60mL, ±360µL	Contact Us >
---------	------	-------	------------------------	---------------------------------

Full specifications ▼

Etcon Analytical and Environmental Systems & Services Limited

Office Address: No 7, Abeokuta Street, Off Murtala Mohammed Way, Ebutte-Metta, Yaba, Lagos, Nigeria.
 Tel: +234-1-4536335, +234-8021422237 | E-mail: info@etconanalytical.com. | Website: www.etconanalytical.com

Catalog number	Unit Size	Precision	Description	Price (NGN)
----------------	-----------	-----------	-------------	-------------

Certifications/Compliance	Individual calibration certificate with serial number
Description	Range: 10-60mL, $\pm 360\mu\text{L}$
Type	Dispenser
Precision	120 μL
Unit Size	Each

