

Thermo Scientific™

Finnpipette™ F1 Good Laboratory Pipetting (GLP) Kits



Related applications:

[Pipettes and Pipette Tips](#)

Thermo Scientific™ Finnpipette™ F1 GLP Kits offer a variety of single channels in a range of volumes that are packaged with compatible Finntips, and a convenient pipette stand to support good laboratory pipetting in your lab. Select the kit that is right for you, and take advantage of the Finnpipette F1 performance and lightweight design to maximize your comfort.

[Contact us for support >](#)

Catalog number	Unit Size	Description	Volume (Metric)	Compatible Tips	Price (NGN)
4700850N	Each	Finnpipette F1 GLP Kit 1: 1-10µL micro, 10-100µL and 100-1000µL	1 to 1000 µL	Finntip Flex tips: 10; 200; 1000	Contact Us >

Full specifications ▼

Catalog number	Unit Size	Description	Volume (Metric)	Compatible Tips	Price (NGN)
----------------	-----------	-------------	-----------------	-----------------	--------------

Compatible Tips	Finntip Flex tips: 10; 200; 1000
Description	Finnpipette F1 GLP Kit 1: 1-10µL micro, 10-100µL and 100-1000µL
Volume (Metric)	1 to 1000 µL
Unit Size	Each

4700860N	Each	Finnpipette F1 GLP Kit 2: 0.2-2µL micro, 2-20µL, 20-200µL and 100-1000µL	0.2 to 1000 µL	Finntip Flex tips: 10; 200; 1000	Contact Us >
----------	------	--	----------------	-------------------------------------	---------------------------------

Full specifications ▼

Compatible Tips	Finntip Flex tips: 10; 200; 1000
Description	Finnpipette F1 GLP Kit 2: 0.2-2µL micro, 2-20µL, 20-200µL and 100-1000µL
Volume (Metric)	0.2 to 1000 µL
Unit Size	Each

4700865N	Each	Finnpipette F1GLP Kit 3: 10-100µL, 100-1000µL and 1-10mL	10 to 10,000 µL	Finntip Flex tips: 200; 1000; FinnTip 10mL	Contact Us >
----------	------	--	-----------------	--	---------------------------------

Catalog number	Unit Size	Description	Volume (Metric)	Compatible Tips	Price (NGN)
----------------	-----------	-------------	-----------------	-----------------	--------------

Full specifications ▼

Compatible Tips	Finntip Flex tips: 200; 1000; FinnTip 10mL
Description	Finnpipette F1GLP Kit 3: 10-100µL, 100-1000µL and 1-10mL
Volume (Metric)	10 to 10,000 µL
Unit Size	Each

4701060N	Each	Finnpipette F1 GLP Kit 4: 2-20µL, 20-200µL and 100-1000µL	2 to 1000 µL	Finntip Flex tips: 200; 1000µL	Contact Us >
----------	------	---	--------------	--------------------------------	---------------------------------

Full specifications ▼

Compatible Tips	Finntip Flex tips: 200; 1000µL
Description	Finnpipette F1 GLP Kit 4: 2-20µL, 20-200µL and 100-1000µL
Volume (Metric)	2 to 1000 µL
Unit Size	Each