

Thermo Scientific™

Finnpipette™ F1 Variable Volume Pipettes



Related applications:

[Pipettes and Pipette Tips](#)

The enhanced volume adjustment on Thermo Scientific™ Finnpipette™ F1 Variable Volume Pipettes is now silent and effortless with textured adjustment knob for a firmer feel when setting volumes. The adjustable finger rest, and lightweight design maximizes comfort without sacrificing performance

[Contact us for support >](#)

Etcon Analytical and Environmental Systems & Services Limited

Office Address: No 7, Abeokuta Street, Off Murtala Mohammed Way, Ebutte-Metta, Yaba, Lagos, Nigeria.

Tel: +234-1-4536335, +234-8021422237 | E-mail: info@etconanalytical.com. | Website: www.etconanalytical.com

Catalog number	Unit Size	Description	Volume (Metric)	Compatible Tips	Price (NGN)
4641020N	Each	Finnpipette F1	0.5 to 5 µL	Finntip 10, Finntip Flex 10	Contact Us >

Full specifications ▼

Compatible Tips	Finntip 10, Finntip Flex 10
Description	Finnpipette F1
Increments	0.01µL
Color	Pink
Volume (Metric)	0.5 to 5 µL
Unit Size	Each

4641010N	Each	Finnpipette F1	0.2 to 2µL	Finntip 10, Finntip Flex 10	Contact Us >
----------	------	----------------	------------	-----------------------------	---------------------------------

Full specifications ▼

Compatible Tips	Finntip 10, Finntip Flex 10
Description	Finnpipette F1

Catalog number	Unit Size	Description	Volume (Metric)	Compatible Tips	Price (NGN)
----------------	-----------	-------------	-----------------	-----------------	--------------

Increments	0.002µL
Color	Pink
Volume (Metric)	0.2 to 2µL
Unit Size	Each

4641030N	Each	Finnpipette F1	1 to 10 µL	Finntip 10, Finntip Flex 10	Contact Us >
----------	------	----------------	------------	-----------------------------	---------------------------------

Full specifications ▼

Compatible Tips	Finntip 10, Finntip Flex 10
Description	Finnpipette F1
Increments	0.02µL
Color	Pink
Volume (Metric)	1 to 10 µL
Unit Size	Each

Catalog number	Unit Size	Description	Volume (Metric)	Compatible Tips	Price (NGN)
4641040N	Each	Finnpipette F1	1 to 10 µL	Finntip 250, Finntip Flex 200	Contact Us >

Full specifications ▼

Compatible Tips	Finntip 250, Finntip Flex 200
Description	Finnpipette F1
Increments	0.02µL
Color	Yellow
Volume (Metric)	1 to 10 µL
Unit Size	Each

4641050N	Each	Finnpipette F1	2 to 20 µL	Finntip 50	Contact Us >
----------	------	----------------	------------	------------	---------------------------------

Full specifications ▼

Compatible Tips	Finntip 50
Description	Finnpipette F1

Catalog number	Unit Size	Description	Volume (Metric)	Compatible Tips	Price (NGN)
----------------	-----------	-------------	-----------------	-----------------	--------------

Increments	0.02µL
Color	Turquoise
Volume (Metric)	2 to 20 µL
Unit Size	Each

4641060N	Each	Finnpipette F1	2 to 20 µL	Finntip 250, Finntip Flex 200	Contact Us >
----------	------	----------------	------------	-------------------------------	---------------------------------

Full specifications ▼

Compatible Tips	Finntip 250, Finntip Flex 200
Description	Finnpipette F1
Increments	0.02µL
Color	Yellow
Volume (Metric)	2 to 20 µL
Unit Size	Each

Catalog number	Unit Size	Description	Volume (Metric)	Compatible Tips	Price (NGN)
4641070N	Each	Finnpipette F1	10 to 100 µL	Finntip 250, Finntip Flex 200	Contact Us >

Full specifications ▼

Compatible Tips	Finntip 250, Finntip Flex 200
Description	Finnpipette F1
Increments	0.2µL
Color	Yellow
Volume (Metric)	10 to 100 µL
Unit Size	Each

4641080N	Each	Finnpipette F1	20 to 200 µL	Finntip 250, Finntip Flex 200	Contact Us >
----------	------	----------------	--------------	-------------------------------	---------------------------------

Full specifications ▼

Compatible Tips	Finntip 250, Finntip Flex 200
Description	Finnpipette F1

Catalog number	Unit Size	Description	Volume (Metric)	Compatible Tips	Price (NGN)
----------------	-----------	-------------	-----------------	-----------------	--------------

Increments	0.2µL
Color	Yellow
Volume (Metric)	20 to 200 µL
Unit Size	Each

4641090N	Each	Finnpipette F1	30 to 300 µL	Finntip 300, FinnTip Flex 300	Contact Us >
----------	------	----------------	--------------	-------------------------------	---------------------------------

Full specifications ▼

Compatible Tips	Finntip 300, FinnTip Flex 300
Description	Finnpipette F1
Increments	1µL
Color	Orange
Volume (Metric)	30 to 300 µL
Unit Size	Each

Catalog number	Unit Size	Description	Volume (Metric)	Compatible Tips	Price (NGN)
4641100N	Each	Finnpipette F1	100 to 1000 μ L	Finntip 1000, FinnTip Flex 1000	Contact Us >

Full specifications ▼

Compatible Tips	Finntip 1000, FinnTip Flex 1000
Description	Finnpipette F1
Increments	1 μ L
Color	Blue
Volume (Metric)	100 to 1000 μ L
Unit Size	Each

4641110N	Each	Finnpipette F1	0.5 to 5 mL	Finntip 5mL	Contact Us >
----------	------	----------------	-------------	-------------	---------------------------------

Full specifications ▼

Compatible Tips	Finntip 5mL
Description	Finnpipette F1

Catalog number	Unit Size	Description	Volume (Metric)	Compatible Tips	Price (NGN)
----------------	-----------	-------------	-----------------	-----------------	--------------

Increments	0.01mL
Color	Green
Volume (Metric)	0.5 to 5 mL
Unit Size	Each

4641120N	Each	Finnpipette F1	1 to 10 mL	Finntip 10mL, Finntip Flex 10mL Ext.	Contact Us >
----------	------	----------------	------------	--------------------------------------	---------------------------------

Full specifications ▼

Compatible Tips	Finntip 10mL, Finntip Flex 10mL Ext.
Description	Finnpipette F1
Increments	0.02mL
Color	Red
Volume (Metric)	1 to 10 mL
Unit Size	Each

Catalog number	Unit Size	Description	Volume (Metric)	Compatible Tips	Price (NGN)
4641130N	Each	Finnpipette F1	5 to 50 µL	FinnTip 50	Contact Us >

Full specifications ▼

Compatible Tips	FinnTip 50
Description	Finnpipette F1
Increments	0.1µL
Color	Turquoise
Volume (Metric)	5 to 50 µL
Unit Size	Each

4641140N	Each	Finnpipette F1	5 to 50 µL	Finntip 250, Finntip Flex 200	Contact Us >
----------	------	----------------	------------	-------------------------------	---------------------------------

Full specifications ▼

Compatible Tips	Finntip 250, Finntip Flex 200
Description	Finnpipette F1

Catalog number	Unit Size	Description	Volume (Metric)	Compatible Tips	Price (NGN)
----------------	-----------	-------------	-----------------	-----------------	--------------

Increments	0.1µL
Color	Yellow
Volume (Metric)	5 to 50 µL
Unit Size	Each

