

Thermo Scientific™

## Finnpipette™ F2 GLP Kits



### Related applications:

[Pipettes and Pipette Tips](#)

### | [SARS-CoV-2 Pathogen Research Solutions](#)

Access everything needed for pipetting in these all-inclusive kits. Thermo Scientific™ Finnpipette™ F2 GLP (Good Laboratory Pipetting) Kits provide durable, high-performance Finnpipette F2 pipetters which offer maximum reliability and comfort. The F2 pipetters have unique features such as AVG (Advanced Volume Gearing mechanism) and a very large display for improved usability and comfort.

[Contact us for support ›](#)

Catalog number	Unit Size	Description	Price (NGN )
4700880	Each	Kit 2: Finnpiettes F2 0.2-2µL micro, 2-20µL, 20-200µL and 100-1000µL	<a href="#">Contact Us &gt;</a>

Full specifications ▼

<b>Compatible Tips</b>	Flex tips: 10, 200, 1000
<b>Description</b>	Kit 2: Finnpiettes F2 0.2-2µL micro, 2-20µL, 20-200µL and 100-1000µL
<b>Color</b>	Multiple
<b>Model</b>	Finnpipette™ F2
<b>Volume (Metric)</b>	0.2 to 1000 µL
<b>Unit Size</b>	Each

4700870	Each	Kit 1: Finnpiettes F2 1-10µL micro, 10-100µL and 100-1000µL	<a href="#">Contact Us &gt;</a>
---------	------	---	---------------------------------

Full specifications ▼

<b>Compatible Tips</b>	Flex tips: 10, 200, 1000
<b>Description</b>	Kit 1: Finnpiettes F2 1-10µL micro, 10-100µL and 100-1000µL

Catalog number	Unit Size	Description	Price (NGN )
----------------	-----------	-------------	--------------

<b>Color</b>	Multiple
<b>Model</b>	Finnpipette™ F2
<b>Volume (Metric)</b>	1 to 1000 µL
<b>Unit Size</b>	Each

4700885	Each	Kit 3: Finnpipettes F2 10-100µL, 100-1000µL and 1-10mL	<a href="#">Contact Us &gt;</a>
---------	------	--	---------------------------------

Full specifications ▼

<b>Compatible Tips</b>	Flex tips: 200, 1000, FT 10mL
<b>Description</b>	Kit 3: Finnpipettes F2 10-100µL, 100-1000µL and 1-10mL
<b>Color</b>	Multiple
<b>Model</b>	Finnpipette™ F2
<b>Volume (Metric)</b>	10 to 10,000 µL
<b>Unit Size</b>	Each

Catalog number	Unit Size	Description	Price (NGN )
4701070	Each	Kit 4: Finnpiettes F2 2-20µL, 20-200µL and 100-1000µL	<a href="#">Contact Us &gt;</a>

#### Full specifications ▼

<b>Compatible Tips</b>	Flex tips: 200, 100-1000µL
<b>Description</b>	Kit 4: Finnpiettes F2 2-20µL, 20-200µL and 100-1000µL
<b>Color</b>	Multiple
<b>Model</b>	Finnpipette™ F2
<b>Volume (Metric)</b>	2 to 1000 µL
<b>Unit Size</b>	Each

#### Product overview

##### Features:

- AVG-volume gearing mechanism for accuracy and precision
- Fully autoclavable for protection
- Soft-touch tip ejector for light tip ejection
- Large display for better vision ergonomics
- Very light pipetting forces

**Etcon Analytical and Environmental Systems & Services Limited**

Office Address: No 7, Abeokuta Street, Off Murtala Mohammed Way, Ebutte-Metta, Yaba, Lagos, Nigeria.  
Tel: +234-1-4536335, +234-8021422237 | E-mail: info@etconanalytical.com. | Website: www.etconanalytical.com

### Kit 1 includes:

---

- Three Thermo Scientific Finnpiquette Pipetters from 1 to 1000 $\mu$ L:
  - Finnpiquette F2
  - 1-10 $\mu$ L
  - 10-100 $\mu$ L
  - 100-1000 $\mu$ L
- Flex tips: 10, 1 x 96; 200, 1 x 96; 1000, 1 x 96
- F-stand
- Good Laboratory Pipetting Guide
- Reagent reservoir demo pack

### Kit 2 includes:

- Four Thermo Scientific Finnpiquette Pipetters from 0.2 to 1000 $\mu$ L:
  - Finnpiquette F2
  - 0.2-2 $\mu$ L
  - 2-20 $\mu$ L univ
  - 20-200 $\mu$ L
  - 100-1000 $\mu$ L
- Flex tips: 10, 12 x 96; 200, 2 x 96; 1000, 1 x 96
- F-stand
- Good Laboratory Pipetting Guide
- Reagent reservoir demo pack

### Kit 3 includes:

- Three Thermo Scientific Finnpiquette Pipetters from 10 to 10,000 $\mu$ L:
  - Finnpiquette F2
  - 10-100 $\mu$ L
  - 100-1000 $\mu$ L
  - 1000-10,000 $\mu$ L
- Flex tips: 200, 1 x 96; 1000, 1 x 96; FT 10mL, 1 x 24 F-stand
- Good Laboratory Pipetting Guide
- Reagent reservoir demo pack

#### Kit 4 includes:

---

- Three Thermo Scientific Finnpiquette Pipetters from 20 to 1000 $\mu$ L:
  - Finnpiquette F2
  - 2-20 $\mu$ L 20-200 $\mu$ L
  - 100-1000 $\mu$ L
- Flex tips: (2) 200 $\mu$ L, 10 x 96; 100-1000 $\mu$ L, 10 x 96
- F-stand
- Good Laboratory Pipetting Guide
- Reagent reservoir demo pack

#### Ordering Information:

The kit contains samples of Finntip Flex, the most advanced traditional tip offering very low attachment and ejection forces. The Good Laboratory Pipetting Guide contains all the information of pipetting from choosing the right type to detailed de-contamination information.

**Warranty:** 5 years with web registration.

