

Thermo Scientific™

Finnpipette™ F2 Variable Volume Pipettes



Related applications:

[Pipettes and Pipette Tips](#)

Thermo Scientific™ Finnpipette™ F2 with a sleek, rugged design is ideal for highly-demanding applications where robustness and durability are essential. Built with tough PVDF components that stand up for harsh chemicals and the damaging effects of UV light. Fully autoclavable.

[Contact us for support >](#)

Catalog number	Unit Size	Description	Volume (Metric)	Compatible Tips	Price (NGN)
4642090	Each	Finnpipette F2	100 to 1000 µL	Flex 1000, Univ. 1000, 1000 Ext	Contact Us >

Full specifications ▼

Compatible Tips	Flex 1000, Univ. 1000, 1000 Ext
Description	Finnpipette F2
Increments	1µL
Color	Blue
Volume (Metric)	100 to 1000 µL
Unit Size	Each

4642010	Each	Finnpipette F2	0.2 to 2 µL	Flex 10, 10, 50	Contact Us >
---------	------	----------------	-------------	-----------------	---------------------------------

Full specifications ▼

Compatible Tips	Flex 10, 10, 50
------------------------	-----------------

Catalog number	Unit Size	Description	Volume (Metric)	Compatible Tips	Price (NGN)
----------------	-----------	-------------	-----------------	-----------------	--------------

Description	Finnpipette F2
Increments	0.002μL
Color	Pink
Volume (Metric)	0.2 to 2 μL
Unit Size	Each

4642020	Each	Finnpipette F2	0.5 to 5 μL	Flex 10, 10, 50	Contact Us >
---------	------	----------------	-------------	-----------------	---------------------------------

Full specifications ▼

Compatible Tips	Flex 10, 10, 50
Description	Finnpipette F2
Increments	0.01μL
Color	Pink
Volume (Metric)	0.5 to 5 μL

Catalog number	Unit Size	Description	Volume (Metric)	Compatible Tips	Price (NGN)
----------------	-----------	-------------	-----------------	-----------------	--------------

Unit Size	Each
------------------	------

4642030	Each	Finnpipette F2	1 to 10 µL	Flex 10, 10, 50	Contact Us >
---------	------	----------------	------------	-----------------	---------------------------------

Full specifications ▼

Compatible Tips	Flex 10, 10, 50
Description	Finnpipette F2
Increments	0.02µL
Color	Pink
Volume (Metric)	1 to 10 µL
Unit Size	Each

4642040	Each	Finnpipette F2	1 to 10 µL	Flex 200, 250 Univ.	Contact Us >
---------	------	----------------	------------	---------------------	---------------------------------

Full specifications ▼

Compatible Tips	Flex 200, 250 Univ.
------------------------	---------------------

Catalog number	Unit Size	Description	Volume (Metric)	Compatible Tips	Price (NGN)
----------------	-----------	-------------	-----------------	-----------------	--------------

Description	Finnpipette F2
Increments	0.02μL
Color	Yellow
Volume (Metric)	1 to 10 μL
Unit Size	Each

4642050	Each	Finnpipette F2	2 to 20 μL	50	Contact Us >
---------	------	----------------	------------	----	---------------------------------

Full specifications ▼

Compatible Tips	50
Description	Finnpipette F2
Increments	0.02μL
Color	Turquoise
Volume (Metric)	2 to 20 μL

Catalog number	Unit Size	Description	Volume (Metric)	Compatible Tips	Price (NGN)
----------------	-----------	-------------	-----------------	-----------------	--------------

Unit Size	Each
------------------	------

4642060	Each	Finnpipette F2	2 to 20 µL	Flex 200, 250 Univ.	Contact Us >
---------	------	----------------	------------	---------------------	---------------------------------

Full specifications ▼

Compatible Tips	Flex 200, 250 Univ.
Description	Finnpipette F2
Increments	0.02µL
Color	Yellow
Volume (Metric)	2 to 20 µL
Unit Size	Each

4642070	Each	Finnpipette F2	10 to 100 µL	Flex 200, 250 Univ., 300, 200 Ext.	Contact Us >
---------	------	----------------	--------------	------------------------------------	---------------------------------

Full specifications ▼

Compatible Tips	Flex 200, 250 Univ., 300, 200 Ext.
------------------------	------------------------------------

Catalog number	Unit Size	Description	Volume (Metric)	Compatible Tips	Price (NGN)
----------------	-----------	-------------	-----------------	-----------------	--------------

Description	Finnpipette F2
Increments	0.2µL
Color	Yellow
Volume (Metric)	10 to 100 µL
Unit Size	Each

4642080	Each	Finnpipette F2	20 to 200 µL	Flex 200, 250 Univ., 300, 200 Ext.	Contact Us >
---------	------	----------------	--------------	------------------------------------	---------------------------------

Full specifications ▼

Compatible Tips	Flex 200, 250 Univ., 300, 200 Ext.
Description	Finnpipette F2
Increments	0.2µL
Color	Yellow
Volume (Metric)	20 to 200 µL

Catalog number	Unit Size	Description	Volume (Metric)	Compatible Tips	Price (NGN)
----------------	-----------	-------------	-----------------	-----------------	--------------

Unit Size	Each
------------------	------

4642100	Each	Finnpipette F2	0.5 to 5 mL	5mL	Contact Us >
---------	------	----------------	-------------	-----	---------------------------------

Full specifications ▼

Compatible Tips	5mL
Description	Finnpipette F2
Increments	0.01mL
Color	Green
Volume (Metric)	0.5 to 5 mL
Unit Size	Each

4642110	Each	Finnpipette F2	1 to 10 mL	10mL, Flex 10 Ext.	Contact Us >
---------	------	----------------	------------	--------------------	---------------------------------

Full specifications ▼

Compatible Tips	10mL, Flex 10 Ext.
------------------------	--------------------

Catalog number	Unit Size	Description	Volume (Metric)	Compatible Tips	Price (NGN)
----------------	-----------	-------------	-----------------	-----------------	--------------

Description	Finnpipette F2
Increments	0.02mL
Color	Red
Volume (Metric)	1 to 10 mL
Unit Size	Each

4642120	Each	Finnpipette F2	5 to 50 µL	Finntip 50µL	Contact Us >
---------	------	----------------	------------	--------------	---------------------------------

Full specifications ▼

Compatible Tips	Finntip 50µL
Description	Finnpipette F2
Increments	0.1µL
Color	Turquoise
Volume (Metric)	5 to 50 µL

Catalog number	Unit Size	Description	Volume (Metric)	Compatible Tips	Price (NGN)
----------------	-----------	-------------	-----------------	-----------------	--------------

Unit Size	Each
------------------	------

4642130	Each	Finnpipette F2	5 to 50 µL	Finntip Flex 200, Finntip Flex 300, Finntip 250 universal, Finntip 200 Extended	Contact Us >
---------	------	----------------	------------	---	---------------------------------

Full specifications ▼

Compatible Tips	Finntip Flex 200, Finntip Flex 300, Finntip 250 universal, Finntip 200 Extended
Description	Finnpipette F2
Increments	0.1µL
Color	Yellow
Volume (Metric)	5 to 50 µL
Unit Size	Each