

Nunc™ Serological Pipettes



Related applications:

Pipettes and Pipette Tips

Stringent quality specifications, and convenient features make Thermo Scientific™ Nunc™ Serological Pipettes a great choice for daily lab work. Bold, clear bi-directional graduations make it easy-to-read fluid volumes and dispensing amounts. Choose from a broad selection of sizes and packaging options.

[Contact us for support ›](#)

Catalog number	Unit Size	Volume (Metric)	Packaging Type	Price (NGN)
170372N	Case of 500	2 mL	Bulk packed, 50/bag, plugged	Contact Us ›

Full specifications ►

170353N	Case of 1000	1 mL	Individually Wrapped, paper/plastic peel packaging, plugged	Contact Us ›
---------	--------------	------	---	------------------------------

Catalog number	Unit Size	Volume (Metric)	Packaging Type	Price (NGN)
Full specifications ►				
170354N	Case of 500	2 mL	Individually Wrapped, paper/plastic peel packaging, plugged	Contact Us >
Full specifications ►				
170355N	Case of 200	5 mL	Individually Wrapped, paper/plastic peel packaging, plugged	Contact Us >
Full specifications ►				
170356N	Case of 200	10 mL	Individually Wrapped, paper/plastic peel packaging, plugged	Contact Us >
Full specifications ►				
170357N	Case of 200	25 mL	Individually Wrapped, paper/plastic peel packaging, plugged	Contact Us >
Full specifications ►				
170358N	Case of 100	50 mL	Individually Wrapped, paper/plastic peel packaging, plugged	Contact Us >
Full specifications ►				

Catalog number	Unit Size	Volume (Metric)	Packaging Type	Price (NGN)
170364N	Case of 1000	1 mL	Individually Wrapped, plastic film packaging, plugged	Contact Us >

Full specifications ►

170365N	Case of 500	2 mL	Individually Wrapped, plastic film packaging, plugged	Contact Us >
---------	-------------	------	---	---------------------------------

Full specifications ►

170366N	Case of 200	5 mL	Individually Wrapped, plastic film packaging, plugged	Contact Us >
---------	-------------	------	---	---------------------------------

Full specifications ►

Show more 10 of 18 | Show more >

Save to list

Product overview

Back to top

Product Details:

- Wide Range of Volumes and Packaging Option
- Color-coded pipettes and packaging for convenient volume selection
- Individually packed pipettes in either paper/plastic or plastic peel for high sterility assurance
- Reduce packaging waste with bulk packed options

Certified Calibration of Accuracy and Precise Graduations

- Ensures dispensing accuracy +/-2% for 1mL and 2 mL and +/-1% for 5 mL, 10 mL, 25 mL, 50 mL at full volume

Etcon Analytical and Environmental Systems & Services Limited

Office Address: No 7, Abeokuta Street, Off Murtala Mohammed Way, Ebutte-Metta, Yaba, Lagos, Nigeria.

Tel: +234-1-4536335, +234-8021422237 | E-mail: info@etconanalytical.com. | Website: www.etconanalytical.com

Clear, Bold, Bi-directional Graduations

- High contrast graduations for easy-to-read volumes
- Reverse graduations aid in determining dispensing volume

Negative Graduations

- Increase pipette capacity, and accommodates additional volumes

Filter Plug

- Prevent aerosol or liquid contamination of the pipetting device while further reducing sample-to-sample contamination risks

Specifications:

- **Material:** Ultra clear, 100% virgin polystyrene for maximum clarity
- **Contaminant Free:** RNase/DNase free, Human DNA free, Non-Cytotoxic, Non-Pyrogenic
- **Sterility:** E-Beam irradiation to SAL of 10^{-6}
- **Shelf Life:** 5 years from Manufacturing Date

Product Testing

- Certificate of Quality