

WELLWASH 1 X 8



	Catalog number	Unit Size	Model	Price (NGN)
	5165000	Each	Wellwash 1x8	Request A Quote

Full specifications ▼

Plate Type	96-Well Plates
Bottles	2L wash bottle, 2L waste bottle
Wash/Waste Bottle Capacity	1
Wash Heads	1x8
Wash Volume	50 to 1000 μ L
Residual Volume	<2 μ L
Dispense Volume	50 to 400 μ L
Accuracy	< 5% @ 300 μ L
Precision	< 3% (CV) @ 300 μ L (96-well plate)

	Catalog number	Unit Size	Model	Price (NGN)
Cell Washing		No		
Prime Volume		5 to 100mL		
Display		LCD		
USB Connections		Yes		
Automation Interface		No		
Width (Metric)		385mm		
Width (English)		15.2 in.		
Depth (Metric)		240mm		
Depth (English)		9.4 in.		
Height (Metric)		345cm		
Height (English)		13.6 in.		
Weight (Metric)		8kg		
Weight (English)		17.6 lb.		

Catalog number	Unit Size	Model	Price (NGN)
Includes	1x2L wash bottle, 1x2L waste bottle, 1x8L wash head, and aerosol cover.		
Hertz	50/60Hz		
Voltage	100-240V		
Certifications/Compliance	CE Marked		
For Use With (Application)	ELISA-Based Assays		
Model	Wellwash 1x8		
Electrical Requirements	100-240V 50/60Hz		
Unit Size	Each		

