



More information on the website
radwag.com/en/info,w1,9XE

100 μ l ÷ 1000 μ l - RW8-106-20-9 Automatic Pipettes, 2 μ l ÷ 20 μ l - RW8-102-20-9 Automatic Pipettes, 10 μ l ÷ 100 μ l - RW8-104-20-9 Automatic Pipettes, 20 μ l ÷ 200 μ l - RW8-105-20-9 Automatic Pipettes, 5 μ l ÷ 50 μ l - RW8-103-20-9 Automatic Pipettes, 0,5 μ l ÷ 10 μ l - RW8-101-20-9 Automatic Pipettes, 500 μ l ÷ 5000 μ l - RW8-107-20-9 Automatic Pipettes



100 μ l ÷ 1000 μ l - RW8-106-20-9 Automatic Pipettes



2 μ l ÷ 20 μ l - RW8-102-20-9 Automatic Pipettes



10 μ l ÷ 100 μ l - RW8-104-20-9 Automatic Pipettes



20 μ l ÷ 200 μ l - RW8-105-20-9 Automatic Pipettes



5 μ l ÷ 50 μ l - RW8-103-20-9 Automatic Pipettes



0,5 μ l ÷ 10 μ l - RW8-101-20-9 Automatic Pipettes



500 μ l ÷ 5000 μ l - RW8-107-20-9 Automatic Pipettes

Datasheet

	100 µl ÷ 1000 µl - RW8-106-20-9 Automatic Pipettes	2 µl ÷ 20 µl - RW8-102-20-9 Automatic Pipettes	10 µl ÷ 100 µl - RW8-104-20-9 Automatic Pipettes
Metrological parameters			
Volume	100 µl ÷ 1000 µl	2 µl ÷ 20 µl	10 µl ÷ 100 µl
Systematic error	2,0% / 0,6%	3,0% / 0,9%	3,0% / 0,8%
Random error	0,7% / 0,2%	2,0% / 0,4%	1,5% / 0,15%

Datasheet

	20 µl ÷ 200 µl - RW8-105-20-9 Automatic Pipettes	5 µl ÷ 50 µl - RW8-103-20-9 Automatic Pipettes	0,5 µl ÷ 10 µl - RW8-101-20-9 Automatic Pipettes
Metrological parameters			
Volume	20 µl ÷ 200 µl	5 µl ÷ 50 µl	0,5 µl ÷ 10 µl
Systematic error	2,0% / 0,6%	2,0% / 0,6%	2,5% / 1,0%
Random error	0,8% / 0,15%	2,0% / 0,3%	1,5% / 0,8%

Datasheet

	500 µl ÷ 5000 µl - RW8-107-20-9 Automatic Pipettes
Metrological parameters	
Volume	500 µl ÷ 5000 µl
Systematic error	2,0% / 0,5%
Random error	0,6% / 0,15%



Compatible with

XA 4Y PLUS Analytical Balance
MYA 4Y Microbalance
Microbalance for Pipette Calibration MYA 4Y.P

XA 4Y.A PLUS Analytical Balance
MYA 4Y PLUS Microbalance
Microbalance for Pipette Calibration MYA 4Y.P PLUS