



More information on the website
radwag.com/en/info,w1,CLW

R-Pipettes



Overview

R-PIPETTES PC software is a RADWAG proprietary solution for determination of errors of pipette volume measurement (batchers) in accordance with the standard ISO 8655. The software supports calibration of fixed volume pipettes (single channel) and adjustable volume pipettes (1, 8, 12 channel).

The program features database of RADWAG pipettes and enables defining custom pipette models.

Compatible with

XA 4Y PLUS Analytical Balance
 MYA 4Y Microbalance
 XA 4Y.M PLUS Microbalance
 Microbalance for Pipette Calibration XA 4Y.M.A.P PLUS
 Microbalance for Pipette Calibration MYA 4Y.P

XA 4Y.M.A PLUS Microbalance
 XA 4Y.A PLUS Analytical Balance
 MYA 4Y PLUS Microbalance
 Microbalance for Pipette Calibration MYA 4Y.P PLUS