



WHY WAIT for test results?



FRENDS™ System



FAST



EASY TO USE



ACCURATE

tPSA

Rapid quantitative assay
for total Prostate Specific Antigen

- Improve patient care
- Increase patient satisfaction
- Increase practice efficiency
- Increase practice profitability


FREND™ System

Accurate, automated testing procedure – simple 3-step operation
No need to batch your sample – run when you are ready

- Controls: Minimum once a month two level liquid controls
- Electronic code chip calibration – no need to run liquid
- Daily electronic system QC control – takes only 1.5 minutes
- Stores 500 Test results



Test procedure




FREND™ PSA Plus

Specification

Analytical measurement range	0.1 ~ 25 ng/mL
Sample type	Plasma or Serum
Sample volume	30 µL
Limit of detection	0.1 ng/mL
Reaction time	Max 4 min.
1 Box	25 tests
Storage condition	35 °F~46 °F

* The specifications are subject to change without prior notice.

Correlation graph



The R² correlation value of FREND™ PSA plus is more than 0.96 against laboratory results.

Ordering information

FREND™ System	Cat# F10
FREND™ PSA Plus	Cat# FRPS 025

Reimbursement (Medicare)

Test	PSA
CPT Code	84153
Reimbursement	\$ 25.09



NanoEnTek



sales@nanoentek.com
www.nanoentekinc.com

NanoEnTek USA Inc.
5627 Stoneridge Drive Suite 304, Pleasanton, CA 94588, USA
Tel : 925-225-0108, 888-988-0108(Toll free) Fax : 925-225-0109