


FDA 510K Cleared



Testosterone

Immediate Consultation Testosterone



FAST



EASY TO USE



ACCURATE



FREND™ System

Quantitative immunoassay analyzer
for Testosterone hormone

- Improve men's healthcare
- Increase customer satisfaction
- Increase practice efficiency
- Maximize your profit
- Test range 20 – 1,500 ng/dL

NanoEnTek

FREND™ System

Accelerating decision-making

Accurate and fully automated FREND™ System brings fast results. With a simple 3-step operation and a wide range of disease markers, the system supports quick decision-making.




Simple and Fast testing



Once a month liquid quality control



User-friendly Icon



6 seconds code chip calibration

Test procedures



STEP 01
Add sample



STEP 02
Insert AP tube



STEP 03
Get Result




Testosterone

Specifications

Assay method	Fluorescent immunoassay
Analytical measurement range	20 ng/dL ~ 1,500 ng/dL
Sample type	Serum or Plasma (Lithium heparin or EDTA)
Sample volume	70 µL
Limit of detection	20 ng/dL
Reaction time	Cartridge reaction: < 3 min.
Pre-incubation	AP (5 minutes)
Package	20 tests/kit
Storage condition	35 ~ 46 °F (2 ~ 8 °C)
Shelf life	18 months

Correlation



Ordering information

FREND™ System	Cat. no: F10
FREND™ Testosterone	Cat. no: FRTEAP 020
CPT Code	84403

