

Vortex Mixer 24



*Versatility
meets
Durability*



OMNI
International
The Homogenizer Company®

Vortex Mixer 24

From vigorous vortexing to gentle mixing, the Vortex Mixer 24 is designed to conveniently handle samples for a variety of applications.

Specifications

Processing Range: 2mL to 50 mL

Power Rating: 120/230 VAC/50/60 Hz

Dimensions: 6.5"D x 4.8"W x 6.5"H

Weight: 8.8 lbs. (4 kg)

Accessories Included: 2

Speed: 600 - 3200 RPM

Orbit: 4mm



ORDER NO. 28-101(115V)
28-102 (220V)



3" Mixing Platform



Mixing Cup

Features

- Mix up to 24 samples in seconds
- Variable speed control
- Hands Free and Touch modes

Optional 2 ml Plus Tube Carriage (sold separately)

- This unique, dual purpose tube carriage allows for vortexing tubes in the center of the device without having to remove the carriage.
- Broad performance range - gentle mixing to high force vortexing.
- Holds 24 - 2 ml tubes

*For best results, use with
Omni International
pre-filled bead kits.*



935C Cobb Place Blvd.
Kennesaw, GA 30144
800.776.4431 • 770.421.0058
www.omni-inc.com



REV.B