Microscope Pen - 330004

Place the pen's angled point perpendicular to the object and look through the lens. Slightly tilt the pen as needed for optimum focus. The object appears inverted.

For best results, use in well-lighted areas such as under a desk lamp. Useful for spot inspections in printing, textile, art, microelectronics, jewelry and other industries.

Also serves as an easy-to-use, inexpensive microscope for science education and in labs.

Clips to a pocket and comes in a protective pouch.

Powerful 50X optics Weight 0.7 oz (20 g) Dim: 51/8" x 1/2" (135x10 mm)

WARRANTY

Sper Scientific warrants this product against defects in materials and workmanship for a period of **five (5) years** from the date of purchase, and agrees to repair or replace any defective unit without charge. If your model has since been discontinued, an equivalent Sper Scientific product will be substituted if available. This warranty does not cover probes, batteries, battery leakage, or damage resulting from accident, tampering, misuse, or abuse of the product. Opening the meter to expose its electronics will break the waterproof seal and void the warranty. To obtain warranty service, ship the unit postage prepaid to:

SPER SCIENTIFIC LTD.

8281 E. Evans Rd., Suite #103 Scottsdale, AZ 85260 (480) 948-4448

The defective unit must be accompanied by a description of the problem and your return address. Register your product online at www.sperwarranty.com within 10 days of purchase.