

MULTI-TX5 Digital

Multi-tube vortex mixer

Versatile, compact and easy to use multi-tube vortex mixer for high sample throughput



Multi-tube Vortex Mixer

The MULTI-TX5 Digital is the versatile and compact vortex mixer for high sample throughput. Vigorous as well as gentle mixing is supported with premium performance making the MULTI-TX5 Digital ideal both for microbiology tasks and R&D or quality control in diverse industries like pharma, cosmetics, food & beverage, and chemicals.

- Stirring speed from 100 up to 2500 rpm thanks to the powerful brushless motor
- Easy to change foam racks for various test tubes
- Structure designed for outstanding comfort, maximized resistance, and minimized footprint
- Easy to operate digital interface with icons
- The built-in timer allows programming of unsupervised operation
- Programmable auto-reverse of the stirring direction
- Exclusive Pulse Mode enabling cyclical vigorous mixing
- The Lock feature prevents unintentional changes



INSTRUMENT - CODE

MULTI-TX5 Digital Vortex

100-240 V / 50-60 Hz

F202A0460



TECHNICAL DATA

USER INTERFACE	Digital
SPEED CONTROL	Electronic
STIRRING SPEED	100 - 2500 rpm
SPEED SETTING ACCURACY	5 rpm
TIMER	From 1 min to 99 h : 59 min
TIMER SETTING STEPS	1 min
AUTOREVERSE	From 5 s to 99 min : 59 s
AUTOREVERSE SETTING STEPS	1s
INTERMITTENT MODE	From 5 s to 99 min : 59 s
INTERMITTENT MODE SETTING STEPS	1s
PULSE MODE	2 s shaking, 1 s pause (cyclic)
LOCK FUNCTION	Yes
MOTOR	BLDC
TYPE OF MOVEMENT	Orbital
ORBITAL DIAMETER	3.6 mm
WORKING IN CONTINUOUS	Max 2 hours
POWER INPUT	40 W
WEIGHT	19 kg; 42 lb
SUPPORT PLATE DIMENSION	294 x 194 mm; 11.6 x 7.6 in
DIMENSIONS (WXHXD)	245 x 380 x 490 mm; 9.6 x 14.9 x 19.3 in



VELP Official Partner

We reserve the right to make technical alterations We do not assume liability for errors in printing, typing or transmission

DESIGNED AND MANUFACTURED IN ITALY



Via Stazione 16 20865 Usmate (MB) Italy Tel. +39 039 628811 velpitalia@velp.com

155 Keyland Court, Bohemia NY 11716 - U.S. Tel. +1 631 573 6002 velpusa@velp.com

USA

CHINA

Xinlong Rd Building 28, Lane 1333 Shanghai city - China Tel. +86 18017557329 velpchina@velp.com



INDIA velpindia@velp.com

ITALY - HQ