

Lovibond[®] Colour Measurement Tintometer[®] Group



Saybolt Color +25

Part Number: 134070



Reference liquid standards are assigned a nominal value and Lovibond Tintometer always endeavors to match this value as closely as possible. On occasion, it may not be possible to match the requested value exactly. The value achieved and the expected performance tolerances will, of course, always be reported.

Saybolt Color +25

500ml Bottle certified to ISO/IEC 17025 with a 12-month shelf life.
Full traceability to internationally recognized standards
Each bottle supplied with full certification including MSDS