

Lovibond[®] Colour Measurement Tintometer[®] Group



EBC (CIE) Colour Filter Set

Part Number: 139400



Conformance filters are assigned a nominal value and the company always endeavours to match this value as closely as possible. On occasion, because of the unique manner of the glass making process, it may not be possible to match the requested values exactly. The values achieved and the performance tolerances expected will, of course, always be reported.

EBC (CIE) Colour Filter Set

A set of five glass conformance filters for checking PFXi instruments in EBC (CIE) scale.

Technical Data

Glass Standard Values	Approximate values
	- 4
	- 9.5
	- 15
	- 18
	- 25

Filters 5

Delivery Scope Text

- EBC 4 filter
- EBC 9.5 filter
- EBC 15 filter
- EBC 18 filter
- EBC 25 filter
- Certificate of Conformity