

Lovibond® Colour Measurement Tintometer® Group



Gardner Colour Filter Set (PFXi995/950/880)

Part Number: 139890



Conformance filters are assigned a nominal value and the company always endeavours to match this value as closely as possible. On occasion, because of the unique manner of the glass making process, it may not be possible to match the requested values exactly. The values achieved and the performance tolerances expected will, of course, always be reported.

Gardner Colour Filter Set (PFXi995/950/880)

A set of four glass conformance filters for checking PFXi instruments in Gardner scale.

Technical Data

Glass Standard Values	Approximate Values
	- 2
	- 8
	- 12
	- 17

Filters 4

Delivery Scope Text

- Gardner 2 filter
- Gardner 8 filter
- Gardner 12 filter
- Gardner 17 filter
- Certificate of conformity