

# Lovibond® Colour Measurement Tintometer® Group



## FAC Colour

Part Number: 343000



Conformance filters are assigned a nominal value and the company always endeavours to match this value as closely as possible. On occasion, because of the unique manner of the glass making process, it may not be possible to match the requested values exactly. The values achieved and the performance tolerances expected will, of course, always be reported.

### FAC Colour

A set of four glass conformance filters for checking the 3000 Comparator in FAC scale

### Technical Data

Glass Standard Values	Approximate Value
	- 7
	- 13
	- 29
	- 39

Filters 4

### Delivery Scope Text

- FAC 7 filter
- FAC 13 filter
- FAC 29 filter
- FAC 39 filter
- Certificate of conformity