

# Lovibond® Colour Measurement Tintometer® Group



## ASTM Color - EC3000

Part Number: 135079



Conformance filters are assigned a nominal value and the company always endeavors to match this value as closely as possible. On occasion, because of the unique manner of the glass making process, it may not be possible to match the requested values exactly. The values achieved and the performance tolerances expected will, of course, always be reported.

### ASTM Color - EC3000

A single glass filter for checking the PFXi instruments in ASTM D1500 scale. This filter was supplied as part of the PFX195/2 accessories kit.

### Technical Data

<b>Glass Standard Values</b>	- 1.1
<b>Filters</b>	1

### Delivery Scope

- Filter
- Certificate of Conformity