

Lovibond[®] Colour Measurement Tintometer[®] Group



Gardner Colour 5

Part Number: 134210



Reference liquid standards are assigned a nominal value and Lovibond Tintometer always endeavours to match this value as closely as possible. On occasion, it may not be possible to match the requested value exactly. The value achieved and the expected performance tolerances will, of course, always be reported.

Gardner Colour 5

500ml Bottle certified to ISO/IEC 17025 with a 12-month shelf life.
Full traceability to internationally recognised standards
Each bottle supplied with full certification including MSDS