

# Haier

## HUMANIZED DESIGN ACCURATE MEASUREMENT

MODEL  
**B15**



ELECTRONIC BLOOD PRESSURE MONITOR

## HUMANIZED DESIGN

### HUMANIZED STYLE

Adopt Nordic style design  
combine home design style with simple and  
practical functional characteristics



## LARGE SCREEN DISPLAY EASY TO READ

The screen of sphygmomanometer reaches 3.5 inches  
and  
the large font design is convenient for users to read



**3.5**  
INCHES

## SCIENTIFIC ALGORITHM ACCURATE MEASUREMENT



The built-in MCU chip bring accurate results  
and comfortable experience

## EASY TO USE EASY TO MEASURE

Humanized design ensures  
that the measurement  
process is simple.

the user is easy to use,  
and easy to read  
the measurement results.



**A**



**B**

## ARRAY MEMORY EASY TO QUERY

It supports array memory,  
and stores two user measurement data  
into different memory groups,  
which is convenient for users to query

## PRODUCT FEATURES

- 1 Memory: 2x120 groups (double users)
- 2 Average function: Last 3 groups average measuring value
- 3 WHO function
- 4 Arrhythmia checking:  
Monitor arrhythmia while measuring the blood pressure.
- 5 Material: ABS+LCD display

### ITEM

Product Size:140\*120\*67mm  
Product Weight:292.5g  
Memory:2\*120

Measurement range  
280mmHg(0-37.30kPa)  
Pulse:40-199 beats/min

Accuracy:Pressure  $\pm 3$ mmHg(0.4kPa) Pulse  $\pm 5\%$  of the reading  
Automatic power-off:60s