

# Lovibond® Colour Measurement Tintometer® Group



## Hazen Colour Filter Set - EC 2000

Part Number: 135159



Conformance filters are assigned a nominal value and the company always endeavours to match this value as closely as possible. On occasion, because of the unique manner of the glass making process, it may not be possible to match the requested values exactly. The values achieved and the performance tolerances expected will, of course, always be reported.

### Hazen Colour Filter Set - EC 2000

A set of five glass conformance filters for checking Pt-Co scale on the PFXi 880, 950 and 995 instruments.

DO NOT USE on the PFXi 195!

### Technical Data

Glass Standard Values	Approximate values
	- 5
	- 20
	- 50
	- 100
	- 300

---

Filters	5
---------	---

### Delivery Scope

- 5, 20, 50, 100, 300 filters
- Certificate of Conformity