

Haier Biomedical

Intelligent Protection of Life Science

CP101 Product Parameters



Product Model	CP101
Flow (L/min)	0-10
Purity (L/min)	0-10, 93±3%
Noise (dB (A))	≤55
Voltage Frequency (V/Hz)	220/50
Power Consumption (VA)	700
Dimension (mm)	400*300*540
Weight (Kg)	20.5
Display Mode	LED
Carton Size (mm)	470*370*590
Packing Weight(Kg)	22.5
Certification	CE


CP Series Product Features



- ✓ Multiple intake filtration
- ✓ Stylish design



- ✓ Warning Function
- ✓ High definition display



- ✓ One-key start
- ✓ Light-touch switch



- ✓ Low-noise design
- ✓ Fast and efficient oxygen production



Qingdao Haier Biomedical Co.,Ltd.

No.280 Feng Yuan Road, High-tech Zone,
Qingdao, 266109, P.R. China
Tel: +86-0532-88935593
Website: www.haiermedical.com



Haier Biomedical International



Haier Biomedical International



@haiermedicalint



Haier Biomedical International



Haier Biomedical International