

# Haier

## HUMANIZED DESIGN ACCURATE MEASUREMENT

MODEL  
U80H



ELECTRONIC BLOOD PRESSURE MONITOR

## HUMANIZED DESIGN HUMANIZED STYLE


Adopt Nordic style design  
combine home design style with simple and  
practical functional characteristics



## LARGE SCREEN DISPLAY EASY TO READ

The screen of sphygmomanometer reaches 3.5 inches  
and  
the large font design is convenient for users to read

3.5  
INCHES



## SCIENTIFIC ALGORITHM ACCURATE MEASUREMENT

The built-in MCU chip bring accurate results  
and comfortable experience

## EASY TO USE EASY TO MEASURE

Humanized design ensures  
that the measurement  
process is simple.

the user is easy to use,  
and easy to read  
the measurement results.



## ARRAY MEMORY EASY TO QUERY

It supports array memory,  
and stores two user measurement data  
into different memory groups,  
which is convenient for users to query

### ITEM

Product Size:103\*103\*65mm

Product Weight:230g

Memory:2\*90

Measurement range

0-299 mmHg(0-37.30kPa)

Pulse:40-199 beats/min

Accuracy:Pressure  $\pm 3$ mmHg(0.4kPa) Pulse  $\pm 5\%$  of the reading

Automatic power-off:180s