

Lovibond® Colour Measurement Tintometer® Group



IP Units Filter Set

Part Number: 137660



Conformance filters are assigned a nominal value and the company always endeavours to match this value as closely as possible. On occasion, because of the unique manner of the glass making process, it may not be possible to match the requested values exactly. The values achieved and the performance tolerances expected will, of course, always be reported.

IP Units Filter Set

A set of four glass conformance filters for checking the PFXi instruments in IP (method B) units

Technical Data

Glass Standard Values	Approximate values
	- 0.75
	- 1.25
	- 2.0
	- 3.5

Filters 4

Delivery Scope Text

- IP 0.75 filter
- IP 1.25 filter
- IP 2.0 filter
- IP 3.5 filter
- Certificate of conformity