

Lovibond® Colour Measurement Tintometer® Group



Honey Colour (Pfund Equivalents) Filter Set

Part Number: 139370



Conformance filters are assigned a nominal value and the company always endeavours to match this value as closely as possible. On occasion, because of the unique manner of the glass making process, it may not be possible to match the requested values exactly. The values achieved and the performance tolerances expected will, of course, always be reported.

Honey Colour (Pfund Equivalents) Filter Set

A set of five glass conformance filters for checking PFXi instruments in Pfund honey scale.

Technical Data

Glass Standard Values	Approximate values
	- 15
	- 30
	- 60
	- 85
	- 100

Filters 5

Delivery Scope Text

- Pfund 15mm filter
- Pfund 30mm filter
- Pfund 60mm filter
- Pfund 85mm filter
- Pfund 100mm filter
- Certificate of conformity