

Lovibond® Colour Measurement Tintometer® Group



Saybolt Colour Filter Set (only PFX/PFXi195)

Part Number: 139390



Conformance filters are assigned a nominal value and the company always endeavours to match this value as closely as possible. On occasion, because of the unique manner of the glass making process, it may not be possible to match the requested values exactly. The values achieved and the performance tolerances expected will, of course, always be reported.

Saybolt Colour Filter Set (only PFX/PFXi195)

A set of five glass conformance filters for checking the PFXi 195 instruments in Saybolt scale. NOT FOR USE with the PFXi 880, 850 or 995.

Technical Data

Glass Standard Values	- -8
	- 0
	- +10
	- +18
	- +25

Filters 5

Delivery Scope Text

- Saybolt -8 filter
- Saybolt 0 filter
- Saybolt +10 filter
- Saybolt +18 filter
- Saybolt +25 filter
- Certificate of conformity