

Lovibond® Colour Measurement Tintometer® Group



Saybolt (only PFX/PFXi880/950/995)

Part Number: 139583



Conformance filters are assigned a nominal value and the company always endeavours to match this value as closely as possible. On occasion, because of the unique manner of the glass making process, it may not be possible to match the requested values exactly. The values achieved and the performance tolerances expected will, of course, always be reported.

Saybolt (only PFX/PFXi880/950/995)

A single glass conformance filter for checking the PFXi 880, 950 or 995 in Saybolt scale.

Technical Data

| Glass Standard Values | Approximate Value: |
|-----------------------|--------------------|
| | - +10 |

| | |
|---------|---|
| Filters | 1 |
|---------|---|

Delivery Scope Text

- Saybolt +10 filter
- Certificate of conformity