

Lovibond® Colour Measurement Tintometer® Group



Klett (blue filter KS-42) Filter Set (PFXi195)

Only for use with the PFXi 195

Part Number: 139860



Conformance filters are assigned a nominal value and the company always endeavours to match this value as closely as possible. On occasion, because of the unique manner of the glass making process, it may not be possible to match the requested values exactly. The values achieved and the performance tolerances expected will, of course, always be reported.

Klett (blue filter KS-42) Filter Set (PFXi195)

A set of five glass conformance filters for checking Klett scale on PFXi instruments.

Technical Data

Glass Standard Values	Approximate Values
	- 21
	- 66
	- 162
	- 318
	- 616

Filters 5

Delivery Scope Text

- Klett 21 filter
- Klett 66 filter
- Klett 162 filter
- Klett 318 filter
- Klett 616 filter
- Certificate of conformity