



Pharmacy Refrigerator

HYC-310

Core Advantages

- Superior alarm system includes high and low temperature alarm, sensor error alarm and power failure alarm
- Air cooling system, temperature uniformity inside the cabinet is $\pm 3^{\circ}\text{C}$, automatic defrosting
- HC refrigerant, energy-saving and environmentally friendly
- Precise temperature control of $2\sim 8^{\circ}\text{C}$

Innovative & Ergonomic Design

- Automatic evaporation of condensed water without manual draining
- Anti-condensation design of front blower
- Double lock design, fitted lock and padlock, safe and convenient
- Adjustable multilayer shelf design

Qingdao Haier Biomedical Co.,Ltd.

No.280 Feng Yuan Road, High-tech Zone,
Qingdao, 266109, P.R. China
Tel: +86-0532-88935593
Website: www.haiermedical.com



Typical Application

The refrigerator is designed for storage of bio products, vaccines, drugs, reagents and is applicable for pharmacies, pharmaceutical companies, hospitals, CDCs and clinics.

Product Advantages



Control System

- Microprocessor control, LED digital display of temperatures, symbol display of working status. The temperature accuracy is 0.1°C.
- Optimized air distribution system designed to guarantee temperature range from 2°C to 8°C.



Refrigeration System

- High quality compressor and fan ensure superior refrigeration performance.
- Air cooling system, the temperature uniformity inside the cabinet is $\pm 3^{\circ}\text{C}$, automatic defrosting.
- HC refrigerant, energy-saving and environmentally friendly.



Safety System

- Equipped with a complete temperature alarm system featuring audible buzzer and visual symbols which is capable of alerting failures due to high and low temperature, sensor error and power failure.
- Cabinet temperature will continue to be displayed after power failure.
- Safety door lock design.
- Features multiple password protection to prevent unauthorized access to ensure security.

Innovative and Ergonomic Design



The front of the unit is designed to blow hot air to prevent condensation.



Multiple adjustable shelves design for better space utilization.



Double lock design, fitted lock and padlock, safe and convenient

Parameters

Model	Voltage (V/Hz)	Power (W)	Temp Range (°C)	Exterior Dimensions (mm)(W*D*H)	Interior Dimensions (mm)(W*D*H)	Capacity (L)	Net Weight/Gross Weight (Kg)	Shelves	Porthole
HYC-310	220	260	2-8	600*615*1,950	560*500*1,290	310	77/90	4	1