

## CALIBRATION AND MAINTENANCE SOLUTIONS



STEIEL Elettronica offers a complete range of standard solutions, for the periodic calibration and maintenance of your sensors.

Item	Description	Part No.
pH4-S	pH 4.01 buffer solution, 90 ml bottle	80090095
pH4-L	pH 4.01 buffer solution, 500 ml bottle	80090093
pH7-S	pH 7.01 buffer solution, 90 ml bottle	80090096
pH7-L	pH 7.01 buffer solution, 500 ml bottle	80090094
pH9-S	pH 9.01 buffer solution, 90 ml bottle	80090097
pH9-L	pH 9.01 buffer solution, 500 ml bottle	80090100
RX220-S	ORP calibration solution at 220 mV (@ 25°C), 90 ml bottle	80190091
RX220-L	ORP calibration solution at 220 mV (@ 25°C), 500 ml bottle	80190090
C0147-S	Standard conductivity solution at 147 µS/cm (@ 25°C), 90 ml bottle	80290156
C0147-L	Standard conductivity solution at 147 µS/cm (@ 25°C), 500 ml bottle	80290155
C01413-S	Standard conductivity solution at 1413 µS/cm (@ 25°C), 90 ml bottle	80290154
C01413-L	Standard conductivity solution at 1413 µS/cm (@ 25°C), 500 ml bottle	80290153
C01288-S	Standard conductivity solution at 12.88 mS/cm (@ 25°C), 90 ml bottle	80290152
C01288-L	Standard conductivity solution at 12.88 mS/cm (@ 25°C), 500 ml	80290151
KCL/250	KCl electrolyte for refillable electrodes, 250 ml bottle	80090087
ST/500	Electrode storage solution, 500 ml bottle	80090091
KRE	Solution kit for cleaning and storage of pH and redox electrode	80099902