




Mini ULT Freezer



Application:
It is suitable for biobanks, vaccination production and transportation users.

-  Long-distance transportation: Satisfies the -80°C ULT transportation application requirements of 24V(DC) vehicle power supply.
-  Dry ice not required: The cooler maintains the temperature at -80°C for a long period time. Superior temperature uniformity inside the cabinet.
-  Multiple applications: The cooler can also be used with 220V power supply for long-term storage.

Qingdao Haier Biomedical Co.,Ltd.

No.280 Feng Yuan Road, High-tech Zone,
Qingdao, 266109, P.R. China
Tel: +86-0532-88935593
Website: www.haiermedical.com



Haier Biomedical
International



Haier Biomedical
International



@haiermedicalint



Haier Biomedical
International



Haier Biomedical
International

Product Advantages



HC Energy Saving

High quality HC high-efficiency compressor with optimized refrigeration system, which can save 50% on power consumption.



One Unit for Multiple Applications

*Vehicle-mounted transportation with 24V(DC) power supply
*Long-term storage with 220V power supply.



Superior Uniformity

Built-in evaporator, fast cool down refrigeration, the temperature is safe, secure and reliable with the uniformity of $\pm 3^{\circ}\text{C}$ at key characteristic points.

Ergonomic Design



Lightweight

The whole unit weighs 25kg, easy to carry.



Superior Insulation Performance

25mm VIP + 70mm thick LBA foam insulation layer design, dual seal strips for better insulation.



Safety Lock Catch

Simple adjustable mechanical lock with lock catch, and it can be equipped with an external lock to guarantee the safety of sample transportation and storage.



Metal Handle

Hidden metal handle, stylish and comfortable, safe and durable, easy to carry.

Parameters

| Model | Power Supply | Temp Range (°C) | Exterior Dimensions (W*D*H mm) | Interior Dimensions (W*D*H mm) | Effective Capacity (L) | Net/Gross Weight (Kg) |
|-----------|---------------------|-----------------|--------------------------------|--------------------------------|------------------------|-----------------------|
| DW-80WZ15 | 220V/50Hz DC 24V | -80--40 | 690*320*485 | 260*180*330 | 15 | 25/30 |