IECHNOLOGY FOR FUTURE



AH- 16p

Portable Autoclave, 16 litre



Compact bench top autoclave, hevey duty with hight safety and simple operation. Multipurpose laboratory use at effective price.



- > This autoclave with high performance, reliability and safety, and easy operation.
- WIDE RANGE USE: It is suitable for wrapped / un-wrapped / liquid / utensils ... etc. Widely used in dental clinics, Hospitals, Clinics, and Laboratories.
- > EASY MAINTENANCE :
 - The sterilizers are made by stainless steel 304, and the heater is made by stainless steel 316L.. So, it is easy to keep clean and maintain.
- > The new standard sterilization technology have been achieved in terms of testing, certification and technical assistance.
- Intelligent design sterilization programs, fully automatic and easy control.
- Automatic control from heating-up sterilization to time-off.
- The entire sterilization process is continuously monitored by a microprocessor control system.
- > The sterilization time starts counting, when the chamber pressure(temperature) arrives the point you set-up.(SPECIAL FUNCTION)



AH-16p t name: (Class-N) Bench-ton Type Autoclave

Madal		ALI 45m
Model		AH-16p
Product name:		(Class-N) Bench-top Type Autoclave
Overall Dimensions:		335 X 430 x 600 mm (W x H x D)
Construction:		Chamber & Door: Stainelss steel 304
Gross weight		15,5 kg
Chamber Dimensions:		230 x 410 mm(Ø x D)
Chamber Capacity:		16 Liters
Max. Sterilization range:		134°C / 2,0 bar
Capacity of water tank:		N/A
Wafer consumption per cycle:		250 cc
Power Supply:		230 V, 50-60 Hz
Power Consumptions:		1400 W, 7 A
Certificate:		CE
Features:		
o Fill Water System:		Manually
o Temperature Display:		Pressure / temoperature Gauge
o Pressure Display:		Pressure / temoperature Gauge
o Progrom(Function) Display:		LED Display
o Indicalor Light:		Power, heat, sterilization, complete
o Program for Select1on:		Sterilization tims select 8, 33 min
o Starilization Temp.Selection		122 °C (1,4 kg/cm²)
o Dry Function		N/A
Safety Device:		
o Pressure Salety Valve:		Setup to 1,7 kg/cm ²
o Over-pressure Protection:		Via Pressure Protection Switch
o Over-heated Protection:		Yes
o Electric Over load Protection:		10 A, No Fuse Breaker X 2
o Door Salety holder:		Yes
o Pressure Door Auto-Lock:		N/A
o Emergency exhaust valve:		Yes
Life time		7 Yeasr
		1x Sterilization Box: (LxWxH): 336 x 174 x 118 mm
		1x Heater cover / Material: Stainless Steel 304
Optional Accessories		1x Tray set : (3x Tray, 1x Tray frame, 1x Holder)
		2x Spring Holder (Pouch Rack) Material: Stainless Steel 304
Order Informations:		
0.120.014	Туре АН-16р	Portable Bench-Top Autoclave, 16 L
0.120.050	AH-box 1	SS Sterilization Box: (LxWxH): 336 x 174 x 118 mm for 16 L size
0.120.053	AH-heater 1	SS Heater cover for 16L size
0.120.056	AH-tray 1	SS Tray set : (3x Tray, 1x Tray frame, 1x Holder) for 16L size
0.120.059	AH-basket 1	SS Sterilization Basket for 16L size
0.120.062	AH-holder	SS Spring Holder (Pouch Rack),
0.120.014	Type AH-16	Bench Top Autoclave, microprocessor control, 16 L
0.120.050	AH-box 1	SS Sterilization Box: (LxWxH): 336 x 174 x 118 mm for 16 L size
		22 2.2228011 BON (2.111.11). 000 X 174 X 110 HIII HOF TO E 3120



SS Box



SS Tray set



SS Pouch rack

Address. S C O - T e c h Auf der Heide 15 D-37351 Dingelstädt Federal Republic of Germany Phon: +49 (0) 36075 / 43439-0 Fax: +49 (0) 36075 / 43439-9 E-mail: sco@sco-tech.com Internet: www.sco-tech.com

Contact